

Product datasheet for **MG217100**

Homer2 (NM_001164086) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Homer2 (NM_001164086) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Homer2
Synonyms: 9330120H11Rik; AW539445; CPD; Vesl-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG217100 representing NM_001164086
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAGAACAGCCCATCTTCACCACGCGAGCGCACGTCTCCAGATTGACCCAGCACCAAGAAGAACT
GGGTGCCGGAAGCAAGCAGGCCGTACCGTTTCTACTTCTATGATGTCACCAGGAACAGCTATCGGAT
CATCAGTGTGGATGGAGCCAAGGTGATCATAAACAGCACTATCACCCGAACATGACTTTCACAAAACG
TCACAGAAGTTCGGCAGTGGGCTGACAGCAGAGCCAACACCGTGTTCGGTTTGGGATTCTCCTCCGAGC
TGCAGCTCACGAAGTTTGCAGAGAAGTTCCAGGAGGTAAGAGAAGCTGCCAGGCTAGCCAGAGACAAGTC
CCAGGAGAAAACCGAGACCTCCAGCAATCATTCCAAGCATCCAGCGTCAATGGCACAGACGACGAAAAG
GCCTCTCACGCGAGCCCAGCCGACACTCACCTCAAGTCTGAGAATGACAAGCTGAAGATCGCGCTGACAC
AGAGTGTGCCAATGTGAAGAAGTGGGAGATGGAGCTGCAGACCCCTGCGGAGAGCAACGCCCGGCTGAC
CACGGCACTGCAGGAGTCGGCGGCCAGCGTGGAGCAGTGGAAAGCGGCAAGTTCATCTGCAGGGACGAG
AATGACAGGCTCCGACGCAAGATCGAGGAGCTGGAAGAACAGTGCAGCGAGATAAACAGGGAGAAGGAGA
AGAACACACAGCTGAAGAGGAGGATCGAGGAGCTGGAGTCAGAGGTCGAGACAAGGAGATGGAGTTGAA
AGATCTCCGAAAACAGAGTGAATCATACTCAGCTCATGTCCGAGTGTGAATATGTCTGAGAAGTTA
GAGGCGGCCGAAAGAGACAATCAAACTTGAAGACAAAGTCCGCTCTCTAAAGACAGACATCGAGGAGA
GTAATACCGACAGCGCCACCTGAAGGGGAGCTGAAGAGCTTCCTTGAGGTGCTGGATGGAAGATCGA
CGACCTCCATGACTTCCGTAGAGGACTCTCCAAGTTAGGCACAGATAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG217100 representing NM_001164086
 Red=Cloning site Green=Tags(s)

MGEQPIFTTRAHVFDIDPSTKKNWVPASKQAVTVSYFYDVTRNSYRIISVDGAKVIINSTITPNMTFTKT
 SQKFGQWADSRANTVFGGLFSSELQLTKFAEKQEVREAAARLARDKSQEKTTSSNHSQASSVNGTDDEK
 ASHASPADTHLSENDLKIALLTQSAANVKKWEMELQTLRESNARLTTALQESAASVEQWKRQFSICRDE
 NDRLRSKIEELEEQCSEINREKEKNTQLKRRIEELESEVRDKEMELKDLRKQSEIIPQLMSECEYVSEKL
 EAAERDNQNLKDYRSLKTDIEESKYRQRHLKGELKSFLEVLGDGKIDDLHDFRRLGSKLGTDN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001164086

ORF Size: 1029 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164086.1](#), [NP_001157558.1](#)

RefSeq Size: 10966 bp

RefSeq ORF: 1032 bp


Locus ID: 26557

UniProt ID: [Q9QWW1](#)

Cytogenetics: 7 D3

Gene Summary: Postsynaptic density scaffolding protein. Binds and cross-links cytoplasmic regions of GRM1, GRM5, ITPR1, DNM3, RYR1, RYR2, SHANK1 and SHANK3. By physically linking GRM1 and GRM5 with ER-associated ITPR1 receptors, it aids the coupling of surface receptors to intracellular calcium release. May also couple GRM1 to PI3 kinase through its interaction with AGAP2 (By similarity). Isoforms can be differently regulated and may play an important role in maintaining the plasticity at glutamatergic synapses (By similarity) Required for normal hearing (PubMed:25816005). Negatively regulates T cell activation by inhibiting the calcineurin-NFAT pathway. Acts by competing with calcineurin/PPP3CA for NFAT protein binding, hence preventing NFAT activation by PPP3CA (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG217100