

Product datasheet for **MG217102**

Klrb1a (NM_010737) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klrb1a (NM_010737) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Klrb1a
Synonyms:	ly-55A; Ly55a; NKR-P1.7; NKR-P1A; Nkrp1-a; Nkrp1a; NKRP12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG217102 representing NM_010737 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACACAGCAAGGGTCTACTTTGGTTTAAAGCCACCCAGGACTCCAGGGGCTTGGCATGAGTCACCCC
CATCTCTTCCCCAGATGCCTGTCGGTGCCACGTTACACAGGTTGGCTCTGAAGCTCAGCTGTGCTGG
GCTCATCTCCTTGTCTGACCCTGATTGGGATGAGTGTTCAGTGCAGTCCTAATAACAAAACCATCA
ATAGAAAATGCTATGTGCTTATTCAAGAGAACCTGAATAAACAACAGATTGTTCAAGCTAAGCTAGAGT
GCCACAAGACTGGCTTTCACACCGAGATAAGTGTCTTTCAGTTTCTCAAGTTTCAACACTTGGGAGGA
AGGTCTAGTTGATTGTGATGGAAAAGGAGCCACTTTGATGCTCATTCAAGACCAAGAAGAACTGAGATTC
CTACTGGACTCAATAAAGGAAAAATACAATTCATTTTGGATTGGACTAAGGTACACATTGCCAGACATGA
ACTGGAAGTGGATAAATGGATCGACTTTAAATTCTGATGTATTAATAAATCACTGGTGACTGAAAATGA
CAGCTGTGCTGCTATCTCAGGAGACAAAGTGACTTTTGGAGAGCTGAATTCAGACAACCGTTGGATTTGC
CAAAGGAACTATACCATGAAACCCTGAGCAACTATGTGGGTTATGGACAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG217102 representing NM_010737
 Red=Cloning site Green=Tags(s)

MDTARVYFGLKPPRTPGAWHESPPSLPPDACRCPRSHRLALKLSCAGLILLVVTLLIGMSVLRVLIQKPS
 IEKCYVLIQENLNKTTDCSAKLECPQDWLSHRDKCFHVSQVSNTWEEGLVDCDGKATMLLIQDQEELRF
 LLDSIKEKYNSFWIGLRYTLPDMNWKWINGSTLNSDLVKITGDTENDSCAAISGDKVTFESCNSDNRWIC
 QKELYHETLSNYVGYGH


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_010737

ORF Size: 699 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010737.3 , NP_034867.3
RefSeq Size:	1284 bp
RefSeq ORF:	702 bp
Locus ID:	17057
UniProt ID:	P27811
Cytogenetics:	6 63.06 cM
Gene Summary:	Plays a stimulatory role on natural killer (NK) cell cytotoxicity.[UniProtKB/Swiss-Prot Function]