

Product datasheet for **MG217122**

Hnrnpf (NM_001166432) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hnrnpf (NM_001166432) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Hnrnpf
Synonyms: 4833420I20Rik; AA407306; Hnrpf
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG217122 representing NM_001166432
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATGCTGGGCCCTGAGGGAGGTGAAGGCTATGTGGTCAAACCTCCGTGGCCTACCTGGTCTGCTCAA
 TTGAGGACGTACAAAACCTCCTCTCCGACTGCACAATTCATGATGGGGTGCAGGTGTTCAATTCATTTA
 TACTAGAGAAGGCAGGCAGAGTGGTGGGCTTTTGTGAACTTGAGTCAGAAGATGATGAAAATTGGCT
 CTGAAAAAAGACAGGGAAAGCATGGACACCGGTATATTGAGGTGTTCAAGTCACACAGAACCGAGATGG
 ATTGGGTGTTGAAGCACAGTGGTCCAAACAGCGCCGACAGTGCCAATGATGGCTTTGTGAGGCTTCGGGG
 ACTCCCATTTGGATGCACAAAGGAAGAAATCGTTCAGTTCCTCTCAGGGTTGAAAATTGTGCCAAACGGG
 ATCACACTACCTGTGGACCCGGAAGCAAGATTACAGGGGAGGCCCTTCGTTCAAGTTTGCCTCACAAAGAGT
 TAGCTGAGAAAGCTTTAGGGAAGCACAAGGAGAGAATAGGGCACAGGTATATTGAAGTGTAAAGAGCAG
 TCAGGAGGAAGTTAGATCATACTCAGATCCACCTCTGAAGTTTATGTCTGTGCAAAGGCCTGGGCCCTTAT
 GACAGGCTGGCACAGCCCGGAGGTACATTGGCATTGTGAAACAGGCAGGTCTGGATAGGATGAGGTCTG
 GTGCCTATAGTGCAGGCTATGGGGCTATGAAGAATACAGTGGCCTCAGTGATGGCTATGGCTTACCAC
 TGACCTGTTTGGGAGAGACCTCAGTATTGTCTCTCAGGAATGTATGACCACAGATATGGAGACAGCGAG
 TTCACAGTGCAGAGCACACCAGCCACTGCGTCCACATGAGAGGGCTGCCCTACAAGCAACGGAGAACG
 ACATTTACAACCTTCTCTCCACTCAACCCTGTGAGAGTTCATATTGAGATTGGTCTGATGGAAGAGT
 GACGGGAGAAGCTGATGTTGAGTTTGCTACTCATGAAGAAGCAGTGGCAGCTATGTCCAAGGACAGGGCC
 AACATGCAGCACAGATACATAGAACTCTTCTGAATTCACAACAGGGGCTAGCAATGGGCTTATAGCA
 GCCAGGTGATGCAGGCGATGGGCGTGTGAGTGCAGCTGCCAGGCAACTTACAGTGGCCTGGAGACCGAGT
 GAGTGGCTGTTACGGGGCCGGCTACAGCGGTGAGAACAGCATGGGCGGATATGAT

ACGGTACGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG217122 representing NM_001166432
Red=Cloning site Green=Tags(s)

MMLGPEGGEGYVVKLRGLPWSCSIEDVQNFLSDCTIHDGVAGVHFITYTREGRQSGEAFVELESEDDVKLA
 LKKDRESMGHRYIEVFKSHRTEMWVWKHSGPNSADSANDGFVRLRGLPFGCTKEEIVQFFSGLEIVPNG
 ITLPVDPEGKITGEAFVQFASQELA EKALGKHKERIGHRYIEVFKSSQEEVRSYSDPPLKFMVSVQRP
 DRPGTARRYIGIVKQAGLDRMRSGAYSAGYGGYEEYSGLSGDYGF TTDLFGRDLSYCLSGMYDHR
 YGDSEFTVQSTTGHCVMRGLPYKATENDIYNFFSPLNPVRVHIEIGPDGRVTGEADVEFATHEE
 AVAAMSKDRANMQHRYIELFLNSTTGASNGAYSSQVMQMGVSAQAATYSGLESQSVSGCYGAGYSGQNSMG
 GYD

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001166432

ORF Size:	1245 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001166432.1 , NP_001159904.1
RefSeq Size:	2288 bp
RefSeq ORF:	1248 bp
Locus ID:	98758
UniProt ID:	Q9Z2X1
Cytogenetics:	6 F1
Gene Summary:	Component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes which provide the substrate for the processing events that pre-mRNAs undergo before becoming functional, translatable mRNAs in the cytoplasm. Plays a role in the regulation of alternative splicing events. Binds G-rich sequences in pre-mRNAs and keeps target RNA in an unfolded state (By similarity).[UniProtKB/Swiss-Prot Function]