

Product datasheet for **MG217263**

Arid3c (NM_001017362) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arid3c (NM_001017362) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Arid3c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG217263 representing NM_001017362 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGCCCTGCAGAGGCAGCAGGCAGCACGGCTGGCCCAAGGGGTGGGGCCCTTGCCCTCCACGCT
TACCGCTTCCACAACCTCCTCTGCTTGGCGCCCGGACCTACAGGCTCCTGAGGGGGCTATAGGGGTGGT
TGGGGCTGAAGAAGAGGAGGGTCTGAAGATGAGGAGGGAGAGACACCTTAGCAGAAGAGGAGACAGCT
GAGCAGAGCCATCCAGGAGCCCGATGTCCCAACTCACCTTCCAGCCAGTCTCCTGGAATCCAACCACATG
AGTGGACCTACGAGGAACAGTTCAAGCAGCTGTATGAAGTGCATGCGGACCCCAAGAGGAAGGAGTTTCT
GGATGACCTATTTAGCTTTCATGCAGAAGAGGGGTACGCCAGTGAACCGGGTACCCATCATGGCCAAACAG
GTGCTGGATCTGTATGCGTTGTTTCGCTTGGTGACAGCCAAGGGAGGCCTGGTGGAAATCATCAACCGAA
AAGTGTGGCGAGAAGTCACACGCGGCCTCAGTTTGCCACCACCATCACCTCTGCCGCTTACAGCTCCG
TACACAGTACATGAAGTATTTGTACCCATACGAGTGCAGACACGGGCACTCAGTCCCCAGGGGAGCTT
CAGGCTGCCATAGACAGCAACAGGCGTGAGGGCCGCCAGGCTTACACCGGGTCCCGCTCTTCAACT
TAGCAGGGCCACACCTCGGGCGCTCCTGGCCCCGCTCGAGTCATGGGCCGCCCGACCGGACCC
GAACTGCCGGGTCCCACCCAGGTTCTGCTTCTGGTTTACCAGCGCACGCTTGCCTCACTGAGCCCG
AGCCAGTAAAGAAAGAGGAGAGTGGAATCCCACCCCTCGTCTGCACTACCTATGGGCTTGGCCTCGG
AAGCTACAAGGGAGAAGCTGGCACCAGGAAGACCTCCGGAGAAGAGGGCTGTGCTAATGGCCCCGTGGA
CTCACCTCGACTGGGGGACCCCGAGTTTCTGCCCCGTGGCAAGGCTCCCTAAGGGAAGAGCGGTTG
GATGGGCCCTCAATCTGGCAGGCACTGGCATCAGCAGTATCAACGTGGCACTAGAGATCAATGGGGTGG
TCTACACTGGTATCCTCTTTGCTCGACGCCAACCTGTGCCAGTTCCTTGGGTCCAACCAATCCTCCACC
TCTACCTCTACAGGGCCCCCTTCCAGTACCTTGCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG217263 representing NM_001017362
 Red=Cloning site Green=Tags(s)

MEALQRQQAARLAQGVGGLAPPRLPLPQPPLLGARTLQAPEGAIGVVGAEEEEGADEEGETPLAEETA
 EQSHPGARCPNSPSSQSPGIQPHEWTYEEQFKQLYELDADPKRKEFLDDLFSFMQKRGTPVNRVPIAKQ
 VLDLYALFRLVTAKGGLVEVINRQVWREVTRGLSLPTTITSAFTLRQYMKYLYPECETRALSSPGEL
 QAAIDSNRREGRRQAYTAVPLFNLAGPTPRGAPGPASSHGPAATPNPCPGPTQGSASGLPAHACAQLSP
 SPVKKEESGIPPPRLALPMGLASEATREKLAPEEPEKRAVLMGPVDSRPLGAPPSFLPRGKAPLREERL
 DGPLNLAGSGISSINVALEINGVYVTGILFARRQVVPASLGPNTNPPPLPSTGPPSSTLP


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001017362

ORF Size:	1227 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001017362.2 , NP_001017362.1
RefSeq Size:	1751 bp
RefSeq ORF:	1230 bp
Locus ID:	550619
UniProt ID:	A6PWW5
Cytogenetics:	4 A5
Gene Summary:	This gene is a member of the ARID (AT-rich interaction domain) family of proteins. The ARID domain is a helix-turn-helix motif-based DNA-binding domain. ARID family members have roles in embryonic patterning, cell lineage gene regulation, cell cycle control, transcriptional regulation and possibly in chromatin structure modification. [provided by RefSeq, Jul 2008]