

Product datasheet for **MG217282**

Ly6g6c (NM_023463) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ly6g6c (NM_023463) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Ly6g6c
Synonyms: 1110003M04Rik; AU016360; G6c; NG24
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG217282 representing NM_023463
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAACACCTCCTGTTGCTCACCTGTCTGCCCTACTCTACTGCTGGGTCTCAGCTGATACTCGATGTC
ACTCCTGCTACAAAGTCCCTGTGCTGGGCTGTGTGGATCGCCAGTCTGCCGCCTGGAGCCGGGCCACAA
ATGCCTGACAACAAACGTGTACCTTGGGAAGATGTGGGTTTTCTAACCTGCGCTGCGGCACACCAGAA
GAGCCTTGTCGGGAGGTCTTCAACGAAACCAACCATAAGCTGGGCTGAACTACAACACCACCTGCTGTG
ACAAGGATAACTGTAACAGCCCGCTCCACGGCCACACCCGCACTGGCCCTCATCTCCCTCACCTCCTT
GGCTGGCCTCGGCTCTGGTTATTGCAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG217282 representing NM_023463
Red=Cloning site Green=Tags(s)
MKHLLLLTLSALLYCWVSADTRCHSCYKVPVLGCVDRQSCRLEPGHKCLTTNVYLGKMWVFSNLRCGTPE
EPCREVFNETNHKLGLNYNTTCCDKDNCNSPAPRPTPALALISLTSLAGLGLWLLH

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN:

NM_023463

ORF Size:

378 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_023463.3</u> , <u>NP_075952.1</u>
RefSeq Size:	960 bp
RefSeq ORF:	381 bp
Locus ID:	68468
UniProt ID:	<u>Q9Z1Q4</u>
Cytogenetics:	17 B1