

## Product datasheet for **MG217407**

### **Cobll1 (NM\_177025) Mouse Tagged ORF Clone**

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                         |
| Product Name:             | Cobll1 (NM_177025) Mouse Tagged ORF Clone                                   |
| Tag:                      | TurboGFP                                                                    |
| Symbol:                   | Cobll1                                                                      |
| Synonyms:                 | 1810047P18Rik; Coblr1; D430044D16Rik                                        |
| Mammalian Cell Selection: | Neomycin                                                                    |
| Vector:                   | pCMV6-AC-GFP (PS100010)                                                     |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                                      |
| ORF Nucleotide Sequence:  | >MG217407 representing NM_177025<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGACCCGACGCTCCCGGATCCCCTACCCCGCAGCGCGCCCGCACCCCGCCATGCAGCCCGCGGGCT  
CCGCAGGCAGGAAAACAAAGGGAAGGCGCCTCTTCTCCGGCCGAGACGAAGGGCACTGATGTCTCTTC  
TGCTGAAGACCTGTAGAATCTACTGCTGTCGTACAGAGCAACAAGACAATATGATAGATAAAGATATT  
GAGCTCTCAGTGGTCTGCCGGGAGATATCCTCAAGTCTACTACTGTTTCATGGCAGTAAGCCAATGATGG  
ACTTGCTCGTGTCTCTGTGCTCAGTATCACTTAAATCCGTCAAGTACACAATTGATTTGCTCTCTGC  
TGAAGAGAACCTTATCAAATTTAAGCCAAACACACCAATAGGAATGCTGGATGTAGAGAAGTAATTTTA  
AAGCCAAATCTCTGGATAAGAAAAAGCCTACACCTATAATCCCAGAGAAAATGTGAGAGTAGTGATCA  
ACTTTAAGAAAACACAGAAGACCATAGTGAGAGTTAGTCCCCACGCACCCCTTCAAGATCTTGCTCCCAT  
CATATGCAGTAAATGTGAGTTTGATCCGTTGCATACCGTGCTGTTGAAGGATTACCAGGCTCAGGAGCCC  
CTCGACTTGACAAAATCTCTTAATGACTTGGGACTAAGAGAGTTATATGCCATGGATATCAGCAGAGAGT  
CCTGCCAAATATCACACAACCCGGACATTGTGAAGGAGAAAGAAAACAAAGGGATTTTCAGTTTTTTCCA  
ACGCAGCAAGAAAAAGCGAGAGCAAACTGCCAGTGCTCCTGCGACCCCACTAGTGAGCAAGCATCGTCCA  
TCTTTCACGAGGTCCAACACCATTTTCAAACCGTACATTTCCAATACCCTGCCCTCGGATGCACCCAAGA  
AGAGACGTGCCCCACTGCCACCGATGCCACATCCCAGGGTGCTGCTCAAGGCCAGGAGAGACGAGCGTC  
TTGTGTGGAGAGATCCACCAGTGTGGACGACACAGATAAGAGTTCCTCTGAAGCGATCATGGTGAGGACA  
GGATCACTGCAGCTCAGCAGCACGTCGATAGGAACTTCATCTCTGAAGAGAACGAAGCGAAAAGCACCTG  
CCCCACCTTCGAAAACCCCCCTGGCTCAGACTGATGAAAGGAACTCTGCAATGGCTCATGGCCTACCTTT  
GGAAGATGGGATTGCTCCAGACAGCATGTTGGAAGTGAAGTCTCCCGAGGGAATGTCCACCCAGAGGGC  
TCCCTTGGTCCCGGTTCTAAGCCAGGAGCAGTGTGCTGTGCCAAGCCGCCGATGAAATCTCTGAGG  
GCCCTGGAAGTCCAGAGGCAGCTGTAGCATCTTGACCTCTGGTGAAGCTCTGATTACAGTCTTGAGGA  
GATAGATGAAAAGGAAGAACTCAGCGAGGCATCTAAAGATCCGGCTGGAAGTATTCTGTGAAGTCGCT



[View online »](#)

GACATTGCTTCAGCATCGACTGATATGAGAATTACAGTGGAGAAGGATCCTGATTACGCCCTCGGCATCA  
 GTGATGGAGAGACCTCTCCAAGCTCCAAGGGAAGCAAGCAGAGAGGGCAAGGGCC  
 GTATCATCCTGTTGTAGGCCATATAGGCAATGAGGACAGAGTATCTGACAGTATAAAGACATGAAGACA  
 CTGGGCCCCAAACCAGGAGAGTGTCAAAATGAAATCATGGTCTGTGCCACAAGCACAGATTATGTGAAAA  
 ACAGACCGGGGAAAATGGAAACCACCATAGAAGGCGAAGGCGAAGCTCTGAAGAAAGCCAGCGACATGGA  
 AACTGATAGATTGTCAGGCTCTCCTGCATGTAGAATGGACAATGTTAAAAGTTTCGAGGGAAAAATCATCTA  
 ACAGCATCGCCAGGACCAGATCAAAAATTGAATCAACCAGGTGTAGAGAAGACGAAAAATGCAGGATGCAG  
 CAATTCAGGCCACCCCTGCCCGCGCACTTTTGATGGGAATCACGAAGTCAGTAATTCATCTGACCCAG  
 GGCTGATGAGACTGTGCAGACTTCGGATGGCAGCATATCAGCTCAGCACTCCTCTGCAAGTCTCCAGGAT  
 TCTGTAATGCCTCAAGAGAATTCAGGAGTCAGGGCACCTCAACTTATGTACAAGATAGACTCCCTGAGA  
 AAGAGCCAGCTTGTACATATGGGAACAATGTCCCTTTGTACCTGTAGATGGAAGTAATAAGAATCCAGC  
 TGCTTCTTACCTAAAGAACTTTCCACTTTATAGGCAAGACTCCAACCCCAAGCCAAAACCGTCAAATGAA  
 ATCAGGAGAGATATACCTAAAATTGGAATGACTACTTATAAAATAGTACCTCCGAAATCCCTGGAAA  
 TGGCGAAAGACTGGGAATCAGAAGCAATGGGGAGGAAAGATGACCAAAAGATGCTCCCTGTAGGGCAGAG  
 GCATACGATTGAGAATATGACAGAACTTCCATGCAAACAGAAGTCCCTGCTACTTCTAAAAGCTCACAG  
 CAGCCTCAGCCTGACCTTAAACCAAGCCAGCTCAGGAACAGAGCGTCAATTTGCACAGGACTTTGAGCT  
 CACCAACTGGTACCGAGACGAACCCTCCCAAAGCACCAGAGTGACTACAGACACTGGCACCATTCTTTT  
 TGCGCCAAATTTGGAAGATATAAACAATATTTTGGAGTCAAAATTTAGATCCCGAGCTTCAAATCCCCAG  
 GCCAAACCAAGTTCATTCTTTTGCAGATGCAAAAGAGAGCATCAGGCCACTATGTGACATCTGCAGCTG  
 CTAAGAGCGTCCACACTGCCCTGGGCTGCCCTAAAGAGCCAACAATCAAGGAGGTGCAGAGGGACCC  
 ACAGCTTTCTCCAGAACAGCATCCTTCTCCTTAAGTGAAGGACTCACTCTGCTCCGCTGCCCAATATT  
 TCAAAAGCTGATGATGATATAATTCAGAAGCCTGCCGAGACCTCTCCTCCTCTGTAGCTCCCAACCAA  
 TGACTCTTCGTGCCGAGACCTCTCCTCCTGTATTTCCCAACCAATGACTCTTCTGTGAGACCTC  
 TCCTCCTCTGTATTTCCCAACCAATGACTCTTCTGCCGAGACCTCTTCTCCTTGTATTTCCCAAA  
 CCAATGACTCTTCGAGCCGAGACCTCTCCTCCTCTGTAGCTGCCAAGCCTGTGGCTCTCCTGGTAGTC  
 AGGGAACATCACTAAACCTGAAGACCTTGAAAACATTTGGGGCCCAAGACCATACTCCAGTTCTGGCCC  
 CTCACCTTTGCCCTTGCTGTGGTGAAAAGGTCACAGTCGTTAGTAAGGCATGCCAGAGTCAGCCAGT  
 GAGGGCTCCTCTGCCCTGCCTCCAGCTGCCACACAGGATGAAAAGACACATACAGTAAACAAGCCTACGG  
 TGGGCTCCCAACATGGTGTGGTGATAAGCAAAATAATCCTGTACAGAATGAGCACAGTTCTCAAGTGCT  
 GACTCCAGCTGATGGCCATCGTTCACCCTTAAGAGACAGAGTTCTTAACATTCCAGAGCTCTGACCCA  
 GAGCAGTACGACAGAGCTTGCTCACTGCAATCCGGTCTGGAGAGGCTGCTGCCAAGTTGAAGAGAGTTA  
 CTGTTCCATCGAATACAATATCTGTGAATGGGAAGTCGGGCTCAGCCAGTCCATGTCCATTGACGCCCA  
 GGACAGTCGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG217407 representing NM\_177025  
Red=Cloning site Green=Tags(s)


MDRSVPDPLPR SAPRTPAMQPAGSAGRKTKGKAPLPPAETKGTDVSSAEDPVESTAVVTEQQDNMIDKDI  
ELSVVLPGDILKSTTVHSGKPMMDLLVFLCAQYHLNPSSHTIDLLSAEENLIKFKPNTPIGMLDVEKVIL  
KPKSLDKKKPTPIPEKTVRVVINFKKTQKTIIVRVSPHAPLQDLAPIICSKCEFDPLHTVLLKDYQAQEP  
LDLTKSLNDLGLREL YAMDISRESCQISHNPDI VKEKENKGIFSFQRSKKKREQTASAPATPLVSKHRP  
SFTRSNISKPYISNTLPSDAPKKRRAPLPPMPTSQGAAQGGERRASCVERSTSVD DTKSSSEAIMVRT  
GSLQLSSTSIGTSSLKRTKRKAPAPPSKTPLAQTDERN SAMAHGLPLEDGIAPDSMLELSSPEGMSTPEG  
SLGPGFLSQEQCAVPKPPDEISEGPGTPEAAVASLTSGVSSDYSLEEIDEKEELSEASKDPAGSISVKSP  
DIASASTDMRITVEKDPDSALGISDGETSPSSKGTQEGRSTEGQGYHPVVGHIGNEDRVSDSIKDMKT  
LGPNQESVQNEIMVCATSTDYVKNRPGKMETTIEGEALKKASDMETDRLSGSPACRMDNVKSSRENHL  
TASPGPDQKLNQPGVEKTKMQDAAIQATPARGTFDGNHEVSNSSDPRADETVQTS DGSISAQHSSASLQD  
SVNASREFRSQGTSTYVQDRLPEKEPACTYGNNVPLSPVDGSNKNPAASYLKNFPLYRQDSNPKPKPSNE  
ITREYIPKIGMTTYKIVPPKSLEMAKDWSEAMGRKDDQKMLPVGQRHTIENMTETSMQTEVPATSKSSQ  
QPQPDLPKPKSSGTERHLHRTLSSPTGTETNPPKAPRVTTDTGTIPFAPNLEDINNILESKFRSRASN PQ  
AKPSSFFLQMQRASGHYVTSAAAKSVHTAPGPAPKEPTIKEVQRDPQLSPEQHPSSL SERTHSAPLPNI  
SKADDDIIQKPAETSPPPVAPKPMTLRAETSPPPVFPKPMTLPAETSPPPVFPKPMTLPAETSLPLVFPK  
PMTLRAETSPPPVAAKPVALPGSQGTSNLKTLKTFGAPRYPSSSGSPFALAVVKRSQSF SKACPESAS  
EGSSALPPAATQDEKTHTVNKPTVGSQHGDKQNNPVQNEHSSQVLTPADGPSFTLKRQSSLTFQSSDP  
EHVRQSLLTAIRSGEAAAKLKRVTVPNTISVNGKSGLSQSMSIDAQDSR

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN: NM\_177025

ORF Size: 3720 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_177025.6</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 5042 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 3723 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 319876                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <a href="#">Q3UMF0</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 2 C1.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

### Product images:



Circular map for MG217407