

Product datasheet for **MG217531**

Pmepa1 (NM_022995) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pmepa1 (NM_022995) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pmepa1
Synonyms:	2210418I02Rik; AW455466; Erg1.2; N4wbp4; Stag1; Tmepai
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG217531 representing NM_022995 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGTCAACGGCACCGCCGCCGCCGCCGCCGGCAGCCCAATGTCTCCTGCGCGTGCAACTGCCAGC
GCTCTTTGTTCCCGAGCATGGAGATCACGGAGCTGGAGTTCGTGCAAATCGTGGTCATCGTGGTAGTGAT
GATGGTGATGGTGGTTATGATTACGTGCCTGCTGAGCCACTACAAGCTGTCAGCCCGCTCCTTCATCAGC
CGACACAGCCAGGCCAGGAGGAGAGACGATGGACTGTCTCGGAAGGATGCCTCTGGCCCTCAGAGAGTA
CGGTGTCAGGTGGAATGCCGGAGCCACAGTCTATGCCCGCCTCGGCCACTGACCGACTCGCTGTGCC
CCCTTCATCCAGCGGAGCCGATTCCAACCCACCTACCCCTACCTGCAGCACGAAATTGCCCTGCCACCC
ACCATCTCACTGTCTGATGGGGAGGAGCCCCACCCCTACCAGGGCCCTGCACCCCTCCAGCTACGGGACC
CTGAGCAACAGCTGGAGCTGAACCGGGAATCTGTGCGCGCACCCCTAACCAGGACCATCTTCGACAGTGA
CCTTATAGACAGCACCATGTCTGGGGGCCCTGTCCCCCAGCAGTAACTCGGGCATCAGCGCCACCTGC
TACAGCAGCGGTGGGCGCATGGAGGGGCCGCCCCACCTACAGCGAGGTCATTGGCCACTACCTGGCT
CCTCCTCCAGCACCAGCAAAGTAACGGGCCATCTCCCTGCTAGAGGGGACCCGGCTCCATCACTCGCA
CATTGCCCACTGGAGAACAAGGAGAAGGAGAAACAGAAAGGTCACCCCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG217531 representing NM_022995
 Red=Cloning site Green=Tags(s)

MGVNGTAAAAGQPNVSCACNCQ RSLFPSMEITELEFVQIVVIVVMMVMVMITCLLSHYKLSARSFIS
 RHSQARRRDDGLSSEGCLWPSESTVSGGMPEPQVYAPRPTDRLAVPPFIQRSRFQPTYPYLQHEIALPP
 TISLSDGEEPPPYQG PCTLQLRDPEQQLNLNRESVRAPPNRTIFDSDLIDSTMLGGPCPPSSNSG ISATC
 YSSGGRMEGPPPTYSEVIGHYPGSSFQHQQSNGPSSLLEGTRLHSHIAPLENKEKEKQKGHPL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_022995

ORF Size: 822 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022995.3 , NP_075371.2
RefSeq Size:	4531 bp
RefSeq ORF:	825 bp
Locus ID:	65112
Cytogenetics:	2 H3