

Product datasheet for **MG217626**

Kbtbd11 (NM_029116) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kbtbd11 (NM_029116) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kbtbd11
Synonyms:	2900016B01Rik; 4930465M17Rik; mKIAA0711
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG217626 representing NM_029116
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGAACTCGGTGGCGCTTTCGTCCTACTCGGGACTGAGCCTAGGACTCCTGGAGAGGACAGCC
 TCCTCTGCCAGCCGAGGAAGAGGGCGCTGCGTCCACAGCGCAGACACCTGCAGTCTCAGCGCCTCATT
 GTGCTTACAGTCCGGGGACGACTCCCGCCGAGTCTCGCGCTCTGCCGCGAGGGCTCCGAAGCCTCT
 CCGCCTTCTCTCCGACGACCTGCGTGTGGTGGAGACGAGTGGGATGTCAGCAGCGCCGCTCCCGAG
 AGTCCCGGAGGAGTGCAGCGCCCGGAAGAGCCCGCTCGCGGAGGATCCCCCTCCGACACGAGCA
 CGCGCTCCCGTGAACCTCGAGTCTTGGACGAGCTTGGGAGCCCGTGCCTGTGCCGCCGGGGTGGG
 TCGGTGCACGGAGAACCAGCTGGTATCGAGGTGGCCGGTCCCGACTACGGGCGCATAAGCGGTGC
 TGGCCGACGCGAGGACTACTCCGCGCGCGCATCCCGGACGTGCTGCGGGTCAAGCGGTGAGCTT
 CACTGCCCTGCGCTACTGCTGGCCGACGATACAGCGGGCGCATGGCGGGCTGCGGCCGACAACGTG
 GCCGAGGTGGTGGCCGGCGCGCCGCTCCAGCTGCCCGGCGTGCAGCAGCGCCACCGAGGCCATGG
 CGCCGACGCTGAGCCTGGACAACCTGCTACGAGGTGCTGAGCGCAGGCAAGCGGCAGCGGCTGACCGAGCT
 GCGCGAGCCGCTACCGCTTATGAGCGACACTACCTGGAGGTGCTGCGCGAGCCCGCCTATTCCGGC
 CGCCTGTCCGGTGCAGAACCGCATCTGCTGCTGCGCGCCGCTGTGTACCGGCCGCGCTTGTCTGCTGG
 CCGCTGCGCTCGGGACAACCGGGAGCGCTTGGCAGTGCACCGCAGAGCCCGTCCGGGGACGAGATC
 GCGCGGGACGCGGCCGTTTATTGCTACCAGCGGAGGCCGGGGAGTGGCGGAGCTGACCCGGCTGCC
 GAGGGCGGCCAGCGCGGCTGCGGCTTCTGCTGCTTCACTACCTTCTCCGCGGGCGGTGGTGG
 CTCTGCGGGTCCGATGGCCGCGCGCCGCTCGGACAGGTGATTGCTACAACCCGTCACCGACTC
 CTGGAGTACTGTGCGCCGCTGCGCCAGGCACGCTCGCAGGTGCAGCTGCTGGCCCTGGACGGCCACTG
 TATGCCGTGGGTGGCAGTGCCTGCTCAGCGTGGAGCGCTACGACCCGCGAGCTGACCGATGGACGGCGG
 TGGCCCGCTTCCCGTGGCGCTTTCGAGTGGCCACGAGGCTGCCACCTGTAACGGAGAGATTTATGT
 GTCCGGGGGCTACTCTTCTACCGCTGCTCAAGTACGACCCCGCGCGAGAGTGGCAGGAGTGCCTCA
 TGCAGCAGCAGCCGAGAGCGCTTGCAGCATGGTGGCCCTCGACGGCTTCTCTATCGCTTTGACCTGT
 GTGGGAGCCGTGGGAAGCTCAGGCTGCAGTGGTTCAGGTGGCGGGTACGCGTTTTCCGCTACCATTG
 CCTGGCCAAGCAATGGAGCAATGTGCTGTACACCTGAGGCCCTGGCGACCCGAGGCTCCAGCCC
 TTTGTTGTGTTGCCCTGGATGGTACCATCTACTGTGTGAGCCGCGCGGACCTGGCGTTTCGTGCTT
 CTAGGACACCGAGGCTGGCAGGATATGGACCTGGTGGCAGCTTTGAGCCGAGCCCTGGGATCCCC
 CCTGGATGTCCGGGGCTACTCTTCCATTTGACTCAACCTGCCTGAGAAGCCGACCGCGGGGAGCAG
 GGTGCAGTC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG217626 representing NM_029116
 Red=Cloning site Green=Tags(s)


MENSVAPFVLYSGTEPRTPGEDSLPLPAEEEGAASTAQTPCSLASLFCSSGDDSPQSRASAAEGSEAS
 PPSLRSDLRVETQWDVSSAASPESPEECARPEEPASPEDPPSRHEHARPVELESDELGEPVPPPGVG
 SVHGEPDLVIEVAGRRLRAHKAVLAARSDYFRARASRDVLRVQVSFTALRLLADAYSGRMAGVRPNDV
 AEVVAGARRLQLPGAQRATEAMAPQLSLDNCYEVLSAGKRQRL TELRDAAYRFMSDHYLEVLRPAVFG
 RLSGAERDLLRRRLCTGRACLLAAALGTTGERSGSRPQSPSGDAESRGDAAVYCYQAEAGEWRELTRLP
 EGAPARGCGLCVLNFNYLFLAGGVAPAGPDGRARPSDQVYCYNPVTDVSWSTVRPLRQARSQVQLLALDGH
 YAVGGECLLSVERYDPRADRWTAVAPLPRGAFVAHEAATCNGEIYVSGSLFYRLLKYDPRRDEWQEQCP
 CSSSRERSADMVALDGFLYRFDLCGSRGEAQAAVGSGGVSVFRYHCLAKQWSQCAVHLRPPGAPAGLQP
 FRCVALDGTIYCVSRAGTWRFVPSQDTEAGSDMGPGGSEFEPEPLGSPLDVRGVLPFVNLNLEPKPDRGEQ
 GAV

TRTRPLE – GFP Tag – V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_029116

ORF Size: 1899 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029116.2](#), [NP_083392.1](#)

RefSeq Size: 6996 bp

RefSeq ORF: 1902 bp

Locus ID: 74901

UniProt ID: [Q8BNW9](#)

Cytogenetics: 8 A1.1