

Product datasheet for **MG217742**

Bean1 (NM_001141925) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Bean1
Synonyms:	Bean
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG217742 representing NM_001141925
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGTCCTTCAAACGCCCTTGCCCTTAGCACGATACAACCGTACCAGCTATTTCTACCCTACCACGTTCT
CCGAGAGCTCTGAGCACAGCCACCTGCTGGTGTCTCCGGTGTGGTGGCGAGCGCTGTCTAGAGTGTGGT
CATCACTCTCTCTGCATCACCATCATTGTGGGAAGCATCAGAAGGGACAGACAGGCCCGGATCCAGCGA
CACCACCATCGTCACCGCCGCCATCATCACCATCACCGCCACCGCAGGCGCCGACACCGAGAATACGTCG
CAGGGGTGACATCCATGAGAGCAGGACCCGTGTGGTCCACCCACCTGCAGCCAGAGGAGGTTTTTC
AGAGGCTATGGAAGGACTGACCTTCACGAGTCTTTGTGCATCTGGAGGACACACACAGCCGTTCTAGC
CCCAGGATGCCCTATGCCTGCAGCCCTGCTGAGGACTGGCCACCGCCCTTGACGCTCAGCTCCGAGGGGG
ATGTGGATGTCACAGTGTCTGGGAGCTGTACCCAGATTCTCCACCAGGCTATGAGGAATGCATGGGACC
AGGGGCCACTCAGCTGTATGTCCCAACGATGCACACCGCCCTACTCCATGACTGACTCCTGTCCCAGG
CTGAATGGCGCCCTTGACTCAGACAGCGGCCAGAGCCGACGCCACCGTCAGCAAGAGCAGAGGACCCAGG
GTCAGAGCAGACTCCACACTGTCTCCATGGACACACTGCCACCTTATGAAGCTGTGTGTGGGACAGGCTC
TCCATCTGACCTGCTGCCATTGCCAGGACCAGAACCATGGCCAAGTAACCCCAGGGCTACCCATCCCA
ACCCAGGCTCCGATGCCAGCCCTGAAAGGATTGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG217742 representing NM_001141925
 Red=Cloning site Green=Tags(s)

MSFKRPCPLARYNRTSYFYPTTFSESSEHSHLLVSPVLVASAVIGVVITLSCITIIIVGSIRRDQRARIQR
 HHHRRHHHHHHRRRRHREYVAGVTSMRAGPVCGPHLQPRGGFSEAMEGLTFTSLCASGGHHSRSS
 PRMPYACSPAEDWPPPLDVSSGDDVDVTLWELYPDSPPGYEECMGPATQLYVPTDAPPPYSMTDSCPR
 LNGALDSDSGQSRSHRQQEQRTQGQSRSLHTVSMIDLPPYEAVCGTGSPSDLLPLGPPEPWPSNSQGSPIP
 TQAMPSPERIV

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001141925

ORF Size: 876 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001141925.1](#), [NP_001135397.1](#)


RefSeq Size: 3191 bp

RefSeq ORF: 585 bp

Locus ID: 65115

Cytogenetics: 8

Product images:



Circular map for MG217742