

Product datasheet for **MG217752**

Ccdc50 (NM_001025615) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc50 (NM_001025615) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ccdc50
Synonyms:	C3orf6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG217752 representing NM_001025615 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGACGTGAGTGTAGATCAGTCCAAGTTGCCGGCGTGAAGGAAGTGTGTAGAGATTTTGCCGTCC
TGGAGGACCACACCCTGGCCATAGCCTGCAGGAACAAGAGATTGAGCATCATTTGGCATCCAACATTCA
GCGGAACCGTCTGGTACAACATGATCTGCAGGTTGCTAAGCAGCTCCAAGAGGAAGACCTCAAAGCCAA
GCTCAGCTCCAGAAGCGCTACAAAGCCCTGAACAACATGATTGTGAAATTGCTCAGGAAATCCAGGAGA
AGCTAACCATTTGAGGCTGAGAGACGACGCATTCAGGAAAAGAAGGATGAGGACATAGCACGTCTTTTGA
AGAGAAGGAGCTACAGGAGGAGAAAAGGAGGAAGAAACACTCCAGAGTTTTCTGGGGCAGTGTTTT
GGAGATAACTACTATCATGAGGATGGAGGAATGAAACCAAGAGGAATAAAAGAAGCTGTCTACTCCAG
CACGAGCAAGCCACAGGGACCAGGAGTGGTATGATGCTGAAATTGCCAGGAAATTGCAAGAAGAAGAACT
TTTGAAATTGCTCGACTTCTAATGGCTGAAGAAAAAAGCTTACAAGAAAGCCAAAGAGCGAGAAAAAG
TCATCTTTGGACAAAAGGAAACATGACCCTGAATGCAAGTTAAAAGCAAAGTCAGCCCACTCAAAGTCAA
AAGAGGGTGATGAAGCACACCGCTCCAAGATTGACAGGCCATCAAGACCACCACCACCTACCATGATGGG
CCTTGAGGACACAGATCCCACCCATTTTACAAACCAGCACAGTACAACATGGCATCTTCCAAAGTCAGAG
TCCTCACAGAAAGGCTTCCATAACAAGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG217752 representing NM_001025615
 Red=Cloning site Green=Tags(s)

MADVSDVQSKLPGVKEVCRDFAVLEDHTLAHSLQEQEIEHHLASNIQRNRLVQHDLQVAKLQEEDLKAQ
 AQLQKRYKALEQHDCEIAQEIQEKL TIEAERRRIQEKKDEDIARLLQEKLQEEKRRKKHTPEFSGGSVF
 GDNYIHEDGGMKPRGIKEAVSTPARASHRDQEWYDAEIARKLQEEELLE IARLLMAEEKKAYKKAKEREK
 SSLDKRKHDPCKLKAHSAHKSKEGDEAHRSKIDRPSRPPPTMMGLETDPTHTFNQHSSTWHLPKSE
 SSQKGFHNKQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001025615

ORF Size: 870 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025615.3](#), [NP_001020786.1](#)

RefSeq Size: 6936 bp

RefSeq ORF: 873 bp


Locus ID: 67501

UniProt ID: [Q810U5](#)

Cytogenetics: 16 18.98 cM

Gene Summary: Involved in EGFR signaling.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG217752