

## Product datasheet for **MG217848**

### Shf (NM\_001013829) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Shf (NM\_001013829) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Shf  
**Synonyms:** RP23-173H17.6  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG217848 representing NM\_001013829  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGCCCTACGAGGCCAAAAGATGATGGCTGAGATCCGGGTTCTAAGGAGACTGCAGCTCAACCC  
TACCTCTCTATGACACACCCTATGAGCCAGAAGATGAGGGGCCAGCCAGAAGGTGAGGGGACTCCCTG  
GCCAGGGAGTCGCGCCTGCCAGAGGATGACGAGAGGCCCTCAGAAGAGTATGACCAGCCCTGGGAGTGG  
AAGAAGGAACGGATTTCAAAGCCTTTGCTGCCAGTTTGAAGGATCGGAGAAGAAGTGCCTGTCACCTG  
GCCGGGAGGAGAAGGGCGGCTACCTCCAGGCTATCTGCTGGGAACCCCAAGTCAGCCAAGCCCTAGG  
CATGGAACCCAGCAGCCCTGGGGAGTGGACAGATCCAGCACTGCCTCTGGAAAACAGGTCTGGTAC  
CACGGGGCCATCAGCCGAACAGATGCCGAGAACCTGCTCAGGCTGTGCAAAGAGGCCAGCTACCTGGTGC  
GCAACAGTGAGACCAGCAAAAATGACTTCTCTCTCCCTCAAGAGCAGTCAGGGCTTCATGCACATGAA  
GCTGTCCCGGACTAAGGAACACAAGTATGTGTTGGGCCAGAACAGCCACCCTTCAGCAGCGTCCCTGAA  
ATCGTGCACTATGCCAGCCGCAAGCTGCCATTAAGGGAGCGGAACACATGTCCCTGCTCTACCCTG  
TGGCCATCCGACTCTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG217848 representing NM\_001013829  
Red=Cloning site Green=Tags(s)

MEPYEAQKMMAEIRGSKETAAQPLPLYDTPYEPEDEGASPEGEGTWPWPRESRLPEDDERPPEEYDQPWEW  
 KKERISKAFAAQFEGSEKNCLSPGREEKGRLLPRLSAGNPKSAKPLGMEPSSPLGEWTDPALPLENQVWY  
 HGAI SRTDAENLLRLCKEASYLVRNSETSKNDFSLSLKSSQGFMMHMKLSRTKEHKYVLGQNSPPFSSVPE  
 IVHHYASRKLPIKGAEHMSLLYPVAIRTL

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001013829

**ORF Size:** 714 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001013829.1</a>
<b>RefSeq Size:</b>	1564 bp
<b>RefSeq ORF:</b>	717 bp
<b>Locus ID:</b>	435684
<b>UniProt ID:</b>	<a href="#">Q8CG80</a>
<b>Cytogenetics:</b>	2 E5
<b>Gene Summary:</b>	Adapter protein which may play a role in the regulation of apoptosis in response to PDGF. [UniProtKB/Swiss-Prot Function]