

## Product datasheet for **MG218000**

### **Mdm1 (NM\_148922) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Mdm1 (NM_148922) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mdm1
Synonyms:	Arrd2; Mdm-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG218000 representing NM\_148922  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGCCGGTGCCTTCAAGGGCTGAGTGAATACCAGAGAACTTCTGTGAAAAAGTCCTATTTGTCAG  
 AGTCTTATAATCCCTCAGTGGGACAAAAGTACTCATGGGCAGGACTTAGATCGGATCAGTTGGGGATCAC  
 GAAAGAACCAGGTTTTATTTCAAAAAGAAGAGTTCCTACCATGACCCTCAGATTTCAAAAATACCTGGAG  
 TGGAACGGAACCGTCAGAAAAGAAGGATACGCTTGTCCCACCAGAACCCAGGCCTTTGGAACGCCAAAGC  
 CACAAGAGGCTGAGCAAGGAGAAGATGCCAATCAAGAAGCAGTTCTCTACTAGAGGCTCCAGGTTCC  
 CAAGAGAACTCGTCTCATTCTGCGGACTCGAGAGCTGAAGGGTTTCAGACACTGTGAAAAAGCACCGAG  
 GGTGTACGAGAAGCCATGCGCCAGTTAGCGCGGATGTGGAGCTGAGACCTCCAGCAAACAACCTCTCT  
 CCCAGAGCATAGATCCCAGTTGGATAGGCATCTTCGTAAGAAAGCTGGATTGGCCGTTGTTCCCACGAA  
 TAATGCCTTGAGAAATCTGAATACCAAAGGCAGTTTGTGGAAAGACTTCTAAAGAAAGCGCTCCAGTG  
 TTTGCATCCAATCAGTTTTCCGTAATAAAAGCCAAATTAATCCACAGTTCCAAGGCAATACATTCACCC  
 ACGAGACTGAATACAAGCGAAATTTCAAGGGTTAACTCCAGTGAAGGAACCAAAGTCAAGAGAGTATTT  
 GAAAGGAAACAGCAGTCTGGAGATGCTGACTCCAGTAAAGAAGGCAGATGAGCCTTTAGACTTAGAAGTA  
 GACATGGCGTCGGAAGACTCAGACCAGTCTGTAAGAAGCCTGCTTCATGGAGACACCAAAGGCTTGAA  
 AAGTGAATTTGAAATATAGAGCAAAGTTCTGAGCCCAGCCAGTATTTCTATAAAGCTGGAGCTTGAC  
 CCGGGTGAAGGAGAACCTGTCAAACCAGTTAAGGAGCTCCGAGAAAAGGCCGAATCTTACAGGAAGCGA  
 GTTCAGGGGACACATTTTCTCGGGACCATCTGAACCAGATTATGTCGGACAGCAACTGCTGTTGGGACG  
 TCTCCTCAGTACAAGCTCGGAAGGCACCGTCACTAGCAACATCCGAGCACTGGATCTTGCTGGAGACCT  
 TACAAACCACAGGACCCCCAGAAAACACCTCCTACCAAAGTGAAGAAAAGAAAAGTTGCCTCGGGAGAG  
 CAGCCCCTGAAAAACTCCACCAGGAGACTGGAGATGCCAGAGCCTGCCGCCTCGGTCCAGGAGGAAGCTGG  
 CTTGGGATGCTGAGGAGAGCACGAAGGAAGACACCAGGAGGAGCCAGGGCGGAGGAGACGGGAGAGA  
 GGAGAGAGGACAGGACAAGCAGACCTGTGCGGTAGAGCTGGAGAAAACCGGACACACAGACACCCAAGGCA  
 GACAGACTGACAGAAGGGTCGGAGACATCTTCTGTTTCTCAGGGAAGGGAGGCAGGCTTCTACACCGA  
 GGCTGAGAGAACTCGGTATCCAGCGGACGCACCATGATCTCACGACGCCAGCTGTTGGTGGCGCAGTCTT  
 AGTGTCTCCATCTAAAGTGAAGCCACCAGGCTCGAGCAGAGGAGGAGAGCGTCTCCCAAGATGGCTTA  
 GAAACTCTGAAGAAAGACATTAAGAAAGGAAAACCCCGTCCCATGTCTCTGTTGACTTCTCCGGCTG  
 CTGGCATGAAGACAGTTGATCCCCTGCCTCTGCGAGAAGACTGTGAAGCCAATGTGCTCAGATTTGCTGA  
 TACTCTTCTGTTTCGAAAATTTTGACCCGTCAGCCCAGCACCCCTGGGAGCTGCCTCCATGTGCCCCG  
 CCTTACTGTATCCGTCACGAGGATCCAGGGCGCTGCGAGACCTGAGTTTTCAGCACAAATGGGAA  
 AACCCAGGACGAACAATTTGCAGCTACATCCACACGATGCCTTTAATGATGAAGATGCAGATAGACTGTC  
 TGAGATCTCTGCTCGCTCTGCAGTTTCCAGCCTCCGGGCTTTCAGACTCTAGCCCGAGCTCAGAAAAGA  
 AAGGAGAAATTTCTGGGCAAGCCA

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:** >MG218000 representing NM\_148922  
 Red=Cloning site Green=Tags(s)

```


MPVRFKGLSEYQRNFLWKKSYLSESYNPSVQKYSWAGLRSDQLGITKEPGFISKRRVPHYHDPQISKYLE
WNGTVRKKDTLVPPEPQAFGTPKPQAEQGEDANQEAVLSLEASRVPKRTRSHSADSRAEGVSDTVEKHQ
GVTRSHAPVSADVELRPSSKQPLSQSIDPRLDRHLRKKAGLAVVPTNNALRNSEYQRQFVWKTSKESAPV
FASNQVFRNKSQIIPQFQGNFTHETHEYKRNFKGLTPVKEPKSREYLKGNSSLEMLTPVKKADEPLDLEV
DMASESDSQSVKKPASWRHQRLGKVNSEYRAKFLSPAQYFYKAGAWTRVKENLSNQVKELREKAESYRKR
VQGTHTFSRDHLNQIMSDSNCCWDVSSVTSSEGTVSSNIRALDLAGDLTNHRTQKHPPTKLEERKVASGE
QPLKNSTRRLEMPEPAASVRRKLAWDAEESTKEDTQEPRAEEDGREERGQDKQTCAVELEKPDTPKA
DRLTEGSETSSVSSGKGGRLPTPRLRELGIQRTHHDLTPAVGGAVLVSPSKVPPGLEQRRRSSQDGL
ETLKKDITKKGKPRPMSLLTSPAAGMKTVDPLPLREDCEANVLRFADTLPVSKILDRQPSTPGQLPPCAP
PYCHPSSRIQGRLRDPEFQHNMGKPRNTNQLLPHPHDAFNDEDADRLSEISARSAVSSLRAFQTLARAQKR
KENFWGKP
  
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_148922

**ORF Size:** 2124 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_148922.3](#), [NP\\_683724.2](#)

**RefSeq Size:** 3131 bp

**RefSeq ORF:** 2127 bp


**Locus ID:** 17245

**UniProt ID:** [Q9D067](#)

**Cytogenetics:** 10 66.65 cM

**Gene Summary:** Microtubule-binding protein that negatively regulates centriole duplication. Binds to and stabilizes microtubules.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG218000