

Product datasheet for **MG218112**

Pum1 (NM_001159606) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pum1 (NM_001159606) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pum1
Synonyms:	AA517475; mKIAA0099; Pumm
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG218112 representing NM_001159606
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGCGTTGCATGTGTCTTGAAGAGAAAAGCAGTGCTTTGGCAGGACTCTTTCAGCCCCACCTGAAAC
 ATCACCTCAAGAACCAGCTAATCCAACATGCCTGTTGTTTTGACATCTGGAACAGGGTCGCAAGCGCA
 GCCACAGCCAGCTGCAAATCAGGCTCTTGCAGCTGGAACACACTCCAGCCCTGTCCAGGATCCATAGGA
 GTTGCAGGCCGTTCCAGGACGACGCTATGGTGGACTACTTCTTTCAGAGGCAGCATGGTGAAGCAGCTTG
 GGGGAGGAGGAAGTGGTGGAGGCGGCTATAATACTAGCAAACATCGATGGCCTACAGGGGATAACATTCA
 TGCAGAACACCAGGTACGCTCTATGGATGAATTGAATCATGATTTTCAAGCACTTGCTCTGGAAGGGAGA
 GCTATGGCGAGCTATTCCAAAGACCAATGCACTTGTGTTGAGCAGTTGACAGCTGCCAGCAGCAGC
 AGTATGCACTTGCAGCAGCCACCAGCCTCACATCGGTTAGCTCCCGCTGCGTTTGTCCCAACCCATA
 CATCATCAGCGCTGCTCCCCAGGGACGGACCCCTACACAGCTGGATTGGCTGCAGCAGCGACACTCGGC
 CCAGCTGTGGTCCCTCACCAATATTATGGAGTACGCCATGGGGTGTCTACCCTGCCAGTCTCTCCAGC
 AGCAGGACAGTGCAGCTGCAGCAGCTACAACTCTGCTACTCAGCAGAGCGCTCCGACGGCCAGCAGGG
 ACAGCAGCAGGTTCTTCGTGGAGGAGCTAGCCAACGTCCTTTGACCCCAATCAGAACAGCAGGGACAA
 CAAACAGATCCCCTTGTAGCAGCTGCAGCAGTGAATTCTGCTCTTGCCCTTTGGACAAGGATTGGCTGCAG
 GTATGCCAGGTTATCCAGTCTTGGCACCTGCTGCTTACTATGACCAAACCTGGTGCCTTGTGGTGAATGC
 AGGGGCAAGAAATGGCCTCGGAGCTCCAGTCCGGCTTGTGCTCCAGCCCAGTCATCATTAGCTCCTCA
 GCTGCACAAGCAGCTGTTGCAGCGGCTGCAGCTCAGCAAACGGAGCAGCTGGTGGTCTTGCCGGAACAA
 CAAATGGACCATTTGACCTTTAGGCAACCCAGCAGCCCAAGCCAGCCTCAGCAGCAGCACTAACAA
 TCTGGCTTCCAGCTCTTCTATGGCAACAACCTCTGAGCAGTAACTCCAGAGCAGCTCCCTCTTCTCT
 CAGGGCTCTGCACAGCCTGCCAACACGTCCTTGGGATTGGAAGCAGCAGTTCCTTGGTCCACCTGG
 GATCGGCCCTTGGAGGTTTGGAAACAGCAGTTGCAAACCTCCAACACTGGCAGTGGCTCCCGCGTGACTC
 CCTGACTGGCAGCAGTGACCTTTATAAGAGGACATCGAGCAGCTTGGCCCCATTGGACACAGTTTTTAT
 AGCAGCCTTAGCTATTCCTCCTCCTGGACCCGTGGGCATGCCTCTCCCTAGTCAGGGACCAGGACATT
 CACAGACACCACCTCCTCCCTCTTTCACATGGATCCTCTTCAAGCTTAAACCTGGGAGGACTCACAAA
 CGGCAGTGAAGATACATCTCTGCTGCTCCCGTGCCGAAGCCAAGTACCGCAGTGAAGCAGTGCCTCC
 AGCCTTTTCAGCCCCAGCAGCACCTTTTTTCTCCTCTCGTTTGGGATGGAATGTCTGATGTTATGC
 CCTCTGGCAGGAGCAGACTTCTGGAAGATTTTCGAAACAACCGGTACCCTAATTTACAACCTGCGGGAGAT
 TGCTGGGCATATAATGGAATTTTCTCAAGACCAGCATGGCTCTAGATTCAATTCAGCTGAAACTAGAGCGT
 GCCACAGCAGCTGAGCGCCAGCTTGTCTTCAATGAAATCCTCCAGGCTGCGTACCAACTTATGGTGGATG
 TGTTTGGAAATACGTCAATTCAGAAGTTTTTGAATTTGGCAGTCATGAACAGAAGCTGGCTCTGGCTGA
 ACGGATTCGAGGCCATGTCTATCATTGGCACTGCAGATGTACGGCTGCCGAGTCATCCAGAAAGCTTTG
 GAGTTTATCCCTCAGACCAGCAGAACGAGATGGTGGGGAGCTAGATGGCCATGTCTGAAGTGTGTGA
 AAGACCAGAACGGCAACCATGTGGTTCAGAAGTGCATTGAGTGTGTACAGCCCCAGTCTTTCAGTTTTAT
 CATCGATGCATTCAGGGGCAGGTGTTTGTCTACGCACCCGTACGGCTGTAGGGTATCCAGAGG
 ATCCTAGAGCACTGTCTCCCTGACCAGACCCTCCCATCCTGGAGGAGCTTACCAGCACACAGAGCAGC
 TTGTGCAGGATCAGTATGGAATTTATGTAATCCAACATGTACTGGAGCATGGTCGCTCCTGAGGATAAAAAG
 CAAAATTTAGCAGAGATCCGAGGCAATGTACTTGTATTGAGTACAGCACAAAGTTTGAAGCAATGTTGTG
 GAGAAGTGTGTTACTCACGCCTCACGTACAGAGCGTGCTGTGCTCATTGATGAGGTGTGCACCATGAACG
 ACGGTCCCACAGTGCCTTATACACCATGATGAAGGATCAGTATGCTAACTATGTGGTCCAGAAGATGAT
 CGATGTGGCTGAGCCGGTACGCGGAAGATTGTCATGCACAAGATCCGGCCACACATTGCGACGCTTCGA
 AAGTACACCTATGGCAAGCACATCCTGGCAAGCTCGAGAAGTACTACATGAAGAATGGCGTGGACTTGG
 GGCCATTGTGGTCCCCTAATGGTATCATC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG218112 representing NM_001159606
Red=Cloning site Green=Tags(s)

```

MSVACVLKRKAVLWQDSFSPHLKHPQEPANPNMPVVL TSGTGSQAQPQPAANQALAAGTHSSPVPGSIG
VAGRSQDDAMVDYFFQRQHGEQLGGGSGGGGYNTSKHRWPTGDNIIAEHQVRSMDLNDHFQALALEGR
AMGELFQRPNALAVQQLTAAQQQYALAAHQPHIGLAPAAFVNPYIIISAAPPGTDPTYAGLAAAATLG
PAVVPHQYYGVTPWGVYPASLFQQQAAAAAATNSATQQSAPQAQQGQQQVLRGGASQRPLTPNQNQGGQ
QTDPLVAAAAVNSALAFGQGLAAGMPGYPLAPAAYYDQTGALVVNAGARNGLGAPVRLVAPAPVISSS
AAQAAVAAAAASANGAAGGLAGTTNGPFRPLGTQQPQPQQPSNNLASSSFYGNSSLSSNSQSSSLFS
QGSAQPANTSLGFGSSSSLGATLGSALGGFGTAVANSNTGSGSRRDSL TGSSDLYKRTSSSLAPIGHSFY
SSLSSSSPGVGMPLPSQGPQHSQTPPPSLSSHGSSSSLNLGGLTNGSGRYISAAPGAEAKYRSASSAS
SLFSPSSTLFSRRLRYGMSDVMPSGRSRLLEDFRNNRYPNLQLREIAGHIMEFSQDQHGSRFIQLKLER
ATAAERQLVFNEILQAAYQLMVDVFGNYVIQKFFEFGSHEQLALAERIRGHVLSLALQMYGCRVIQKAL
EFIPSDQQNEMVRELDGHVLCVKDQNGNHVVQKCIQVQPSLQFIIDAFKQGVFALSTHPYGCRIQR
ILEHCLPDQTLPILEELHQHTEQLVQDQYGNVVIQHVLEHGRPEDKSKI VAEIRGNVLLSQHKFASN VV
EKCVTHASRTERAVLIDEVCTMNDGPHSALYTMKDKQYANYVVQK MIDVAEPGQRKIVMHKIRPHIATLR
KYTYGKHILAKLEKYYMKNGVLDLGPICGPPNGII
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001159606

ORF Size: 2832 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001159606.2
RefSeq Size:	4635 bp
RefSeq ORF:	2835 bp
Locus ID:	80912
Cytogenetics:	4 D2.2

Gene Summary:

Sequence-specific RNA-binding protein that acts as a post-transcriptional repressor by binding the 3' UTR of mRNA targets. Binds to an RNA consensus sequence, the Pumilio Response Element (PRE), 5'-UGUANAUA-3', that is related to the Nanos Response Element (NRE). Mediates post-transcriptional repression of transcripts via different mechanisms: acts via direct recruitment of the CCR4-POP2-NOT deadenylase leading to translational inhibition and mRNA degradation. Also mediates deadenylation-independent repression by promoting accessibility of miRNAs. Following growth factor stimulation, phosphorylated and binds to the 3' UTR of CDKN1B/p27 mRNA, inducing a local conformational change that exposes miRNA-binding sites, promoting association of miR-221 and miR-222, efficient suppression of CDKN1B/p27 expression, and rapid entry to the cell cycle (By similarity). Acts as a post-transcriptional repressor of E2F3 mRNAs by binding to its 3' UTR and facilitating miRNA regulation (By similarity). Represses a program of genes necessary to maintain genomic stability such as key mitotic, DNA repair and DNA replication factors. Its ability to repress those target mRNAs is regulated by the lncRNA NORAD (non-coding RNA activated by DNA damage) which, due to its high abundance and multitude of PUMILIO binding sites, is able to sequester a significant fraction of PUM1 and PUM2 in the cytoplasm (By similarity). Involved in neuronal functions by regulating ATXN1 mRNA levels: acts by binding to the 3' UTR of ATXN1 transcripts, leading to their down-regulation independently of the miRNA machinery (PubMed:25768905). In testis, acts as a post-transcriptional regulator of spermatogenesis by binding to the 3' UTR of mRNAs coding for regulators of p53/TP53 (PubMed:22342750). Involved in embryonic stem cell renewal by facilitating the exit from the ground state: acts by targeting mRNAs coding for naive pluripotency transcription factors and accelerates their down-regulation at the onset of differentiation (PubMed:24412312). Binds specifically to miRNA MIR199A precursor, with PUM2, regulates miRNA MIR199A expression at a postranscriptional level (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG218112