

## Product datasheet for **MG218201**

### Ccdc90b (NM\_025515) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccdc90b (NM\_025515) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ccdc90b  
**Synonyms:** 2310015N07Rik; C86558  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG218201 representing NM\_025515  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGAGGAGCCGCTGGATTTGGCGGTTCTGCGCCCTGACGGCGGCGGAATCCGCTGGACTTCAACTCCTC  
 ATGGGCGTCTGTACCCGCCCTGCGGAGAGGGTTCTTGACCACCACAACAAAGTCAGATTATGATAGAAG  
 ACCAGTGGATATAACTCCTCTAGAACAAAGAAAATTAACTTTTGATACCCATGCGTTGGTGCAGGACTTG  
 GAAACTCACGGATTTGACAAAACACAAGCACAGACAATTGTATCAGTGCTAAGTACGCTATCAAACGTCA  
 GCCTGGATACTATCTACAAGGAGATGGTGACTAAAGCTCAGCAGGAAATAACAGTCCAACAGCTAATGGC  
 TCAATTTGGATTCCATCAGAAAAGACATGGTCATCCTAGAGAAAAGTGAATTTGCAAATCTGAGAGCAGAG  
 AATGAGAAAATGAAAATTGAACTAGATCAAGTTAAGCAGCAGCTGACTAATGAAACAGTCGAATCAGAG  
 CGGACAATAAGCTGGACATTAACCTGGAAGGAGCAGAGTGACAGACATGTTACAGATCAAGAAAAGCA  
 ACTTATAGAAGCAACCAATGAATTTGCCAAAAGGATACAAAACAAAAGTATTATTTTCAGAGACCAGT  
 AATAAAATTGACTGAAATCGCTTCTTAAAAACTGATGGAGTCGAGCAAACCTTGAACAATTTCGTT  
 ATCTTGCAGCCTCAGTGTCACTTGCTTGCGATAGCATTGGGGTTTTATAGTTCTGGAAGGAGAAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_025515.3](#), [NP\\_079791.2](#)

**RefSeq Size:** 1615 bp


**RefSeq ORF:** 771 bp

**Locus ID:** 66365

**UniProt ID:** [Q8C3X2](#)

**Cytogenetics:** 7 E1

**Product images:**



Circular map for MG218201