

## Product datasheet for **MG218278**

### Lrrn4cl (NM\_001013019) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Lrrn4cl (NM\_001013019) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Lrrn4cl  
**Synonyms:** 1110067112Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG218278 representing NM\_001013019  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTTGGGCTCTTTCCCTTCTGTGGCTGGCAGCCATGACCACCTCCTTGGTTCCCAACCTCAGATCT  
TGACCCTGGAAGACTACCAGGAAGGGGAAGAGGATGATGTGACAGTAGCTACACCTTCCTTAGCTGTCCG  
TTGCGACTATGACCGTTGCCGCCACCTGCAGGTGTCTGCCAGGAGCTGCAGAAGTTGGCCAGTAGCC  
TGCCTGTGCCAGGGCTCTCCAGGAAGATCAACAGCCAGAGCCTCCTCGCTGGGAGAAGTGCAAATAA  
TGGCTGAAGAAGGCTACGCAGTGGTTCAGTGGTGTCTCCCTTCTCCAGTCAGCCACTACTGGCTTCT  
GCTTTGGGAAAGCAACGGGGCTCCACAGAAGAGTGCCCTCTCAATGCTACAGTTGGAAGAGCAGAGCTG  
AAGGGACTAAAGCCTGGGGTTGCTTATGTCCTTTGCGTGGTGGCTGCTAATGACGCAGGTGAGAGCAATG  
TTCTGGGGCAGAAGTCGAGGGTCTGAGAAGTGGACTGGCCCTTCTTTGGGCCCTGTCGCAAGTTTAT  
CATGCCGCTAAGCCTGTTACCCTGGTCTATGCAGCCGTGGGAGTGGGCACAGCCTTAGCTCTGCTGAGC  
TGTGCAGCCCTGGTTTGGCATTCTGCCTTCGTGAGCGATGGGGTTGCCCCGACGTCAAGGTATGGCC  
AAGCGTCAGAAGCTCTC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG218278 representing NM\_001013019  
 Red=Cloning site Green=Tags(s)

MLGSLLLWLAAMTTSLVSQPQILTLEDYQEGEEDDVTVATPSLAVRCDYDRCRHLQVSCQELQKVGPPVA  
 CLCPGLSREDQQPEPPRLGEVQIMAEEGYAVVHWCAPFSPVSHYWLLLWESNGAPQKSAPLNATVRRRAEL  
 KGLKPGVAYVLCVVAANDAGESNVPGAIEVEGPNWTGPSFGPCRKFIMPPKPVTLVYAAVGVGTALALLS  
 CAALVWHFCLRERWGCPRRQGMQAQASEAL

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001013019

**ORF Size:** 717 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001013019.3</a>
<b>RefSeq Size:</b>	2583 bp
<b>RefSeq ORF:</b>	720 bp
<b>Locus ID:</b>	68852
<b>UniProt ID:</b>	<a href="#">Q3TYX2</a>
<b>Cytogenetics:</b>	19 A