

Product datasheet for **MG218296**

Slirp (NM_026958) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Slirp (NM_026958) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Slirp
Synonyms: 1810035L17Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG218296 representing NM_026958
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCTCCGCCATAAAGGGGCTAAGCGCACTGCGTTCAGTACTGGTCGGCCTATTGCGTTTGTCAG
GAAAAATTCCTTGGACCGCGGCAGCGAGTGAGCTGAGAGAACATTTTGCACAGTTTGCCATGTAAGAAG
GTGCACTGTACCTTTGACAAAGAGACTGGCTTTCACAGAGGCATGGGTTGGTTTCAGTTTCTCTCAG
GAAGAATTCAGAATGCACTACAACAAGAACATCATATTATTGATGGAGTAAAAATCCACGTTCAAGCTC
AAAGAGCAAAGGCTTTGCACGGAGCTCAAACATCTGATGAAGAAAGATTCTGAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG218296 representing NM_026958
Red=Cloning site Green=Tags(s)
MAASAIKGLSALRSSTGRPIAFVVRKIPWTAASELREHFAQFGHVRRTVPFDKETGFHRGMGWQFSSQ
EELQNALQQEHHIIDGVKIHVQAQRAKALHGAQTSDERFLR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_026958

ORF Size: 336 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026958.3</u> , <u>NP_081234.2</u>
RefSeq Size:	427 bp
RefSeq ORF:	339 bp
Locus ID:	380773
UniProt ID:	<u>Q9D8T7</u>
Cytogenetics:	12 D2
Gene Summary:	RNA-binding protein that acts as a nuclear receptor corepressor. Probably acts by binding the SRA RNA, and repressing the SRA-mediated nuclear receptor coactivation. Binds the STR7 loop of SRA RNA. Also able to repress glucocorticoid (GR), androgen (AR), thyroid (TR) and VDR-mediated transactivation (By similarity).[UniProtKB/Swiss-Prot Function]