

## Product datasheet for **MG218411**

### Ociad2 (NM\_026950) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ociad2 (NM\_026950) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ociad2  
**Synonyms:** 1810027I20Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG218411 representing NM\_026950  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCCTCAGTGTCCACTCATGGGAACCAAGAAAAAGCCCCACTTGCCGCCACTGAGCAAGCAGAGCC  
TGTTGTTTTGCCCAAATCAAACACTGCACATCCACAGAGGAGAGATTGCCAAGATTATCAGAGAGTGTC  
AGAAGAGAGTTTCTGGAAGAGAGCCCTGCCTTTTTCTCTTATAAGCATGCTTGTCACCCAGGGACTCGTC  
CACCAAGTTATTTAGCCGCTAATCCAAGATTTGGATCATTGCCTAAAGTTGCATTGGCTGGTCTCTTGG  
GATTTGGCCTTGGGAAGGCTTCTATATAAGAGTGTGCCAGAGTAAATTCATTCTTTGAAGATCAGCT  
CCGTGGTCTGGCTTTGGTCTGAGCATAACAGGCACTGCCTGCTTACCTGTGAGGACTGTAAAACCAGG  
CGTGGATTAAGCGAGAAGGCAGGGTCACAACCATCTGCTTCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG218411 representing NM\_026950  
Red=Cloning site Green=Tags(s)

MASVSTHGNQEKSPHLPPLSKQSLLLFCPKSKLHIHRGEIAKIIRECQEESFWKRALPFSLISMLVTQGLV  
HQGYLAANPRFGSLPKVALAGLLGFGLGKASYIRVCQSKFHSFEDQLRGAGFGPEHNRHCLLTCEDCKTR  
RGLSEKAGSQPSAS

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_026950.4, NP_081226.1</u>
<b>RefSeq Size:</b>	2201 bp
<b>RefSeq ORF:</b>	465 bp
<b>Locus ID:</b>	433904
<b>UniProt ID:</b>	<u>Q9D8W7</u>
<b>Cytogenetics:</b>	5 C3.2