

Product datasheet for **MG218580**

Mtcp1 (NM_001039373) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mtcp1 (NM_001039373) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Mtcp1
Synonyms: P13MTCP1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG218580 representing NM_001039373
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAAGAGAGGATGTGGGGCTCCACCCGATCATCTCTGGGTTACCAAGAAGGAGTCTACCGCGACG
 AATATCAGCGCACATGGGTAGCTGTGGTGAAGAGGAGACAAGTTTCTAAAGGCTCGAGTCCAGCAAGT
 TCAGGTTCTTTAGGTGACGCACTAAGCCAAGTCACCTTCTTACCTCCCACTACCTCTCATGTGGCAA
 CTCTACCCAGAGGAGCGCTACATGGATAACAACTCTCGCTTGTGGCAGATTGAGCATCATTTAATGGTCA
 GGGGAGTACAGGAACACTCCTTAACTTTGCCTGATGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG218580 representing NM_001039373
 Red=Cloning site Green=Tags(s)

MAREVDGAPPDHLVWHQEGVYRDEYQRTWVAVVEEETSFLKARVQVQVPLGDATKPSHLLTSQLPLMWQ
 LYPEERYMDNNSRLWQIQHLMVRGVQELLLKLLPDD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001039373

ORF Size: 321 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001039373.5](#), [NP_001034462.1](#)

RefSeq Size:	1594 bp
RefSeq ORF:	324 bp
Locus ID:	17763
UniProt ID:	Q60945
Cytogenetics:	X A7.3
Gene Summary:	Enhances the phosphorylation and activation of AKT1 and AKT2.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG218580