

Product datasheet for **MG218759**

Rad52 (NM_001166381) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Rad52 (NM_001166381) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Rad52 |
| Synonyms: | Rad52yh |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG218759 representing NM_001166381 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGGCCTGAAGAAGCAGTCCACAGAGGGTGTGACAACCATCCTCCCTTTGTTGGTGGGAAGTCTG
TGCTCCTCTTTGGACAGAGCCAGTATACAGCGGATGAATACCAGGCCATCCAGAAAGCTCTGAGACAGAG
ACTGGGTCAGAGTACATTAGCAGCCGCATGGCTGGAGGAGGTCAGAAGGTGTGTATATTGAAGGTCAT
CGGTAATTAATCTGGCCAATGAGATGTTTGGTTACAATGGCTGGGCACACTCCATCACCCAGCAAAATG
TGGATTTTGTGACCTCAACAATGGCAAGTTCTACGTGGGAGTCTGTGCATTTGTAAGGTGCAGTTAAA
GGATGGTTCCTATCATGAGGACGTGGGCTATGGAGTTAGTGAGGGCCTCCGTTCAAAGGCCTTGTCACTG
GAGAAGGCCAGGAAGGAGGCTGTGACTGATGGGCTGAAGCGAGCACTCAGGAGTTTTGGGAATGCACTTG
GAACTGTATTCTGGACAAAGACTATCTGAGGTCACTAAATAAGCTTCCACGACAGCTCCCTCTTGATGT
GGATTTAACTAAAACAAAGAGAGAAGATTTTGAACCATCTGTAGAACAGGCAAGATACAACAGCTGCCGA
CAGAAATGAAGCATTGGGACTCCCCAAACCACAGGAAGTGACTTCCCCTTGCGATCAAGCCCCCACATG
ACTCGAACATTAAGCTTCAGGGGCTAAGGACATCAGCAGCTCCTGCAGTCTGGCCGCCACCTTGGAGAG
TGACGCCACTCACAGAGGAAGCTCCGGAAGCTCCGGCAGAAGCAGCTCCAGCAGCAGTTCGGGAACAG
ATGGAGACTCGCCGGCAGAGCCATGCACCTGCCGAAGAGGTGGCAGCCAAGCATGCGGCAGTACTTCCAG
CCCTCCAAAACACAGCACCCCTGTAAGTGCAGCCTCAGAATCCTCCAGGAGAAAGTCGTCTTTCCAGA
TAACCTTGAAGAGAACCTTGAATGTGGGACCTTACTCCAGACTTAGAGGACATCATTAAGCCCTTGTGT
AGAGCAGAACCAGCCAACTTCTGCCACTCGAACCTTCAACAACCAGGACAGCGTCCCACATATCCATT
GCCATCAGAAACCACAAGAAAAGCCTGGACCTGGGCACCTGCAGACCTGCAACACCAACCAGCATGTTCT
AGGTAGCAGAGAAGACTCTGAACCTCATAGGAAGGCCAGGACCTGAAGAAAAGGAACTAGATCCATCC

ACGGTACGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >MG218759 representing NM_001166381
Red=Cloning site Green=Tags(s)

MAGPEEAVHRGCDNHPPFVGGKSVLLFGQSQYTADYQAIQKALRQRLGPEYISSRMAGGGQKVCYIEGH
RVINLANEMFGYNGWAHSITQQNVDFVLDLNNKGFYVGVCAFVKVQLKDGSYHEDVGYGVSEGLRSKALSL
EKARKEAVTDGLKRALRSFGNALGNCILDKDYLRSLNKLPRQLPLDLDLTKTKREDFEPSVEQARYNSCR
QNEALGLPKPQEVTSPCRSSPPHDSNIKLGAKDISSSCSLAATLESDATHQRKLRKLRQKQLQQQFREQ
METRRQSHAPAEVAAKHA AVL PAPPKHSTPVTAASELLQEKVVFDPNLEENLEMMWDLTPDLEDI IKPLC
RAEPAQTSATRTFNNQDSVPHIHCHQKPKQEPGPGHLQTCNTNQHVLSGREDSEPHRKSQDLKRRKLDPS

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2248_b05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001166381


ORF Size: 1260 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|---|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001166381.1 , NP_001159853.1 |
| RefSeq Size: | 1698 bp |
| RefSeq ORF: | 1263 bp |
| Locus ID: | 19365 |
| UniProt ID: | P43352 |
| Cytogenetics: | 6 56.86 cM |
| Gene Summary: | Involved in double-stranded break repair. Plays a central role in genetic recombination and DNA repair by promoting the annealing of complementary single-stranded DNA and by stimulation of the RAD51 recombinase (By similarity).[UniProtKB/Swiss-Prot Function] |

Product images:



Circular map for MG218759