

## Product datasheet for **MG218786**

### Vgll1 (NM\_133251) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Vgll1 (NM\_133251) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Vgll1  
**Synonyms:** Tondu  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG218786 representing NM\_133251  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGGAAATGAAGAAGAAATGGTGTCTGGATACCTGAAAGCAGTCAAAAACCGGTAAAGACAGAGTGGA  
ACGCCGGTTCTGTCATTTTCACCTACTTCGAAGGGGACATCAACAGTATGGTGGATGAACATTTCTCCAG  
AGCTCTGCGTAACCTCAAGAGGCCCCACAATGCAGTTCCTCAAGTCACAGTGACAATTTGACCCTCAAG  
AACGGGGTGAATAAAACACGGACGGACACCAGCCACACGGATAGAGAGCGCAAAGTCGAATCCAGTCAAA  
CTGAAAATTCAGGAGAACTGAAAGACGTGCCTGTGGAGGAGAAGGAACCAATAAGAGCATGCCCCCGAA  
TCAGAGGCATCTCTTCTTGGACAATTGCACAGCCAGCAGCATCTCTGAAAATGGTGCCAGAAGCAGC  
AGCTTGGATGAATATGGTCTTAAGGCTATGGATCAGCACTCACTGTCCATGCCAAGACCCCTCTGCAC  
TTCTCTGGAGCCGTGGCACTTCTCTCCCTAGCGAGGCCAGGCTTCATAGCGCCTGCCTACTTCCCTGT  
CTTTCCCGACAGGCATCTGACTCCAGAGGTCTACCGTGTCTTTCCCGACAGGCATCTGACTCCAGAGGTC  
TACCGTGTCTTCTGACAGGCATCTGACTCCAGAGGTCTACCGTGTGGAGATGTGGGCCCTTCAAC  
ATTTAGTCCAGCAAGATAGATACCAAACCATCTCTGGAACCTGCTGCCAGGGAGAAGTGCAGCCCTGC  
CAAGATAGCTGGAAGCACAGGATCGTCTATGAACCTGCCTCCCTACCTGTCCACTACAAGAAGAAAATC  
TGTGCTCGTGGAGCTGCCAGTGCCAGTCTTGATGATGAAAAGAAGCCCATCTCAAAGCGAAGAAGAGATG  
CCTACTATTAT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG218786 representing NM\_133251  
Red=Cloning site Green=Tags(s)

MEEMKNGVWIPESQKPVKTEWNAGSVIFTYFEGDINSMVDEHFSRALRNLKRPHKCSSSSSHDNLTLK  
 NGVNKTRTDTSHDRERKVESSQTENSSELKDVPEEKGTNKSMPNQRHLSSWTIAQPAASPENGARSS  
 SLDEYGLKAMDQHSLMPKTPSALPLEPWHFSSLARPGFIAPAYFPVFPDRHLTPEVYRVFPDRHLTPEV  
 YRVFPDRHLTPEVYPDGRCGPLQHLVQDDRYQNHLEPAARENCSPAIA GSTGSLMNLPPYPVHYKKKI  
 CARGAASASLDDERSPSPKRRRDAYYY

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_133251

**ORF Size:** 921 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_133251.1](#), [NP\\_573514.1](#)

**RefSeq Size:** 1187 bp

**RefSeq ORF:** 924 bp


**Locus ID:** 170828

**UniProt ID:** [Q99NC0](#)

**Cytogenetics:** X A6

**Gene Summary:** May act as a specific coactivator for the mammalian TEFs.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG218786