

## Product datasheet for **MG218877**

### Scfd2 (NM\_001114660) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Scfd2 (NM_001114660) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Scfd2
Synonyms:	E430013M20Rik; STXBP1L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG218877 representing NM\_001114660  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGAGCGCCGCGGCATACTGGCCTTTGCCAGCAAGGCTGGGAGCAGGTGCTGGCCAAAGTGAAGTGGT  
 CCGTGGTCTACCTGGACGCTGCCTGCGCCGAGAGCCTGCACTGGAGCTGCGGATCCAGCCGGCTTCTGGA  
 GCGGTTGAAGGGTCTGCGTGTACCTGCGGGAATTTGAACCCAGGCCATTGGCGGGGAGCCAAGCAG  
 CCCCAGCGGTTTTCTGTTTTGAGCAGCCCGCTGAAAGGCCGATCGTAGACACTCTACAGAGTATCATT  
 GCCGCAGTCACTTTCAGCACTGTGTGGTGGTCACAGCCGTGAGCCACGCAGTCCACCTCACTGCTAACCA  
 CGTACCGGCCGCTGCAGCGGCCGAGCTCGAGGGGCAGCAGCCGGTGTTCGAACAGCTGGAGGAGAAGCTC  
 TGTGAGTGGATGGCAACGAGAATAACGCGCCGAGGTGCTCCACGTCCCCTGTTTCTCGCCCCTGTTG  
 CCTCCCATCTCGTTTTACTCCAGCCTTGTACCTCTTCCCCTGCTACCTCAGGATGCTCCACGCTCT  
 TAACAGTCCAGACCCGACAAGAGGAGGCTGAGCAGCTTGGGTGAGGTGGATGCCACTGCCCTAACCCCA  
 GAGCTGCTGCTGTACATCAGATGCCTCGTGTGAGGCTCAGTTCTCTGTGTGAACATCTAGGGGTACGAG  
 AGGAGTGTTCGCCGTGGGTCCCCTCAGTCGGGTATCGCTACAGATCTGGCAAATTATGCCCCCGCCAA  
 GAACAGGAAAAAGACTGCTACAGGCAGGGCTTCTGTGGTCTTTGTGGATAGAACTCTGGATCTCACAGGA  
 GCGGTTGGACACCACGGAGACAACCTTAGTAGAGAAGATCATGTCAGTGTTCACAGCTTCCAGGCCACA  
 CCCACGACGTGATGGTCAACATGGCCGAGCTCACTGCTGTCCAGGCTGTGGAGGAAAACAGAACGTGAT  
 TGCACCAGGCTGCCTTGCCTTCAATGAGGCCCTCAGCCAAGGCCCTGTGGGAAGCCTTACTGAACTCC  
 AAGCACAAGAGGAGTGTGGAAGTTCGAGACATCTAGTGAAGCAGCCAGCAGAGAGAACCTGCCCA  
 TCAAGATGAGTATGGGAAGAGTCACTCCAGGGCAGCTCATGTCCATATTTCAACTTTCAAGAACAACCT  
 CAGAGCTCTTAGAAATCACTGTGGCTTCTACAGCTTGGGATGGCCACAGTTCAAACCTCTGAAGCACCCA  
 CAGACTGCCAAGTGGATAAATTCTGGCATTGAAAGGCTCCTTCTACAGAGCCTTGGGATTCAACCA  
 TGGCTGGTGTGCTAAATCAGCTGCTGCCAATGATTAAGTCTCAAGCCAGAGAACCAGTACGACCTCAA  
 CCCTGAGGAACTGCTGATCCTGCTCATATACATTTACTCTGTCCCTGGAGACGTACGCTGGACAGAGAC  
 CTGGGAGACGTGGAAGAAAAGGTGAAGAAAGCATTGGCCCATGCTCTCTGAGGAGTCTGAGCTGTCCC  
 CTTTGTCTCAGAAAATCACAGGCTGTGACTCTGCAGTTGACCTGACGTTGCCAAATCCCAAATTTGCTGT  
 GAATGATGTCTTCATGGCTCTCCGGGAGATCGCTGGAGCTAGGAATCTCATGAGACAGTTTAAAGTCTGTG  
 TATGTTCTGGAATAACACCCATCAGGCATCCTATAAGCCGCTGTTGAAGCAGGTTGTGGAGGAAATAT  
 TTAATCCTGAGAAGTCAATCCTATTGACATTGAACATATGTCCTCAGGCCTCACTGATCTCCTGAAAAC  
 TGGATTACAGATGTTTATGAAGGTGAGCCGGCCCATCCAGTGACCATCCCCTCTTGATCCTCTTTGTG  
 GTGGGCGCGCTCACAGTCGCTGAAGCCAAAATGGTCAAAGACCTTGTAGCATCTCTGAAGCCAGGAACCC  
 AGGTGATGGTGTGTCCACAAGACTCCTGAAGCCGCTTAACATTTCCCGAGCTGCTCTTGAACCTGACCC  
 GCTACACCCAGACCTCGGCTTC

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG218877 representing NM\_001114660  
 Red=Cloning site Green=Tags(s)


MSAAGILAF AQQGW EQVLAKVKWSVVYLDAACAESLHWSGSSRLLEAVKGPACHLREFEPQAIGGGAKQ  
 PRAVFVLS SPLKGRIVDTLQSIICRSHFQHCVVVTA VSHAVHLTANHVPAAAAAELEGQQPVFEQLEEKL  
 CEWMGNENYTAEVLHVPLFLAPVASHLAF TPAFATLFP LLPQDVHALNSARPDKRRLSSLGEVDATALT  
 P ELLLYIRCLV SGLSSLCEHLGVREECF AVGPLSRV IATDLANYAPAKNRKKTATGRASVVFVDR TLDLTG  
 AVGHHGDNLVEKIMS VLPQLPGH THDVMV NMAELTAVQAVEENQ NVIAPGCLAPSNEASAKALWEALLNS  
 KHKEAVMEVRRHLVEAASRENLP IKMSMGRVTPGQLMSYIQLFKNNLRALRNHCGLLQLGMATVQTLKHP  
 QTAKWDNFLAFERLLLQSLGDSTMAGV LNQLLPMIKSSSQRTSDDL NPEELLILLIYIYSVPGDVTLDRD  
 LGDVEEKVKALAHV LSEESL SPLLQKITGCDSAVDL TPKSQI AVNDVFMALREIAGARNLMRQFKSV  
 YVPGNNT HQASYKPL LKQVVEE I FNPEKSDPIDIEHMSSGLTDLLKTGFSMFMKVS RPHPSDHPLLILFV  
 VGGVTVAEAKMVKDLVASL KPGTQVMV LSTRLLKPLNIPELLFATDRLHPDLGF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM\_001114660

ORF Size: 2052 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001114660.2](#), [NP\\_001108132.1](#)

**RefSeq Size:** 3687 bp

**RefSeq ORF:** 2055 bp


**Locus ID:** 212986

**UniProt ID:** [Q8BTY8](#)

**Cytogenetics:** 5 C3.3

**Gene Summary:** May be involved in protein transport.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG218877