

Product datasheet for **MG219036**

Tmem30c (NM_027651) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem30c (NM_027651) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tmem30c
Synonyms:	4933401B01Rik; 4933409A18Rik; Cdc50c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG219036 representing NM_027651 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGATGATGCCCCAGTATGACCTGTCCAGACTGCCAGAAAACACAGCATTGAAACAGCAAACGTTGC
CCACCCAGCAGCTAAATCTATCTGCCTCGGTGGTCTCTATCTTTTCATCACAGGAGGATTCTGCCT
CTCTATAGGCATCATCCTTCTGTGTCTGCAAAGAGCACCAAGAAAATTGAGATAAATTACAAAAGACG
TGTGCGAATTGCGCTCAGCTCCGGGAAAATAGTTCTAATTTTGACAAAGCATGCAACTGTTCTTCTCCTT
TTTACCTCCGGAGAAAATGGAGGGTATGTTTATATGTACTACAACTGTATGGTTTCTATCAGAACTC
CTATCAATATATCTTATCAAGAAGCAATAGCCAATTGGTGGGTAAGGATATTTGGGATACTACGAACTGT
GATCCGTTCCAAGTGTCCACAATGACACTCCCATCATTCTTGTGGTGAATTGCCAACAGCATATTCA
ATGATACCATAACTCTTTCATACAACCTGAATTCATCTACACAGATAGAAGTCCCAATGCTAAAGAGTGG
GCTTACATGGTGGACAGATAAGTATGTCAAATTCGGAATCCAAGATCCAGTAACCTCACAGTACATTT
GCAGGAAGTTCGAAGCCTCTGCACTGGGCCAAACCTATCTATGAACTGGATCTAGATGACCCAGGGAACA
ATGGCTTTCTCAATGAGGACTTCATTGTGTGGATGCGGACGGCTGCCTTTCCACTTTCAAAAACTGTA
CCGCCGACTGAAGCGGTGCACGCCTTCGCTGAAGGCTTACCTGCTGGTAACTACAGTTTGAGCATCTCC
TACAACCTTCCAGTAACCATGTTCCAAGGAGAAAAGTCCATTGTCCTTTCCACTGACATGGATCGGGG
GAGGAGGCCTTTTCTTGGGGCTGACTTACACAGTGACGGGGCACTGACGCTGTTGGCCTCTTTTGCCAT
CCTGACGATTCATTGATGCTGAAAAGGAGTAACTGAATTTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG219036 representing NM_027651
 Red=Cloning site Green=Tags(s)

MEMMPQYDL SRLPENTALKQQTLP TQQLNL SASVLSIFFITGGFCLSIGIILLLSAKSTKKIEINYTKT
 CANCAQLRENSNFDKACNCSLPFYLPEKMEGDVYMYKLYGFYQNL YQYILSRNSQLVGKDIWDTTNC
 DPFQVSHNDTPIIPCGAIANSIFNDTITLSYNLNSSTQIEVPM LKSGLTWWTDKYVKFRNPRSSNFTSTF
 AGSSKPLHWAKPIYELDLDDPGNNGFLNEDFIVWMRTAAFP TFKKLYRRLKRVHAF AEGLPAGNYLSIS
 YNFPVTMFQGEKSIVLSTLTWIGGGGLFLGLTYVTGAL TLLASFAILTIHMLKRSKLNFL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_027651

ORF Size: 1026 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027651.1 , NP_081927.1
RefSeq Size:	1382 bp
RefSeq ORF:	1029 bp
Locus ID:	71027
UniProt ID:	Q9D4D7
Cytogenetics:	16 C1.1