

## Product datasheet for **MG219078**

### **Armc8 (NM\_028768) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Armc8 (NM_028768) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Armc8
Synonyms:	1200015K23Rik; Gid5; HSPC056
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>MG219078 representing NM\_028768  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCGTGCCTGTTGGAGACCCCAATCCGCATGAGCGTCTCTCGGAGGTAACAGCTAGCAGTCCGCACT  
 ATGTTGACAGGCTATTTGACCTGATCCCAGAAAAGTTCTCCAAGGTGTCATAGACATGAAAAATGCTGT  
 AATTGGAAAACAATAAACAGAAAGCCAACCTCATCGTTCTAGGAGCGGTTCCAAGATTGTTGTATTTGCTT  
 CAACAAGAAACCTCAAGTACAGAGCTGAAAAGTGAATGTGCAGTGGTGTAGGGAGTCTTGCGATGGGCA  
 CAGAAAACAATGTCAAGTCTCTACTAGACTGCCATATTATCCCTGCCTTACTACAAGGACTACTGTCCCC  
 AGACCTGAAGTTCATTGAAGCTTGCCCTCCGATGCCTCCGGACAATCTTACCAGCCCTGCTACTCCAGAG  
 GAGCTACTGTATACAGATGCTACGGTGATACCACACCTCATGGCACTGCTGAGCAGGTCCCGCTACACCC  
 AGGAGTATATCTGCCAGATCTTCTCACACTGTTGCAAAGGGCCAGATCATCAAACAATTTTGTAAACCA  
 TGGGGCAGTGCAGAACATTGCTCATCTCCTAACCTCACCTCTTACAAAGTTCGAATGCAAGCACTGAAA  
 TGCTTCTCTGTTTTAGCCTTTGAGAATCCCAAGTATCGATGACCCTGGTAAACGTTTTGGTTGATGGGG  
 AGTTGTTACCACAGATATTTGTGAAGATGTTACAGAGGGATAAGCCTATTGAGATGCAGCTCACGTCAGC  
 AAAATGTTAACTTACATGTGCAGAGCTGGAGCAATTCGCACAGATGACAGCTGTATTGTATAAAACT  
 TTGCCCTGTTTGGTTTGAATGTGCAGTAAGGAGAGGTTACTAGAGGAGAGAGTTGAAGGAGCTGAGACCC  
 TCGCCTATCTGATTGAGCCAGATGTTGAGCTGCAGAGAATCGCTAGCATAACCGACCACCTCATTGCCAT  
 GCTGGCTGACTATTTCAAGTATCCCAGCTCAGTGCAGGCCATCACTGACATTAAGGCTTGATCATGAC  
 TAAAAACATGCCCATGAGCTCCGCCAGGCTGCGTTCAAGCTCTATGCCTCTCTGGAGCAAATGATGAAG  
 ACATCCGGAAGAAGATCATTGAGACTGAAACTATGATGGACCGAATTGTGACTGGCTTGCTGAGTCTAG  
 TGTCAAGGTGCGGTTAGCTGCTGTCAGATGTTTGCACAGTTTATCCAGATCTGTACAACAGCTTCGAACC  
 AGTTTCCAGGATCATGCGGTGTGAAAACCTTTAATGAAGTTTTACAAAATGCACCCGATGAAATCCTAG  
 TAGTAGCATCTTCCATGTTATGTAACCTTCTTCTGAATTTTCTCCAAGTAAAGACCAATTTTGAATC  
 AGGAGCTGTCGAGCTGCTCTGTGGGTTGACACAGAGTAAAAACCTGCTCTACGAGTGAATGGGATTTGG  
 GCTTTAATGAATATGGCATTTCAGGCGGAACAAAAATAAAAGCAGATATTTACGAAGTCTGAGTACTG  
 AACAGTTATTCGGTTATTGTCAGATTCTGATATGAATGTGCTGATGAAGACATTGGGGCTTCTTAGAAA  
 TCTCCTCTCCACTCGGCCTCATATAGATAAAAATCATGAGTACTCATGGGAAGCAGATTATGCAAGCTGTC  
 ACTCTGATACTGGAAGGGGAGCATAGCATTGAAGTCAAAGAGCAGACCCTGTGCATCCTAGCCAACATTG  
 CAGATGGGACAACAGCAAAGAGCTTATTATGACCAATGATGACATCCTGCAGAAAATCAAATATTACAT  
 GGGTCATTACATGTCAAAGTGCAGCTTGTGCCATGTTTTGTATATCAAACCTCATATGGAACGAAGAG  
 GAAGGTTACAGGAACGCCAGGATAAATTACGAGACATGGGGATTGTAGATATTCTACACAAACTGAGTC  
 AGTCGGCAGACTCAAACCTCTGTGACAAGGCAAAGACGGCGCTGCAGCAGTACCTGGCG

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:** >MG219078 representing NM\_028768  
 Red=Cloning site Green=Tags(s)


MACLLETPIRMSVLSEVTASSRHYVDRLFDPDPQKVLQGVIDMKNNAVIGNNKQKANLIVLGAVPRLLYLL  
 QQETSSSTELKTECAVVLGSLAMGTENNVKSLLDCHIIPALLQGLLSPDLKFEACLRLCLRTIFTSPVTP  
 ELLYTDATVIPHLMALLSRORYTQEYICQIFSHCCKGPDHQTILFNHGAVQNIHLLTSPSYKVRMQALK  
 CFSVLAFENPQVSMTLVNVLVDGELLPQIFVKMLQRDKPIEMQLTSAKCLTYMCRAGAIRTDDSCIVLKT  
 LPCLVRMCSKERLLEERVEGAETLAYLIEPDVELQRIASITDHLIAMLADYFKYPSSVSAITDIKRLDHD  
 LKHAHEL RQA AFKLYASLGANDEDIRKKI IETETMMDRIVTGLSESSVKVRLA AVRCLHSLSRVQQLRT  
 SFQDHAVWKPLMKVLQNAPDEILVVASSMLCNLLLEFSPSKEPILESGAVELLCGLTQSENPALRVNGIW  
 ALMNMAFQAEQKIKADILRSLSTEQLFRLLSDSDMNVLMKTLGLLRNLLSTRPHIDKIMSTHGKQIMQAV  
 TLILEGEHSIEVKEQTLCILANIADGTTAKELIMTNDLILQKIKIYYMGHSHVKLQLAAMFCISNLWNEE  
 EGSQERQDKLRDMGIVDILHKLQSADSNLCDKAKTALQQYLA

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



## Plasmid Map:



ACCN: NM\_028768

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_028768.3](#), [NP\\_083044.2](#)

RefSeq Size: 4628 bp  
RefSeq ORF: 2022 bp  
Locus ID: 74125  
UniProt ID: [Q9DBR3](#)  
Cytogenetics: 9 E3.3  
Gene Summary: Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1.[UniProtKB/Swiss-Prot Function]