

## Product datasheet for **MG219174**

### Gpr160 (NM\_027965) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gpr160 (NM\_027965) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Gpr160  
**Synonyms:** 1700025D19Rik; GPCR150  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG219174 representing NM\_027965  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACAGCTCTTTCTTCAAAGAAGTCTCTCTTCAATACCAGCTCCATCAGTCTCCGAGCTCCTGGAGG  
CCTCCTGTCTGTTGTTCTTGATCATCTGGGAAAGTGTATTAACATCCTCCTTCTAAGAGTGAGAAG  
GGGAGACGCCGCTGGACTTTGATGGAGTACTTCTGCTTCTCACTCGCACTTGTCGATCTCCTGCTGTTG  
GTAACATTTCCATCCTGACCTACTTCAGGGACTTTGTAGTTCTAGGTATTAGGTTTACCAGATACCACA  
TCTGTCTGCTCACACAAATTATTTCTTCACTTATGGCTTCTTGCAATACCCAGTCTGCTACTGGCTTG  
TATAGACTACTGGTGTAGTCTGTCCAGAGCCAGCAAGCAGTCATCTAGGTGGCAGAAATTTATTTTC  
TTGACAGTCATTTAACTTGATCTCAGTCCTTGCTTACGTTTTGGTTGACCCGGCCATCTCAGTCAGCC  
TAAAGGCACACAGGGGCTATGTGTACCAGTGCCCGCCTATGTCAGCACTCAGAGTCACTGGCTCTCACT  
GTCTATGCTGATGGTTTTATTTGTGGCTTTTCTGATTTTCATGGCAGGAAGTTGTTGCCTTGCTTCAGGCT  
ATGAGGATAGCATCCTATAAGAGCAAGGCCGCTCTACTTTCTTTCCCACTCCACTGCGGTTACGCC  
TGAGCTGCAGAGAAGCACTCTTGCCAGGCTCATCGTGTGCTTTCTTGGACCTGGTTTCCCTTTGTAGC  
ACTTCAGTCCCTCATCTCTCCCTCAGAGTTTTCAGATCCCGGCATACATAGAAATGAACGTCCCCTGGCTG  
TACTTTGTCAACAGTTTCTCATCGCTGCAGTTTATTGGTTTAACTGCCACAAGCTTGATTTAAGAGACA  
GCTCGTTACCTGTGGATCCTTTTCATCAACTGAAATGCTGCTTTGTTCCAGTCCATAGACTTAAGCAAGT  
GGAGAGGCCCATGTCCATAGTCATCTGT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG219174 representing NM\_027965  
 Red=Cloning site Green=Tags(s)


MTALSSKNCSLQYLHQSPQLLEASCLLFLIILGKVLNILLRVRRGDARWTLMEYFCFSLALVDLLLL  
 VNISILTYFRDFVVLGIRFTRYHICLLTQIISFTYGFHLHYPVCSLACIDYWCSL SRASKQSSRWQKLLYF  
 LTVILTWISVLAYVLVDPAISVSLKAHRGYVYQCPAYVSTQSHWLSLSMLMVL FVAFLISWQEVVALLQA  
 MRIASYKSKAALYFPFPLHCGYALSCREALLPRLIVCFLGTWFFPVALQVLILSLRVQIPAYIEMNVPWL  
 YFVNSFLIAAVYWFNCHKLDLRDSSLPVDPFINWKCCFVPVHRLKQVERPMSIVIC

TRTRPLE - GFP Tag - V


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_027965

**ORF Size:** 1008 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_027965.1</a> , <a href="#">NP_082241.1</a>
<b>RefSeq Size:</b>	2018 bp
<b>RefSeq ORF:</b>	1011 bp
<b>Locus ID:</b>	71862
<b>UniProt ID:</b>	<a href="#">Q3U3F9</a>
<b>Cytogenetics:</b>	3 A3
<b>Gene Summary:</b>	Orphan receptor.[UniProtKB/Swiss-Prot Function]