

Product datasheet for **MG219239**

Acsl6 (NM_001033597) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Acsl6 (NM_001033597) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Acsl6
Synonyms:	A330035H04Rik; AW050338; FacI6; LACS; Lacsl; mKIAA0837
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG219239 representing NM_001033597
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTGACCTTCTTCTGGTGTGCGGGGTTCTCTCTGGCTGTTTCCGAGATTGCTCTCACTTCTGG
AGAAGATGCAGACCCAGGAGATCCTGAGGATCCTGCGGCTGCCGAGCTATCGGACTTGGCCAGTTTTT
CCGCAGCCTCTCAGCTACCACCCTCGTGAGTGTGGGTGCACTGGCTGCCGTCCTTGCCTACTGGCTCACT
CACCGGCCAAAGGCCTTGAACACCATGCAACCTCCTGAAGCAGTCGGAAGAAGTGAGGACGGTGGT
GAGCCCGCGATCTGTGATTGGGGTTGCACTCAATTGCTTACCCATTACTATGACGATGCCCGACCAT
GTACCAGGTGTTCCGCCGTGGGCTTAGCATCTCAGGGAATGGACCCTGTCTTGGCTTCCGGAAGCCAGAG
CAACCTTACCAGTGGCTGTCTACCAAGAGGTGGCCAAAAGGGCTGAATTTCTGGGGTGGGGCTTCTCC
AGCAGACTGTAAAGTGGGCACAGAGCAGTTTCGTTGGTGTTCACAAAATCGGCCGAGTGGATCAT
TGCCGAGCTGGCCTGCTACACGTATTCCATGGTGGTGGTACCGCTCTACGACACCCTGGGCCAGGGTCT
ATCAGCTACATCAATACTGCGGACATTTGCAGGTAATCGTCGATAAACCCACAAGGCAACACTTC
TGCTGGAACACGTGGAGAGGAAGGAGACTCCGGGCCTCAAGCTGGTCATCCTCATGGAGCCGTTTGGGA
CGCCCTGAGAGAGAGAGAAAGAAGTGCGGGTGGACATCAAGTCCATGCAGGCTATAGAGGACTGTGGC
CGAGAGAATCATCATGCCCCGTGCCCCACGGCTGATGACCTCTCATCGTGTGTTTCAAAAGTGGTA
CAACAGGGAACCCAAAGGTGCAATGCTCACCCACGGGAACGTGGTGGCCGATTTCTCGGGCTTCTGAA
AGTGACAGAGAGTCAGTGGGCCCCGACCTGTGCGGACGTCACCTTTTCTACTTGCCTTTAGCACACATG
TTTGAGCGCATGGTGCAGTCTGTTGTCTACTGCCACGGAGGCGCGTGGGCTTCTCCAGGGAGACATCC
GCCTCCTCTCAGATGACATGAAGGCTCTCCGCCCTACCATCTTCCCTGTGGTCCCACGGCTGTAATCG
GATGTATGACAAGATCTTCCACCAGGCAGACACCTCACTAAAGCGCTGGCTCCTGGAGTTTGCAGCAAAG
CGCAAGCAGGCAGAGGTCCGGAGCGGAATCATCAGAAACAATAGTATCTGGGATGAACTCTTCTTAAATA
AGATTCAGGCCAGTCTTGGTGGGCATGTGAGGATGATTGTCACTGGAGCAGCCCCGCGTCACCAACGGT
TCTGGGATTCTACGAGCAGCTCTGGGGTCCAGGTCTATGAAGTTATGGGCAAACTGAATGCACAGCT
GGGTGTACCTTCAACGCCAGGGGACTGGACATCAGGGCATGTAGGGCACCTTGCCTGCAACCACA
TCAAGCTGGTCGATGCAGAGGAACCACTACTGGACCTGCAAAGGAGAAGGAGAGATATGTGTAAAGG
ACCAAATGTGTTCAAAGGCTACTTAAAAGACGAGGACAGGACAAAGGAGGCCCTGGACAGCGACGGCTGG
CTTCACACTGGAGACATTGGGAAATGGCTGCCGAGGGAACACTCAAATCATTGATCGAAAAAGCACA
TATTTAAACTTGCTCAGGGGGAATATGTTGCGCCAGAGAAGATCGAGAACATCTACATCCGGAGTGAGCC
TGTGGCACAATCTACGTCCACGGGACAGCTTAAAGGCCTTTTTGGTTGGCATTGTCTGCCTGACCCCT
GAAGTCATGCCTTCTGGGCTCAGAAGAAAGGAATCGAAGGGACCTATCAGGAACTCTGCATGAAAAAGG
AATTGAAGAAAGCCATTCTGGATGACATGGTGTGCTGGGAAAGAAAGCGGACTGCATTCTTTGAACA
GGTTAAAGCCATTTACATCCATTGTGACATGTTCTCTGTTCAAAATGGTCTGCTGACACCAACACTAAAG
GCTAAGAGACCGGAGCTGAGAGAGTACTTCAAAAAGCAAATAGAAGAGCTTTACTTAGTCTCCGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG219239 representing NM_001033597
 Red=Cloning site Green=Tags(s)


MLTFFLVSGGSLWFLAEIALSLLEKMQTQEILRILRLELPELSDLGQFFRSLSATTLVSVGALAAVLAYWLT
 HRPKALQPPCNLLKQSEEVEDGGGARRSVIGGCTQLLTHYYDDARTMYQVFRRLSISGNGPCLGFRKPE
 QPYQWLSYQEVAKRAEFLGSGLLQHDCKVGTQFVGVFAQNRPEWIIAELACYTYSMVVPLYDTLGPGS
 ISYIINTADICTVIVDKPHKATLLEHVERKETPGLKLVLMEPFEDALRERGGKCGVDIKSMQAIEDCG
 RENHHAPVPPRPDDLIVCFSTGTTGNPKGAMLTGHNVVADFSGFLKVTESQWAPTCADVHFSYLPLAHM
 FERMVQSVVYCHGGRVGFQGDIRLLSDDMKALRPTIFPVVPRLLNRMVDKIFHQADTSLKRWLLEFAAK
 RKQAEVRSIIRNNSIWDELFFNKIQASLGGHVRMIVTGAAPASPTVLGFLRAALGCQVYEGYGTCTECTA
 GCTFTTPGDWTSGHVGAPLPCNHIKLVDAEELNYWTKGEGEICVKGNVFKGYLKDEDRTKEALDSDGW
 LHTGDIGKWLPEGTLKIIDRKKHIFKLAQGEYVAPEKIENIYIRSEPVQIYVHGDLSKAFVLGIVVDPDP
 EVMPSWAQKKGIEGTQYELCMKKELKKAILEDVMLGKESGLHSFEQVKAIYIHCDMFSVQNGLLTPTLK
 AKRPELREYFKQIEELYLVSV


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001033597

ORF Size: 2166 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033597.1](#), [NP_001028769.1](#)

RefSeq Size: 2592 bp

RefSeq ORF: 2169 bp

Locus ID: 216739

Cytogenetics: 11 32.13 cM

Gene Summary: Activation of long-chain fatty acids for both synthesis of cellular lipids, and degradation via beta-oxidation. Plays an important role in fatty acid metabolism in brain and the acyl-CoAs produced may be utilized exclusively for the synthesis of the brain lipid (By similarity).
[UniProtKB/Swiss-Prot Function]