

Product datasheet for **MG219399**

Cd72 (NM_007654) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cd72 (NM_007654) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cd72
Synonyms:	CD72c; Ly-19; Ly-32; Ly-m19; Lyb-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG219399 representing NM_007654 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGACGCTATCACGTATGCAGACCTGCGCTTTGTGAAAGTGCCCTGAAGAACAGCGCATCTAACC
ATCTAGGACAGGACTGTGAGGCCTATGAAGATGGGAACTCACCTACGAGAATGTGCAAGTGTCTCCAGT
CCCAGGAGGGCCACCAGGCTTGGCTTCCCCTGCACTAGCGGACAAAGCAGGGTTCGGGTGAGCAACCA
ACTGCGACCTGGAGCTCTGTGAACCTGCTGCTCTCAGGCAGATTCCCCGCTGTCCTACAGTCTGCTTGC
AATACTTCTTGGCTTGGCCTTCTCGTGTCTGTCTGATGTTAGGGTGGCTGTCATCTGCCTGGGAGTTTCG
CTATCTGCAGGTGCTCGGCAGTTCAGGAGGGGACCAGGATTTGGGAAGCCACCAATAGCAGCCTGCAG
CAGCAGCTCAGGAGAAGATAAGTCAGCTGGGGCAGAAGGAGGTGGAGCTTCAGAAGGCTCGGAAAGAGC
TGATCTCGAGCCAGGACACATTACAGGAGAAGCAGAGGACTCACGAGGACGCTGAGCAGCAACTACAAGC
CTGCCAGGCTGAGAGAGCGAAGACCAAGGAGAACCTGAAAACCTGAGGAGGAGCGGAGGAGGGACCTGGAC
CAGAGGTTGACAAGCACGCGGGAGACACTGAGGCGCTTCTTCTGATTTCATCAGACACCTGCTGTCCAT
GCGGATGGATTCCATATCAGGAAAGGTGCTTTTACATCTCACATACCCTCGGAAGTCTGGAGGAGAGCCA
AAAATACTGCACATCTCTGTCTCCAACTGGCAGCATTTCGATGAACCTTCTAAGTATTACTATGAAGTT
TCTCTGCCAGCGCTTAGAGGAGTTGCTAGATCGTTCGAAGTCATATTGGATACAGATGAGCAAGAAGT
GGAGGCAGGACTCTGACTCTCAAAGCCGACATTGTGTCAGGATAAAAACATATTACCAGAAGTGGGAAAG
AACAATTTCCAAGTGTGCAGAGCTTACCCTGCATTTGTGAGTCGGAGGCTTTCAGTTTCTGATGGG
ATCAATCTGAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

ORF Size:	1062 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_007654.2 , NP_031680.2
RefSeq Size:	1482 bp
RefSeq ORF:	1065 bp
Locus ID:	12517
UniProt ID:	P21855
Cytogenetics:	4 23.04 cM
Gene Summary:	Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]