

Product datasheet for **MG219509**

Col25a1 (NM_029838) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Col25a1 (NM_029838) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Col25a1
Synonyms:	2700062B08Rik; CLAC-P
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG219509 representing NM_029838
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTGGTGAAGAAGCTTGCAGGAAAGGAGGGGGACGAGAGTCTGGATCAGAAGATCCGCGCCCTTGG
GACAGCGTTGTGCCGGCACCATGCCCTCGTGACGGCCCTGGCGACCCTTTGTCAAGTGGTTGCTGTGGC
TTTCTGTTTTTATCTTGGGGTAAAAACCAACGACCTCCAGGCGAGGATTGCTGCTTTGAATCTGCTAAA
GGGACCCCTTCTTCCATCCGCTGTCTGACACCGTGGATGAGCTGAAGGCAATGGTTCAGGAGAAAGTGG
AGCGTCTCTTGGCTCAGAAATCTACGAACATATGGCTAAAATAAGAACGGCCAGGGAGGCACCTTTAGA
GTGCAACTGCCAGCAGGTCTCCAGGAAACGAGGGAAGAGAGGCCGAAGAGGAGAATCTGGTCTCTCT
GGTCAGCTGGTCTCAGGGCCCTCTGGTCCAAAAGGTGATAAGGGAGAACAAGGTGATCAGGGACCTC
GGATGGTGTTCCTAAAATCAATCACGGCTTTCTCTGCTGATCAGCAGCTCATTAAACGCCGGCTGAT
TAAGGGTGACCAAGGACAGGCAGGGCCTCCAGGACCTCCAGGCCCTCTGGTCCAAGAGGCCACCTGGG
GACACAGGAAAGGACGGCCCCGAGGAATGCCAGGAGTACCTGGTGAACAGGAAAACAGGAGAAACAAG
GCTTGATGGGACCTCTGGGGCCTCCAGGACAAAAGGGTCCATTGGAGCACCTGGGACCCAGGCATGGA
TGGGCAAAAGGGTGAACCTGGATCACCTGGAGCAGCCGGGAGAGTGGACTACCAGGACCTAAGGGAGAA
CCTGGAAGAAGGAGAAAAGGGAGATGCTGGAGAAAATGGTCCCAAAGGTGATACAGGAGAAAAGGGTG
ACCCTGGATCATCTGCTGCAGGAATTAAGGGAGAACCTGGAGAATCTGGCCGCCCGGGGAGAAAGGGTGA
ACCAGGGCTGCCTGGCTGCCTGGACTTCCGGGAATAAGGGAGAACAGGCTTCATTGGTCTCAAGGA
AACGAGGAGATCCTGGAGCCCCGGGGCAAGGGGAAGCAAGGTGAATCAGGAGCTAGAGGCCGAAGGG
GTCAAAGGGTATCGTGGAGACAAAGGAGACTCTGGCGCTCTGGGACCACGGGGTCCACCTGGACAAAAG
GGGATCCAGGAGCCACAGAGATCATAGACTACAATGGCAACCTCCATGAGGCCTTACAGAGAATTACCA
CCTTAACTGTACGGGCCCCCTGGACCTCCTGGACCTCAAGGACTACAAGGGCCAAAGGGTGAAGCAAGG
CTCTCCAGGAATCCCGGAGTTGATGGAGAACAGGGACTCAAAGGCTCCAAGGGAGACATGGGGACCCA
GGTGTGCCAGGTGAAAAGGAGGACTGGGACTTCTGGATTGCCGGGTGCCAATGGAGTAAAAGGAGAGA
AAGGAGACACCGTTTGGCAGGTCTCAGGGCCTTCTATCATAGGCCACCAGGCCCTCCAGGTCCCCA
TGGCCACCTGGTCCATGGGGCCCATGGACTCCTGGACCAAGGGAGCATCTGGCTTAGACGGAAG
CCAGGATCCCGGGTGCAGATGGTCTATAGGACCCACGGCCCTGCAGGACCAAGGAGAAAAGGAGGAG
AGAAAGGAGCTATGGGAGAGCCTGGACCCAGAGGGCCCTATGGGCTGCCTGGCAAAGATGGAGAACCTGG
TCTTGATGGCTTCCCTGGTCTCGAGGCGAGAAGGGTGAACCTGGGAGAAAAGGGAGAAAAGGGATTCCGT
GGCGTTAAGGGGGAAAAGGGGGAGCCAGGCCAGCTGGCCTGGATGGGCTGGATGCTCCTTGCCAATTGG
GACCTGATGGGTTACCTATGCCTGGCTGCTGGCAAAAG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG219509 representing NM_029838
 Red=Cloning site Green=Tags(s)

```


MLVKKLAGKGGGREGSGEDPRPLGQRCAGTMPSTALATLLSVVAVAFCFYLGVKTNDLQARIAALESAK
GTPSFHPLSDTVDELKAMVQEKVERLLAQKSYEHMAKIRTAREAPLECNCPAGPPGKRGKRGRREGSGPP
GQPGPQGPPGPKGDKGEQDQGPRMVFVKINHGFLSADQQLIKRRLIKGDQAGPPGPPGPPGPRGPPG
DTGKDGPRGMPGVPGEKPGEQQLMGPLGPPGQKGSIGAPGTPGMDGQKGEPEGSPGAAGQSGLPGPKGE
PGKEGEKGDAGENPKGDTGEKGDPSAAGIKGEPGESGRPGQKGEPLPGLPGLPGIKGEPGFIGPQG
EPGLPGLPGTKGDRGEAGPPGRGERGDPGAPGPKGQGESGARGPKGSKGDRGDKGDSGALGPRGPPGQK
GDPGATEIIDYNGNLHEALQRITTLTVTGGPPGPPGQGLQGPKGEQGSPIGVDGEQGLKGSKGMGDP
GVPGEKGLGLPGLPGANGVKGEKGTGLPGPQGPSIIIGPPGPPGPHGPPGPMGPHGLPGPKGASGLDGK
PGSRGADGPIGPHGAPGPKGERGEKGAMGEPGPRGPYGLPGKDGEPGLDGFPPGRGEKGDLEKGEKGF
GVKGEKGEPEGQPLDGLDAPCQLGPDGLPMPGCWQK
  
```


TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_029838

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029838.4](#), [NP_084114.2](#)

RefSeq Size: 2637 bp

RefSeq ORF: 2001 bp

Locus ID: 77018

UniProt ID: [Q99MQ5](#)

Cytogenetics: 3 G3

Gene Summary: Inhibits fibrillization of amyloid-beta peptide during the elongation phase. Has also been shown to assemble amyloid fibrils into protease-resistant aggregates. Binds heparin (By similarity).[UniProtKB/Swiss-Prot Function]