

Product datasheet for MG219625

Gml (NM_001177524) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gml (NM_001177524) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Gml
Synonyms: EG625599; HemT
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG219625 representing NM_001177524
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGATGCTCCCCTTCTTCTGTCAATTCTCATGGGTCTCCCATGGGTGGATACCTCTATCAACAATACCA
 GTGGATTGGTATTCGACAACACTACTGGTGGCTTGGACAATGCCATTGAGCCTAGATGGACTCCTCAGAT
 GAGATGCCATTCATGCCAGGAGCATAACACCTTCTATTGTCCACATATCCACTATTGTGACATGGATATC
 AGACGCTGTTTGACAGTTGCCATACGTGTGAACCTTCGGCTGCTCTATGTCTACAAGCACTGCACCAAGG
 ACTGCACATTTATTTATAGAAGGCATGTGCCTCCAGAACTCCCTAGGGTGCTTAAAGATGTCAAGAGCTT
 CTACTTTGTCTTGTGCTGTGGGAGTGTGCTTTGCAATGAGGGAGGACCCAGAAATATGGAGAGGGACTTG
 TTAGGTGAAACTTCCATAGAGGAAGAGGTAATTGCTAGAGCTGTGAGTCTGGGTGGGTCAACCTTCTCC
 TTTGTCTTGCCTTGATCCTCTCCAGCATCATACTGACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG219625 representing NM_001177524
 Red=Cloning site Green=Tags(s)


MMLPFFLSILMGLPWVDT SINNTSGLVFDNTTGGLDNAIEPRWTPQMRCHSCQEHNFTYCPHIHYCDMDI
 RRCLTVAIRVNLRLLYVYKHCTKDCTFIYRRHVPELPRVLKDVKSFYFVLCGSVVCNEGGPRNMRDL
 LGETSIEEEVIARAVSLGWVNL LCLALILSSIILT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001177524

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001177524.1](#), [NP_001170995.1](#)


RefSeq Size: 895 bp

RefSeq ORF: 531 bp

Locus ID: 625599

Cytogenetics: 15 D3

Product images:



Circular map for MG219625