

Product datasheet for MG219664

Gm4027 (NM_001177564) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm4027 (NM_001177564) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Gm4027
Synonyms: 100042776
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG219664 representing NM_001177564
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGCCAAAGAATAAAGGCAAAGGAGGCCAAAACAGGCGCAGAGGTAAAAATGAAATGAATCTGAGAAA
 GAGAGTTGGTGTAAAGAATATGGGCAAGAGTATGCTCAGGTGACCAAAATGCTGGGATGTGGATGGT
 GGAAGCAATGTGTTTGTGTTGTAAGGAGGCTGTGCCACATAAGAGGAAAGTTGAGAAAAAGGTTTG
 ATAAATACCTCAGACATTATATTGATTGGTCTACGAGACTATCAAGATAACAGGGCTGATATAATCTTAA
 AGTACAATCTAGATGAAGCAAGAAGTCTGAAGGCCTATGGAGAACTCCAGAACATGCCAAAATCAATGA
 AATGGACACATTTGGTGTGGGATGATGAAGAAATCGTGTGTTGATGATTGGAGAGGATGATGAAGAC
 ATTGATGACATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG219664 representing NM_001177564
 Red=Cloning site Green=Tags(s)

MPKNKGKGGKNRRRGKNENESEKRELVFKEYGQEYAQVTKMLGCGWLEAMCFDGVRRRLCHIRGKLRKKVW
 INTSDIILIGLRDYQDNRADIILKYNLDEARSLKAYGELPEHAKINEMDTFGAGDDEEIVFDDIGEDDED
 IDDI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001177564

ORF Size: 432 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001177564.1](#), [NP_001171035.1](#)


RefSeq Size: 998 bp

RefSeq ORF: 435 bp

Locus ID: 100042776

Cytogenetics: 12

Product images:



Circular map for MG219664