

Product datasheet for **MG219795**

Pth2 (NM_053256) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pth2 (NM_053256) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Pth2
Synonyms: Tifp; Tifp39; TIP; Tip39
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG219795 representing NM_053256
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGACCTGCCAGATGTCCAGGAGCCCCGAGAGCGGCTGTTGCTGCTTTTGCTGCTGCTACTGCTTG
TGCCCTGGGGCACTGGCCCTGCCTCAGGTGTTGCCCTGCCCTCGCTGGTGTGTTTCAGCCTCCGCGCCCC
GGGTCGTGCCTGGGCGGGCTTGGGTAGCCCCCTGTCTCGGCGCAGCCTGGCGCTAGCTGACGACGCGGCC
TTTCGGGAGCGCGCGCCTGCTGGCCGCCCTGGAGCGCCGCCGCTGGCTGGACTTTACATGCAGAAGC
TGTGCTACTGGACGCGCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG219795 representing NM_053256
Red=Cloning site Green=Tags(s)
METCQMSRSPRERLLLLLLLLLVPWGTGPASGVALPLAGVFSLRAPGRAWAGLSPLSRRSLALADDAA
FRERARLLAALERRRWLDSYMQKLLLLDAP

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053256.2</u> , <u>NP_444486.1</u>
RefSeq Size:	569 bp
RefSeq ORF:	303 bp
Locus ID:	114640
UniProt ID:	<u>Q91W27</u>
Cytogenetics:	7 B3
Gene Summary:	This gene encodes the precursor of a peptide hormone that shares sequence similarity with the parathyroid hormone. This gene is expressed in various regions of the brain where it plays a role in the release of pituitary hormones, anxiety and nociception. The encoded precursor protein is proteolytically processed to generate the biologically active neuropeptide. Mice lacking the encoded protein display increased fear and anxiety after exposure to stressful events, and decreased sensitivity to pain. [provided by RefSeq, Aug 2015]