


Product datasheet for **MG219999**

Gpr35 (NM_022320) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr35 (NM_022320) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gpr35
Synonyms:	KPG_007
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



Cloning sites used for ORF Shuttling:

CTATAGGGCGCGCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGCGCCAGATCT

EcoRI *BamHI* *KpnI* RBS Kozac Consensus *SgfI* *AscI*

CAAGCTTAACTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- --- ---

HindIII *NheI* *RsrII* *MluI* *NotI* *XhoI* GFP Tag

T R T R P L E M E S D - - -

--- --- GAA GAA AGA GTT TAA ACGGCCGCCCGGAGCT

- - - E E R V stop

PmeI *FseI*

ACCN: NM_022320


ORF Size: 921 bp



View online »

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022320.4 , NP_071715.3
RefSeq Size:	4233 bp
RefSeq ORF:	924 bp
Locus ID:	64095
UniProt ID:	Q9ES90
Cytogenetics:	1 D
Gene Summary:	Acts as a receptor for kynurenic acid, an intermediate in the tryptophan metabolic pathway. The activity of this receptor is mediated by G-proteins that elicit calcium mobilization and inositol phosphate production through G(qi/o) proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG219999