

## Product datasheet for **MG220115**

### **Rbbp9 (NM\_015754) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rbbp9 (NM\_015754) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Rbbp9  
**Synonyms:** Bog  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG220115 representing NM\_015754  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGTCCCCAACAAAGGCAGTGATTGTTCCCGGAACGGAGCGGGGATGTGGCCACCCACGGCTGGT  
 ACGGCTGGGTGAAAAAGGGGCTGGAGCAGATTCCTGGTTCCAGTGTGGCTAAAAACATGCCTGACCC  
 AATTACAGCGGAGAGAGCATCTGGCTGCCCTTCATGGAGACAGAGCTGCACTGTGACGAGAAGACCATC  
 ATCATAGGCCACAGTCCGGGGCCATCGCAGCCATGAGGTATGCAGAGACACATCAGGTATATGCTCTCG  
 TATTGGTGTCTGCATACACATCAGACTTGGGAGATGAAAATGAGCGTGCAAGTGGTACTTCAGCCGCC  
 CTGGCAGTGGGAGAAGAACAAGCCAAGTCCCTCACATTGTACAGTTTGGCTCTACTGATGACCCCTTC  
 CTTCCCTGGAAGGAACAACAAGAAGTAGCAGATAGGCTGGACGCCAAATTGTACAAATTCAGTACCGTG  
 GTCACCTTCAGAACACAGAGTTCATGAAGTATTAGTGTGGTGAAGTCTATGTTGAAAGGACCAGAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG220115 representing NM\_015754  
 Red=Cloning site Green=Tags(s)

MASPNKAVIVPGNGGGDVATHGWYGVVKKGLEQIPGFQCLAKNMPDPITARESIWLPFMETELHCDEKTI  
 IIGHSSGAIAMRYAETHQVYALVLSAYTSDLGDENERASGYFSRPWQWEKNKANCPHIVQFGSTDDPF  
 LPWKEQQEVADRDLAKLYKFTDRGHFQNTFHELISVVKSMKLGPE


**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online >](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_015754

**ORF Size:** 558 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_015754.2</a></u> , <u><a href="#">NP_056569.2</a></u>
<b>RefSeq Size:</b>	2183 bp
<b>RefSeq ORF:</b>	561 bp
<b>Locus ID:</b>	26450
<b>UniProt ID:</b>	<u><a href="#">O88851</a></u>
<b>Cytogenetics:</b>	2 G1
<b>Gene Summary:</b>	May play a role in the transformation process due to its capacity to confer resistance to the growth-inhibitory effects of TGF-beta1 through interaction with retinoblastoma and the subsequent displacement of E2F-1.[UniProtKB/Swiss-Prot Function]