

Product datasheet for **MG220190**

Npas4 (NM_153553) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Npas4 (NM_153553) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Npas4
Synonyms:	LE-PAS; Nxf
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG220190 representing NM_153553
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTACCGATCCACCAAGGGCGCCTCCAAGGGCGCCGACCAGATCAACGCCGAGATTCGGAACCTCA
 AGGAGCTGCTGCCGTTGGCTGAAGCGGACAAGGTCGGCTGTCTACCTGCACATCATGAGTCTTGCCCTG
 CATCTACACTCGCAAGGGTGTCTCTTTGCTGGAGGCACTCCTTTGGCTGGCCCCACCGGCTTCTCTCT
 GCTCAAGAGCTTGAAGACATTGTGGCAGCACTACCTGGATTTCTCCTTGTATTACAGCTGAGGGGAAGT
 TGCTATACCTGTCGGAGAGTGTGAGCGAGCATCTGGCCACTCTATGGTGGACCTGGTTGCCAGGGCGA
 CAGTATCTACGATATCATTGACCCTGCTGACCATCTCACTGTGCGCCAGCAGCTCACCATGCCCTCTGCT
 CTGGATGCTGATCGCTTTTCCGTTGTCGATTCAACACCTCCAAGTCCCTCCGGCCAGAGTTCAGGAA
 ACAAAGTGGTGTATTTCGAGGTCGATTCCATGCTCACCCACTGGGGCCTACTGGGCAGGAAACCTGT
 GTTCACCGCTTTCTGCGCCCACTGGAGCCAAGACCCCGCCCTGGCCCCGGCCCTGGCCCTGGCCCTGGT
 CCTGCTTCTCTTCTCCTGGCCATGTTCCAGAGCCGGCATGCTAAGGACCTAGCCCTACTGGACGTTTCTG
 AAAGTGTCTTAATCTACCTGGGCTTTGAGCGCAGCGAACTGCTCTGTAATCATGGTATGGACTGCTACA
 CCCCAGGACCTGGCCCAAGCTTCTTCTCAACTACCGCCTGTTGGCTGAAAGTGGAGATATTCAGGCT
 GAAATGGTGGTGAGACTTCAAGCCAAGCATGGAGGCTGGACATGGATTTACTGCATGCTATACTCAGAAG
 GTCCAGAAGGCCCTATTACTGCCAATAACTACCCTATCAGTGACACGGAAGCCTGGAGCCTCCGCCAGCA
 GCTAAACTCTGAAGACACCCAGGAGCCTATGTCCTAGGAACCCAGCTGTGCTACCCCTATTCTCTGAG
 AATGTCTTCTCCAGGAGCAATGCTCTAATCCACTCTTACACCATCCCTGGGGACTCCTAGAAGTGCCA
 GCTTCCCAGGGCTCCTGAACTAGGTGTGATCTCAACACCAGAAGAGCTTCCCAACCTCCAAAGAGAT
 GGACTTCAGTTACCTGCCATTCCCTGCTAGGCTGAGCCTTCCCTCCAAGCAGACCTGAGCAAGGATTTG
 GTGTGACTCCACCTTACACACCCACAGCCAGGAGGCTGTGCCTTCTCTTACGCTCCATGAACCTT
 TCCAGACTCACTTGGCCCTCCGTCCAGCTCTCTCCAAGAACAGCTGACACCAAGTACAGTGACTTTCTC
 TGAACAGTTGACACCCAGCAGTGTACCTTCCAGACCCACTAACAGTTCACACAGGACAGTTGACA
 GAAAGCTCAGCCAGAAGCTTTGAAGACCAGTTGACTCCATGCACCTTCTCCTTCCCTGACCAGCTACTTC
 CCAGCACTGCCACATTCCAGAGCCTCTGGGCAGCCCCGCCATGAGCAGCTGACTCCTCCAGCACAGC
 ATTCCAGGCTCATCTGAACAGCCCCAGCCAAACCTTCCAGAGCAACTGAGCCCCAATCCTACCAAGACT
 TACTTCGCCCAGGAGGGATGCAGTTTTCTCTATGAGAAGTTGCCCAAGTCTAGCAGCCCTGGTAATG
 GGGACTGTACACTCCTGGCCCTAGCTCAGCTCCGGGGCCCTCTCTGTGGATGTCCCCCTGGTGCCCGA
 AGGCTGCTCACACCTGAGGCTCTCCAGTCAAGCAAAGTTTCTTCCACTACACAGAGAAAGAGCAAAAT
 GAGATAGATCGTCTCATTAGCAGATCAGCCAGTTGGCTCAGGGCGTGGACAGGCCCTTCTCAGCTGAGG
 CTGGCACTGGGGGGCTGGAGCCACTTGGAGGGCTGGAGCCCTGAACCTAACCTGTCCCTGTCAGGGGC
 TGGACCCCTGTGCTTAGCCTGGATCTTAAACCTGGAATGCCAGGAGCTGGACTTCTGGTTGACCTT
 GATAATTTATTCCTGGAAGAGACGCCAGTGAAGACATCTTATGGATCTTTCTACTCCAGACCCCAATG
 GGAATGGGGTTCAGGGGATCCTGAGGCAGAGGTCAGGAGGGACCCTGTCACCTTGAACAACCTGTC
 CCCAGAAGATCACAGCTTCTGGAGGACTTGGCCACCTATGAAACCGCCTTTGAGACAGGTGTCTCAACA
 TTCCCCTACGAAGGGTTTGTGATGAGTTGCATCAACTCCAGAGCCAAGTTCAAGACAGCTTCCATGAAG
 ATGGAAGTGGAGGGGAACCAACGTTT

ACCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >MG220190 representing NM_153553
 Red=Cloning site Green=Tags(s)

```

MYRSTKGASKARRDQINAEIRNLKELLPLAEADKVRLSYLHIMSLACIYTRKGVFFAGGTPLAGPTGLLS
AQELEDIVAALPGFLLVFTAEGKLLYLSSESVSEHLGHSMVDLVAQGDSIYDIIDPADHLTVRQQLTMP
LDADRLFRCRFNTSKSLRRQSSGNKLVIRGRFHAHPPGAYWAGNPVFTAFCAPLEPRPRPGPGPGPG
PASLFLAMFQSRHAKDLALLDVSESVLIYLGFERSELLCKSWYGLLHPEDLAQASSQHYRLLAESGDIQA
EMVVRQLQAKHGGWTWIYCMYSEGPEGPITANNYPI SDTEAWSLRQQLNSEDQAAYVVLGTPAVLPSFSE
NVFSQEQC SNPLFTP SLGT PRSASF PRAPELGVISTPEELPQPSKELDFSYLPF PARPEPSLQADLSKDL
VCTPPYTPHQPGGCAFLFSLHEPFQTHLPPSSSLQEQLTPSTVTFSEQLTPSSATFPDPLTSSLQGQLT
ESSARSFEDQLTPCTSSFPDQLLPSTATFPEPLGSPAHEQLTPPSTAFQAHLNPSQTFPEQLSPNPTKT
YFAQEGCSFLYEKLPPSPSPGNGDCTLLALALQLRGPLSVDVPLVPEGLLTPEASPVKQSFHYTEKEQN
EIDRLIQQISQLAQGVDRPFSAEAGTGGLLEPLGGLEPLNPNL SL SGAGPPVLSLDLKPWKQCQLDFL
DNL FLEETPVEDI FMDLSTDPNGEWGSDPEAEVPGGTLSPCNNLSPEDHSFLEDLATYETAFETGVST
FPYEGFADELHQLQSQVQDSFHEDGSGGEPTF
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfi-MluI

Plasmid Map:


ACCN: NM_153553

ORF Size: 2406 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_153553.4](#), [NP_705781.1](#)

RefSeq Size: 3292 bp

RefSeq ORF: 2409 bp

Locus ID: 225872

UniProt ID: [Q8BGD7](#)

Cytogenetics: 19 A

Gene Summary:

Transcription factor expressed in neurons of the brain that regulates the excitatory-inhibitory balance within neural circuits and is required for contextual memory in the hippocampus (PubMed:18815592, PubMed:22194569, PubMed:23029555, PubMed:24201284, PubMed:24855953). Plays a key role in the structural and functional plasticity of neurons (PubMed:23172225). Acts as an early-response transcription factor in both excitatory and inhibitory neurons, where it induces distinct but overlapping sets of late-response genes in these two types of neurons, allowing the synapses that form on inhibitory and excitatory neurons to be modified by neuronal activity in a manner specific to their function within a circuit, thereby facilitating appropriate circuit responses to sensory experience (PubMed:24201284, PubMed:24855953). In excitatory neurons, activates transcription of BDNF, which in turn controls the number of GABA-releasing synapses that form on excitatory neurons, thereby promoting an increased number of inhibitory synapses on excitatory neurons (PubMed:18815592, PubMed:22194569, PubMed:24201284). In inhibitory neurons, regulates a distinct set of target genes that serve to increase excitatory input onto somatostatin neurons, probably resulting in enhanced feedback inhibition within cortical circuits (PubMed:24855953). The excitatory and inhibitory balance in neurons affects a number of processes, such as short-term and long-term memory, acquisition of experience, fear memory, response to stress and social behavior (PubMed:18815592, PubMed:22194569, PubMed:23029555, PubMed:24201284, PubMed:27238022). Acts as a regulator of dendritic spine development in olfactory bulb granule cells in a sensory-experience-dependent manner by regulating expression of MDM2 (PubMed:25088421). Efficient DNA binding requires dimerization with another bHLH protein, such as ARNT, ARNT2 or BMAL1 (PubMed:14701734, PubMed:15363889, PubMed:19284974). Can activate the CME (CNS midline enhancer) element (PubMed:14701734, PubMed:15363889).[UniProtKB/Swiss-Prot Function]