

Product datasheet for **MG220191**

Mrpl24 (NM_026591) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrpl24 (NM_026591) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mrpl24
Synonyms:	2010005E08Rik; 2810470K06Rik; 6720473G22Rik; AA407670
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG220191 representing NM_026591 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGGCTCTCAGCCCTGCTGGCCTTGGCATCCAAAGCCACTTCCTCCCCCTTCTACCGCTATGGAATGA
GCCGTCCAGGTTCCATAGCAGACAAAAGGAAGAATCCTCCATGGAGCAGGCGGCCAGTAGTGGTAGA
GCCATCTCTGATGAAGACTGGCACCTATTCTGTGGAGACATGGTAGAAATCTTAGAAGGCAAGGATGCT
GGGAAGCAAGGCAAAGTAGTCCAAGTCGTCGGCAGCGAAACTGGGTTGCTCTGGAGGGTTGAACACGC
ATTACCGATACATTGGCAGAACCAAGGACCACCGAGGGACCATGATCGCTAGTGAAGCCCCCTTGCTTCA
TCACCAGGTCAAGCTAGTGGATCCTGTGGACAGGAAACCCACTGAGATCCAGTGGAGATTTACTGAGGCA
GGAGAACGTGTTGAGTCTCTACAAGATCTGGGAGAATTATCCCCAAACCTGAATTTCTAGAGCAGATG
GCATTGTCCTGAAACATGGACTGATGGCCCCAAGGACACATCAGTGGAAAGATGCTCTTGAAGAAGCTTA
TGTGCCCCGACTAAAGACATTAGAAGAAGATGTGATGGAGGCCATGGGGATCCAGGAGACTCGAAGATTC
AAGAAAGTTTATTGGTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG220191 representing NM_026591
Red=Cloning site Green=Tags(s)

MRLSALLALASKATSSPFYRYGMSRPGSIADKRKNPPWSRRRPVVVEPISDEDWHLFCGDMVEILEGKDA
 GKQGVVQVVRQRNWWVLEGLNTHYRYIGRTKDRGTMIASEAPLLHHQVKLVDPVDRKPTEIQWRFTEA
 GERVVRVTRSGRIIPKPEFPRADGIVPETWTDGPKDTSVEDALERTYVPRLKTLEEDVMEAMGIQETRRF
 KKVVYWY

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_026591

ORF Size: 648 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026591.3](#), [NP_080867.2](#)

RefSeq Size: 949 bp


RefSeq ORF: 651 bp

Locus ID: 67707

UniProt ID: [Q9CQ06](#)

Cytogenetics: 3 F1

Product images:



Circular map for MG220191