

## Product datasheet for **MG220235**

### Lrrc69 (NM\_028499) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc69 (NM_028499) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Lrrc69
Synonyms:	1700034K16Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG220235 representing NM_028499 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTGAAAGATTATTGGTAACAGCACTGAAGGGTGGCAAGAACAATAAGATCCTGACTTTGAATGGCA  
AGAGGATTACCAAGATGCCCTCGACACTGGAAAACTACCTAATCTGAAGACTCTAGACCTTCAAATAA  
TTCCATCTCAAAGTATGCCAGAGCTAAGGACCTTGACCCAGTTGACGTTACTAACTGGGAAACAAC  
CATTTACAAGAAGTTCCTGAAGAGATAAAATATCTTACCTCCCTGAAGAATCTTCATTTATTTGAAACA  
GGATTTGCAGAATTGCTCCTGGAGTTTTAATGGCCTACATCGCTTGATCATGCTTAATCTGAATGACAA  
CAGACTTACGAGTCTTCCCAAGAAATTGGCAGGTTAAGGAGTCTTACCTACTTAAGCCTCAACCGTAAC  
AATCTCACCGTGATTTCCAAAGAGCTTTGTTCCCTTGAGCATCTTCTGAACTGCATCTAACTACAATC  
AGATCGTGTACATCCCTGAGGAGATCAAGTTCTGAAAGAACCTTCAGCAGCTCTTCTGGTCAGGAACAA  
CATTGAAGAGTTGCCAGAGGAAATTTGTCATCTTGAAAACTAAGAGTCTGGACATAGCTGGAAATGTA  
ATCCAGATATTTCCGGCAGGCTTTCAGAATCTGAGGCTGACAGAGTCTACTGTGAGGAAAACCCGCTCT  
TCCTGAAGAGGCCATTTTTGCTGTACAACCTAAAGACCTCTGGACTCTACGGGAAATAGCAGCAAGATT  
TGTCCTCTCTGTTAGAAGAAAACGACCCTTTGATTATGAATGCATAGAGAGTTACCCAGAGGTCAAG  
GACAACTCTCAAAGCAAAAAGTGTCAATATGCAGGAAACCTTTTTAACTGAGTGGCTGGAGTGCG  
TTTACTTTGTGGCTCCATCGAAGAAGTGAAGATAAGCAGAAACCTAAAGCTGATTCTGTCAGACATC  
AGTCTGTTCTATCAATGCTTTGACCAGCGTGATCCTGATGTCTTTGGAATTGCTCAAGAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


**Protein Sequence:** >MG220235 representing NM\_028499  
 Red=Cloning site Green=Tags(s)

MAERLLVTALKGGKNTKILTLNGKRITKMPSTLEKLPNLKTLDLQNNISISKVCPPELRTLQTLTLLNLGNN  
 HLQEVPEEIKYL TSLKNLHLFGNRIACRIAPGVFNGLHRLIMLNLNDNRL TSLPQEIGRLRSLTYLSLNRN  
 NLTVIPKELCSLEHLSELHLNQNQIVYIPEEIKFLKNLQQLFLVRNNIEELPEEICHLEKLRVLDIAGNV  
 IQIFPAGFQNLRLTEFYCEGNPLFLKRPFFAVQPKDLWTLREIAARFVLSLLEENDPLIMNVIESYPEVK  
 DKLSKAKKCSICRKPFLTEWLECVYFVAPSKNWKISRNLKLIIPVQTSVCSYQCFDQRDPDVFVIAQE


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_028499

**ORF Size:** 1041 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_028499.2</a> , <a href="#">NP_082775.1</a>
<b>RefSeq Size:</b>	1210 bp
<b>RefSeq ORF:</b>	1044 bp
<b>Locus ID:</b>	73314
<b>UniProt ID:</b>	<a href="#">Q9D9Q0</a>
<b>Cytogenetics:</b>	4 A1