

Product datasheet for **MG220678**

Cd300c (NM_199225) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cd300c (NM_199225) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Cd300c
Synonyms: Clm6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG220678 representing NM_199225
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATCCAGAGTAATAAGATTGTGGCTGCCTTCAGCTGTGCTACTCTCCTTGGTTCCAGGCCATTTCC
CTGTGCGTGGTCCCAGCACTGTACAGGCACTGTGGGGAGTCACTCAGTGTGCCTGTCAGTATGAGAA
GAAATTGAAAATAAAAAAAAACTGGTGCAAATGGAATCTAATGTACTGTGTAAGATATTGTCAAG
ACCAGTGCCTCAGAAGAAGCTAGGAATGGCCGAGTGTCCATCAGGGACCATCCAGACAACCTCACCTTCA
CAGTGACCTTGGAGAACCTCACCTGGAGGATGCAGGCACCTACATGTGTATGGTGGATATAGGATTCTT
TTATGATGCCTACTTACAGATCGATAAGTCCTCAAGTTGAGGTGTTCGTGGTCCCAGGTAAGCCCTCT
GTCAAAGTTTCACGCCCTGGGAACGGTATCAACATTCTGGCTTCCCCACATCCTCAGCGGTGCACACTC
AGCCCAACGTGACCACAGACGACAAATTCCTGCTCCCAGCCCTGAGCTTCGGTCTCTTCTGAGCAGCCC
CCACTTCTGGATCCTGGTGTCTCTGAAGTTGCCCTGTTCTGAGCATGCTTGGTGCCTTGTCTGGGTG
AACAGGCTCAGAGGTGCTCTGGGGGAAGCAGCGCTTGGCCCTGTTATGAGAACCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG220678 representing NM_199225
Red=Cloning site Green=Tags(s)

MNPRVIRLWLPSAVLLSLVPGHFPVRGPSTVTGTVGESLSVSCQYEKLLKTKKKIWKWKSINVLCCKDIVK
 TSASEEARNGRVSIRDHPDNLTFVTLENLTLEDAGTYMCMVDIGFFYDAYLQIDKSFKVEVFVVPKPP
 VKGSRPNGINILASPTSSAVHTQPNVTTDDTIPAPSELRSLSSPHFWILVSLKLPFLSMLGALLWV
 NRPQRCSGGSSAWPCYENQ

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja1703_h05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_199225

ORF Size: 687 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199225.1](#), [NP_954695.1](#)

RefSeq Size: 690 bp


RefSeq ORF: 690 bp

Locus ID: 387565

UniProt ID: [A2A7V7](#)

Cytogenetics: 11 E2

Product images:



Circular map for MG220678