

Product datasheet for **MG220869**

Fkbp10 (NM_001163481) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|-----------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | Fkbp10 (NM_001163481) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Fkbp10 |
| Synonyms: | AI325255; FKBP-10; FKBP-65; Fkbp-rs1; Fkbp1-rs; Fkbp6; FKBP65; Fkbprp |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>MG220869 representing NM_001163481
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTTCTTGTGGGGTCTCCAGCCACACCTCCATCGGCTCCGCATACTGCCGTTGCTGTGCTTCTAC
 AGACCTTGAGAGGGGACTGGGCCGTGCCAGCCCGGCCGAGCCCTTGGAAGATGTGGTCATCGAGAG
 ATACCACATCCCTCGGGCCTGTCCCGAGAAGTGCAGATGGGGATTTTGTGCGTTACCACTACAATGGC
 ACTTTGAAAGACGGGAAAAAGTTTGACTCCAGCTATGACCGTAGCACCTGGTGGCCATCGTTGTGGCG
 TAGGCCGCTCATCACCGGCATGGACGGGGTCTCATGGGCATGTGTGTCAACGAGCGCCGCCCTCAT
 TGTGCCTCCACCTGGGCTACGGCAGCATCGGTGTGGCGGGCTCATCCCCCTGATGCCACCCTCTAT
 TTTGACGTGGTCTGCTGGACGTGTGGAACAAAGCAGACCGGTGCAGTCAACTATCCTCTGCGCCCTC
 CCTACTGCCCCGAATGGTGCAGAACAGTGACTTTGTGCGCTATCACTACAATGGCACTCTGCTGGATGG
 CACTGCCTTTGACAACAGCTACAGTAGGGGAGGCACTTATGACACCTACATCGGCTCTGGTTGGCTGATC
 AAAGGCATGGACAGGGGCTGCTGGGCATGTGCCCTGGAGAGAAAAGGAAGATCATTATCCCTCCCTTCC
 TGGCTTATGGGGAGAAAGGCTATGGGACTGTGATACCCCGCAGGCCCTCCCTGGTCTTCTATGTCCTGCT
 GCTGGATGTCCACAACCGAAGGACACGGTCCAGCTGGAGACGCTGGAGCTGCCCCAGGGCTGTGTCGG
 CGAGCTGTGGCGGGGACTTCATGCGTTACCACTACAATGGCTCTTTGATGGATGGTACCCTCTTTGATT
 CCAGCTACTCCGAAACACACCTACAATACCTATGTCGGGCAGGGTTACATCATCCCTGGGATGGACCA
 GGGGCTGCAAGGCGCATGCATAGGGGAGCGAAGGAGGATTAAGTGTGCCCCCTCACCTTGCCTACGGGGAG
 AATGGGACAGGAGACAAGATCCCTGGCTCGGCTGTGCTCATCTTCGATGTGCACGTCATCGACTCCACA
 ACCCTTCGGACCCTGTGAAATCAAGACGCTGTCCCGCCCTCTGAGAAGTGAATGAGACATCCAAAAT
 CGGGGACTTCATTCGCTACCACTACAACCTGTTCTCTGCTGGACGGCACCAGGCTCTTCTCGTCCACGAC
 TATGAGGCCCTCAAGAGATCACCTCGGAGCCAAACAAAGTGATCGAAGGTCTGGACAGGGGCTGCAGG
 GCATGTGTGTAGGAGAGAGGAGGAGCAGCTCATTGTCCCCACACCTGGCCACGGGGAGAATGGAGCCCG
 GGGTGTCCCTGGCAGTGTGTGCTATTTGAGGTGGAGCTGGTATCCCGAGAGGATGGCTGCCACACA
 GGCTACCTGTTTGTGTTGACCAGGATCCTTCCACTAGCCTTTTGAAGACATGGATCTCAATAAAGATG
 GAGAGGTGCCCCAGAAGAGTTCTTCTTCCATCAAGGCTCAAGTGAATGAAGGCAAAGGACGCCTCAT
 GCCTGGGCAAGACCCAGACAAAACCATCAGTGACATGTTTCAAGACCAGGACCGGAACCAGGATGGCAAG
 ATCACAGCTGAGGAGCTCAAGCTGAAGTCAAGTGAAGGACCAGGAGCGGGTCCATGAGGAGCTC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG220869 representing NM_001163481
 Red=Cloning site Green=Tags(s)


MFLVGSSTLHRLRILPLLLLLQTLERGLGRASPAGAPLEDVVIERYHIPRACPREVQMGDFVRYHYNG
 TFEDGKKFDSSYDRSTLVAIVVGVGRLITGMDRGLMGMVCNERRRILVPPHLGYGSIGVAGLIPPDATLY
 FDVLLLDVWVKADTVQSTILLRPPYCPRMVQNSDFVRYHYNGTLLDGTAFDNSYSRGGTYDYIYSGWLI
 KGMDQGLLGMCPGEKRKIIIPFFLAYGEKGYGTVIPPQASLVFYVLLLDVHNPKDTVQLETLELPQGCVR
 RAVAGDFMRYHYNGSLMDGTLFDSSYSRNHTYNTYVQGYIIPGMDQGLGACIGERRRITVPPHAYGE
 NGTGDKIPGSAVLIFDVHVIDFHNPSDPVEIKTLRPPENCNETSKIGDFIRYHYNCSLLDGTRLFSSHD
 YEAPQEITLGANKVIEGLDRGLQGMCVGERRQLVPPHLAHGENGARGVPGSAVLLFEVELVSREDGLPT
 GYLFVWYQDPSTSLFEDMDLNKDGVEPPEEFSSF IKAQVNEGKGRMPGQDPDKTISDMFQNDQDRNQGK
 ITAEELKLKSDQERVHEEL

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001163481

ORF Size: 1740 bp

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001163481.1 , NP_001156953.1 |
| RefSeq Size: | 2616 bp |
| RefSeq ORF: | 1743 bp |
| Locus ID: | 14230 |
| Cytogenetics: | 11 63.47 cM |
| Gene Summary: | PPlases accelerate the folding of proteins during protein synthesis.[UniProtKB/Swiss-Prot Function] |