

Product datasheet for **MG220950**

Yipf4 (NM_026417) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Yipf4 (NM_026417) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Yipf4
Synonyms: 2310034L04Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG220950 representing NM_026417
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGCCTCCGGGCCGCCCGGCCTATGCGCCTGCCAACGGGGACTTCACCTTTGTGTCTTCGGCGG
 ACGCGGAAGATCTCAGCGGTTCTATAGCAGCTCCAGATGTCAAATTAACCTCGGTGTAAGCGGAGATTT
 CATCAAAGAATCTACAGCTACTACATTTCTGAGACAAAGAGGCTATGGCTGGCTTTTGAAGTTGAAGAT
 GAGGATCCTGAAGATAACAAGCCACTTTTGAAGAATTGGACATCGATCTAAAGGACATTTACTACAAGA
 TCCGCTGTGTTCTGATGCCAATGCCATCGCTGGGTTTTAATAGACAGGTGGTGGAGACAATCCTGACTT
 CTGGGGTCCACTGGCTGTTGTTCTTTCTTTCCATGATATCATTATATGGACAGTTTAGGGTGGTCTCA
 TGGATTATAACCATTTGGATATTTGGTTCACTGACAATTTTCTTACTGGCTAGAGTTCTTGGTGGAGAAG
 TTGCATATGGCCAAGTCTTGGAGTTATAGGATATTCATTACTCCTCATTGTAATAGCCCCTATACT
 TTTGGTGGTTGGATCTTTTGAATGGTATCTACACTTATAAACTGTTTGGTGTGTTCTGGGCTGCCTAC
 AGTGCTGCTTCATTACTCGTAGGTGAAGAATCAAGACCAAAAAGCCCCTCCTCATTACCCCATCTTTT
 TATTGTACATTTATTTTTGTCTTATATACTGGATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG220950 representing NM_026417
Red=Cloning site Green=Tags(s)


MQPPGPPPAYAPANGDFTFVSSADAEDLSGSIAAPDVKLNLVSGDFIKESTATTFLRQRGYGWLLEVED
 EDPEDNKPLLEELDIDLKDIYYKIRCVLMPMPSLGFNRQVVRDNPDFWGPLAVVLFSSMISLYGQFRVVS
 WIITIWIFGSLTIFLLARVLGGEVAYGQVLGVIGYSLLPLIVIAPIILLVGSFEMVSTLIKLFVFWAAAY
 SAASLLVGEEFKKKPLLIYPIFLLYIYFLSLYTG

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_026417

ORF Size: 738 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026417.4 , NP_080693.2
RefSeq Size:	2112 bp
RefSeq ORF:	741 bp
Locus ID:	67864
UniProt ID:	Q8C407
Cytogenetics:	17 E2
Gene Summary:	Involved in the maintenance of the Golgi structure.[UniProtKB/Swiss-Prot Function]