

Product datasheet for **MG221082**

Dbn1 (NM_001177372) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dbn1 (NM_001177372) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Dbn1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG221082 representing NM_001177372
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCCGGCGTCAGCTTCAGCGGCCACCGCTGGAGCTGTGGCGGCGTACGAGGAGGTGATCCGGGAGG
AGAGCGCAGCCGACTGGGCTCTGTACACATACGAGGATGGCTCAGATGACCTCAAGCTTGCAGCGTCAGG
AGAAGGGGGCTTGCAGGAGCTTCCGGCCACTTCGAGAACCAGAAAGTATGATGATGGTTTCTGCAGCGTC
AAGGACTCCCAAGCTGCCCTGCCAAAATATGTGCTCATCAACTGGGTTGGTGAGGATGTGCCTGATGCC
GAAAATGTGCTTGCAGCAGTGTGTGGCCAAGGTGGTGAATTCTCCAGGGTGTGATGTCATTGTGAA
TGCCAGCAGTGTGAAGACATCGATGCTGGTGCCATTGGGCAGCGGCTCTCCAATGGACTGGCACGGCTC
TCCAGCCCAGTATTGCACCGCTGCGCCTTCGGGAGGATGAAAATGCTGAACCGTGGGTACCACCTACC
AGAAGACGGATGCAGCAGTGGAGATGAAGCGGATTAACCGTGAGCAGTTTTGGGAGCAGGCCAAGAAGGA
GGAAGAGCTGCGGAAGGAGGAGGAGCGGAAGAAGGCTCTGGACGCCAGGCTCAGGTTTGAACAGGAACGG
ATGGAGCAGGAGCGGCAGGAGCAGGAAGAAGCTGAGCGGCGCTACCGGGAGCGGGAGCAGCAGATTGAGG
AGCACAGGAGGAAACAGCAGAGTCTGGAAGCTGAAGAAGCCAAGAGGAGGTTAAAGGAGCAGTCTATCTT
TGGTGACCAGCGGGATGAAGAGGAAGAGTCCCAGATGAAGAAGTCGGAGTCAGAGGTGGAGGAGCGGCT
GCCATCATTGCCAGCGGCTGATAACCCACGGGAGTTCTTCAGACAGCAGGAACGAGTGGCATCGGCT
CTGGTGGCAGCTGTGACGCGCTGCGCCTGCACCTTCAACCACCGACCAGGCAGCCACCTGGACAGCCA
CCGGAGGATGGCACCCTCTATTCCCACCGGAGCCCATCTGATTCCAGCACAGCCTTACCCCATC
GCTGAGCAGATCGAGAGGGCCCTGGATGAGGTCACATCCTCGCAGCCTCCACCTCCACCTCCACCTC
CACCAACTCAAGAGGCCAGGAGACTACCCCAAGCCTGGATGAAGAGCTCAGCAAGGAGGCCAAAGTAAC
AGCAGCTCCTGAGGTCTGGGCTGGCTGTGCGGCAGAGCCCCCTCAGGCACAGGAACCTCCCCTGTTGCAA
AGCAGCCCCCTGGAGGACTCGATGTGCACAGAATCTCCAGAGCAGGCTGCCCTGGCTGCCCTGCGGAGC
CTGCTGCCTCTGTACCTCAGTAGCTGATGTCCATGCAGCTGACACCATTGAGACCACCTGCCACTAC
TGACACCCTATTGCCAACACGTACCCCTGCCGCTGCCAGCCTCATTGATCTATGGCCTGGCAACGGG
GAAGAGGCTCAACACTCAGGCTGAACCCAGGGTGCCACACCACCTCAGGTGCTGAGGCTCCCTGG
CAGAGGTGCCCTGCTGAATGAGGCGCTCAGGAGCCGCTGCCGCGGTAGGCGAAGGCTGTGCTAACCT
TCTTAATTTTATGAGCTGCCAGAACCTCCAGCCACCTTCTGTGACCCAGAGGAGGAAGTAGGAGAAACG
CTGGCTGCCTCCAGGCTCTAATATGCCCTCAGCTCTAGAGGAGGTAGATCAGGTGCTGGAGCAGGAGC
TGGAGCCAGAACCTCACCTGCTGACCAATGGAGAGACCTCAAAAGGAGGGGACCCAGCAGGCCAGCGA
AGGATACTTCAGTCAGTCACAGGAGGAAGATTCGCCAATCAGAAGAGCCATGTGCAAAGGTTCCGCT
CCTGTATTTTACAACAAGCCTCCAGAAATCGACATCACCTGCTGGGATGCAGACCAGTTCTGAAGAGG
AAGAGGGCTTCGAGGGTGGTGAT

ACCGGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG221082 representing NM_001177372
Red=Cloning site Green=Tags(s)

```
MAGVSFSGHRLELLAAYEEVIREESAADWALTYEDGSDDLKLAASGEGGLQELSGHFENQKVMYGFCSV
KDSQAALPKYVLINWVGEDVPDARKCACASHVAKVAEFQGVVDIVNASSVEDIDAGAIGQRLSNGLARL
SSPVLHRLRLREDENAEPVGTYYQKTDAAVEMKRINREQFWEQAKKEEELRKEEERKKALDARLRFEQER
MEQERQEERERRREREREQQIEEHRRKQQSLEAEEAKRRLKEQSIFGDQRDEEEESQMKKSESEVEEAA
AIIAQRPDNPREFFRQQERVASASGGSCDAPAPFNFHRPGSHLDSHRRMAPTPIPTRSPSDSSTASTPI
AEQIERALDEVTSSQPPPPPPPPPTQEAQETTPSLDEELSKAKVTAAPVWAGCAAEPQAQEPPLLQ
SSPLEDSMCTESPEQAALAAPAEPAAASVTSVADVHAADTIETTTATDITIANNVTPAAASLIDLWPGNG
EEASTLQAEPRVTPPSGAEASLAEVPLLNEAAQEPLPPVGECCANLLNFDELPEPPATFCDEPEEVGET
LAASQVLTPMSALEEVDQVLEQELEPEPHLLTNGETTQKEGTQQASEGYFSQSQEEFFAQSEEPKAVPP
PVFYNKPEIDITCWDADPVPEEEEGFEGGD
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001177372

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177372.1](#), [NP_001170843.1](#)

RefSeq Size: 2909 bp

RefSeq ORF: 1986 bp


Locus ID: 56320

UniProt ID: [Q9QXS6](#)

Cytogenetics: 13 30.06 cM

Gene Summary: Actin cytoskeleton-organizing protein that plays a role in the formation of cell projections (By similarity). Required for actin polymerization at immunological synapses (IS) and for the recruitment of the chemokine receptor CXCR4 to IS (By similarity). Plays a role in dendritic spine morphogenesis and organization, including the localization of the dopamine receptor DRD1 to the dendritic spines (PubMed:25865831). Involved in memory-related synaptic plasticity in the hippocampus (PubMed:25865831).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG221082