

Product datasheet for **MG221226**

Ldlrad1 (NM_001081272) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ldlrad1 (NM_001081272) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ldlrad1
Synonyms:	OTTMUSG00000008594
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG221226 representing NM_001081272 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGTGGGCTGTGAGAGAGCAGAATTTTGCCAGCCTGCTTCTGGTTTCTTGCCAGGCAGATGCCG
AGAGTAACATCACCCTAGAAAGAAAGCCGCCCTGGAGATGAAGGCTGTGGTCCACTCTGCTGCTCAG
CCGTGGAGCCTGCCTCTCGGTCGGCCTGCTGCTCCTGCTGGCCATGCTCGCTGCCTTGCTGTGGCC
ACCATCCTCGGACATCCACCCCGCACACCAGAGACTCACTCCTGCGTGACTAGCCAACAGGACAGGGT
TCTTGTGCCACGACCGAAGGCACTGCATCCAGCCATGCTGTGTGCGATGGCATCCGCACCTGTCTCA
TGGTGAGGACGAGGCTGAGAGCCTGTGCCGAGATGTGCCCCAGAGCCTCCCAGCTTCTCTGACCACC
TGTGGAGACCCAGACTCCTGGATCTACTCAGACCAGAAATGTGATGGCGTCAACAACCTGTGGAGACTGT
CAGATGAGCTGAGCCAGTGACTACGTGCCCGCCTGTGGCCCTGGGTGGTGGCGCTGCTCGCCAACCTGT
GTTCAAGTACTGCAGCTGTGCCGAGGGATCTGTGCCGGACAGCGTGCAGCACTGCTCTGACTGGTCA
GATGAGTACTCCTGTCTGGACCC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG221226 representing NM_001081272
 Red=Cloning site Green=Tags(s)

MSVGCERAEFCACFWFLGQADAESNITRRKARPGDEGCGPLCCSRRGACLSVGLLLLLAMLAALLAVA
 TILGHPPTPEHSCVTLANRTGFLCHDRRHCI PAHAVCDGIRTCPHGEDEAESLCRDVPSLPSFLLTT
 CGDPDSWIYSDQKCDGVNNGDCSDELSPVTTCPPCGPGWWRCSPTVFKYCSCVPRDLCRDSVQHCSDWS
 DEYSCPGP


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001081272

ORF Size: 654 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001081272.2
RefSeq Size:	947 bp
RefSeq ORF:	657 bp
Locus ID:	546840
Cytogenetics:	4 C7