

## Product datasheet for **MG221336**

### **Klk12 (NM\_027097) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Klk12 (NM\_027097) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Klk12  
**Synonyms:** 2310008B01Rik; AI324852; KLK-L5  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG221336 representing NM\_027097  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGAGGTTCAGCATTCTCTTGCTCCTGTGTGCTGTGGGTCTCAGTCAGGCAGACAGGGAGAAGATCTATA  
 ATGGTGTGAGTGTGTGAAGAACTCACAGCCTTGGCAGGTGGGGCTGTTCCATGGCAAATACCTGCGCTG  
 TGGCGGTGCTCTGTGGACCGCAAGTGGTCTTACGGCTGCTCACTGCCGTGACAAGTACGTGGTGCGC  
 CTTGGGGAACACAGCCTACCAAGCTGGACTGGACAGAACAGTTGAGACACACCACCTTCTCCATAACAC  
 ACCCCAGTACCAGGGCGCCTATCAGAACCATGAGCATGATCTTCGACTCCTGCGGCTGAACAGACCCAT  
 CCACCTGACCCGTGCTGTCCGGCCCGTGGCCCTGCCAGTTCTGTGTAACACAGGGGCCATGTGTCAT  
 GTCTCAGGATGGGTAACAACAAGCCATGGGACCCATTCAGACCCGCTCCAGTGCCTCAACCTCT  
 CCACTGTCTCAAATGAGACGTGCCGGGCTGTGTTTCTGGAAGAGTGACGGAGAATATGTTGTGTGCAGG  
 TGGAGAAGCCGGGAAGGATGCCTGTGAGGGTACTCTGGAGGCCCCCTGGTGTGTGGAGGGGTTCTTCAA  
 GGTCTGTTTCTGGGATCTGTTGGCCTTGTGGTCAAAAAGGTATTCCAGGAGTCTACACGAAGGTCT  
 GCAAATACACGGACTGGATCAGAATAGTCATCAGGAATAAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG221336 representing NM\_027097  
 Red=Cloning site Green=Tags(s)

MRFSILLLLCAVGLSQADREKIYNGVECVKNSQPWQVGLFHGKYLRCGGVLVDRKWVLTAAHCRDKYVVR  
 LGEHSLTKLDWTEQLRHITFSITHPSYQGAYQNHEHDLRLRLNRP IHL TRAVRPVALPSSCVTTGAMCH  
 VSWGTTNKPWDPFPDRLQCLNLSTVSNETCRAVFPGRVTENMLCAGGEAGKDACQDSSGGPLVCGGVLQ  
 GLVSWGVSVPCCGQKIGIPGVYTKVCKYTDWIRIVIRNN

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_027097

**ORF Size:** 741 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_027097.1</a> , <a href="#">NP_081373.1</a>
<b>RefSeq Size:</b>	930 bp
<b>RefSeq ORF:</b>	744 bp
<b>Locus ID:</b>	69511
<b>Cytogenetics:</b>	7 28.26 cM