

Product datasheet for **MG221349**

Rnf150 (NM_177378) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf150 (NM_177378) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rnf150
Synonyms:	A630007N06Rik; C030044C12Rik; Greul5; mKIAA1214
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG221349 representing NM_177378 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAATGTCTCTCATCCAAGCGTGCCGAGTCTGGCTCTCTCAACATGGCTGCTTTCCTTTTGTTCG
TGCATCTGCTCTGCCTGGACTTCACCGTGGCGGAGAAGGAGGAATGGTACACCGCCTTCGTGAACATCAC
CTATTTGGAGCCCGAGCCCGGGCCGCGGTGGCGGCAGCGGTGGCGCGCGGAGCTGCACACGGAGAAG
AGCGAGTGCAGGCGATACGGGGAGCACTCGCCCAAGCAGGACGCGCGCGGAGGTGGTATGGCCAGCT
CGGCCCAAGACCGCTGGCTTGCACCCCAACACCAAGTTCGCCGCGCCGCCACGGCAAGCACTGGAT
AGCCCTCATCCCCAAGGCAACTGCACGTACAGGATAAGATCCGGAACGCGTTCCTGCAGAACGCTCC
GCTGTGGTTCATCTTAACGTGGGCTCCAACACAAACGAGACCATCACCATGCCCCACGCAGGTGTGGAAG
ACATCGTGGCTATAATGATTCTGAGCCCAAGGGAAGGAGATAGTGAGCCTGCTGGAAAGGAACATCAC
CGTGACGATGTACATCACCATCGGCACCCGGAACTTGCAAAAATATGTGAGCCGCACGTCGGTTGTGTTT
GTCTCCATCTCCTTATTGTCTGATGATCATTCTCTGGCATGGCTGGTCTTTTATTATATCCAGAGGT
TTCGGTATGCAATGCCAGGGACAGGAACAGCGCCGACTGGGAGATGCAGCAAAGAAAGCCATCAGCAA
GCTCCAGGTGAGGACATCAGGAAGGGTGACAAGGAAACAGAGTCCGATTTTGATAACTGTGCAGTCTGC
ATAGAGGGGTACAAGCCCAACGATGTGGTCAGGATCCTGCCCTGCCGCGCATCTTTCCACAAGTCCCTGTG
TTGACCCCTGGCTTCTAGACCATCGCACCTGCCCATGTGCAAGATGAACATTCTTAAGGCCCTAGGGAT
CCCGCCCAATGCTGACTGCATGGATGACTTGCCCATTTGATTTTGGAGGCTCTCTAGGAGGCCACCAACC
AATCAAATCACAGGAGCCAGCGACACAACGTCAATGAGAGCTCTGTCACTCTGGACCCTGCTGTGCGGA
CAGTGGGAGCCCTGCAGGTGGTCCAGGACCCAGACCCTGCTCCTCAGGAGGGAGAAGCCATCTTTACAAC
TAACAGTGGTCAAGAGCCAGCACTGAGCAGTACTCTGACATCTCATTAAATTATGGCACTGGAAGTGGGC
CTGTCTGATGTAGAATTTCCACTGACCAAGACTGCGAGGAGGTGAAATCT

AGCGGACCGACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

ORF Size:	1311 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_177378.4 , NP_796352.2
RefSeq Size:	9685 bp
RefSeq ORF:	1314 bp
Locus ID:	330812
UniProt ID:	Q5DTZ6
Cytogenetics:	8 C2