

## Product datasheet for **MG221489**

### H4c2 (NM\_178193) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** H4c2 (NM\_178193) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** H4c2  
**Synonyms:** H4c1; H4c3; H4c4; H4c6; H4c8; H4c9; H4c11; H4c12; H4c14; H4f16; Hist1h; Hist1h4b  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG221489 representing NM\_178193  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTCTGGTCGCGGCAAGGGAGGAAAAGGCCTGGGCAAAGGTGGCGCTAAGCGCCATCGCAAAGTGTCTGC  
 GCGACAACATCCAGGGCATCACCAAGCCCGCCATCCGCCGCTGGCCCGCGGGGGAGTGAAGCGCAT  
 CTCCGGCCTCATCTACGAGGAGACCCGGTGTGCTGAAGGTGTTCTGGAGAACGTGATCCGCGACGCC  
 GTCACCTACACGGAGCACGCCAAGCGAAGACCGTCACCGCCATGGACGTGGTCTACGCGCTCAAGCGCC  
 AGGGCCGCACCCTCTACGATTCTGGGGT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG221489 representing NM\_178193  
 Red=Cloning site Green=Tags(s)  
 MSGRGKGGKGLGKGGAKRHRKVLRDNIQGITKPAIRRLARRGGVKRISGLIYEETRGVLKVFLENVIRDA  
 VTYTEHAKRKTVTAMDVVYALKRQRTLYGFGG

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)



---

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_178193.2</a> , <a href="#">NP_835500.1</a>
<b>RefSeq Size:</b>	312 bp
<b>RefSeq ORF:</b>	312 bp
<b>Locus ID:</b>	326620
<b>UniProt ID:</b>	<a href="#">P62806</a>
<b>Cytogenetics:</b>	13 A3.1
<b>Gene Summary:</b>	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H4 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. [provided by RefSeq, Aug 2015]