

## Product datasheet for **MG221530**

### **Kctd17 (NM\_001081367) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kctd17 (NM\_001081367) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Kctd17  
**Synonyms:** 2900008M13Rik; AA414907; AW742389; N28155  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG221530 representing NM\_001081367  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGACAACGCGGCCGGCGATGAGGATGGAGGCCGGGAGGCAGCGCCACCGGTGGGGCGGGCGGCC  
GCCCCGGGGCGGCTGGGGCAAGTGGGTGCGGCTGAACGTGGGGGCACGGTGTTCTGACCACCCGCCA  
GACGCTGTGTCGGGAGCAGAAGTCCTTCTCAGCCGCTGTGTCAGGGCGAAGAGCTTCAATCGGACAGG  
GATGAGACCGGGCCTACCTCATTGACCGTGACCCACCTACTTCGGGCCATCCTGAATTTCTCAGGC  
ATGAAAAATTGGTCTGGACAAGGACATGGCCGAGGAGGGGTCTGGAGGAAGCAGAGTTCTACAACAT  
TGGTCTCTGATTGGAATCATCAAAGACAGGATGGAGGAGAAAGACTACACAGTGGCCCAGGTGCCACCT  
AAGCATGTATACCGAGTCTGCAGTGCCAAGAGGAGGAGCTCACACAGATGGTATCCACTATGTCCGACG  
GCTGGCGCTTCGAGCAGCTGGTGAACATCGGCTCCTCTACAACATATGGGAGTGAGGACCAGGCCGAGTT  
CCTGTGTGGTGTCTAAGGAGCTCCACAGTCCCCACACGGCCTGAGCTCCGAGTCCACTCGCAAGGCC  
AAGAGCACGGACGAGCAGCTGGAGGAGCAGCGCGGCAGGAGGTGGAGGAGTGGAGGTGGCCCAGGTGC  
AGGTGGAGGCGGATGCCAGGAGAAAGCTATTACAAGCCAGAGGCACCCGGATGTGAGCCCCAGATCAC  
CTCCAGGGACTTGGGGTCCGATCTGAAATCCTTTATTTTTGTACCACGGGTGGGCCCTGGCCCAGG  
AGGAAGACAAGCCCTCTACCTTGTCAGTGCTGCCTGCATCTGCCGGGC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG221530 representing NM\_001081367  
 Red=Cloning site Green=Tags(s)

MQTTRPAMRMEAGEAAPPVGGGRPGGGWGKWRNLNVGGTVFLTTRQTLCREQKSFLSRLCQGEELQSDR  
 DETGAYLIDRDPTYFGPILNFLRHGKLVLDKDMAEEGVLEAEFYNIIGPLIRIIKDRMEEKDYTVAVQVPP  
 KHVYRVLQCQEEELTQMVSTMSDGWRFEQLVNIIGSSYNYGSEDQAEFLCVVSKELHSSPHGLSSESTRKA  
 KSTDEQLLEEQRREQVEVEVAQVQVEADAQEKAITSQRHPDVSPQITSRDLGRFSEILYFCTTGWAPGPR  
 RKTSP LPCQCLHLPG

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_001081367

**ORF Size:** 888 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001081367.1</a> , <a href="#">NP_001074836.2</a>
<b>RefSeq Size:</b>	1574 bp
<b>RefSeq ORF:</b>	891 bp
<b>Locus ID:</b>	72844
<b>Cytogenetics:</b>	15 E1