

Product datasheet for **MG221576**

Stx7 (NM_016797) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Stx7 (NM_016797) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Stx7
Synonyms:	AI315064; AI317144; AW107239; Syn7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG221576 representing NM_016797 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTACACTCCGGGGATTGGTGGGGACTCTGCGCAGCTAGCACAACGCATCTCCTCTAACATCCAAA
AGATCACACAATGCTCCGTGAAATACAAAGGACTCTGAATCAACTTGGAAACACCTCAAGACTCACCGGA
ACTGAGGCAACAGCTGCAGCAGAAGCAGCAGTATACTAACCAACTTGCAGAAAGAAACCGATAAGTACATT
AAAGAGTTTGGCTCTCTGCCACCACCCCAAGTGAGCAGCGTCAGAGAAAGATACAGAAAGACCGATTAG
TGGCAGAGTTCACCACATCGCTGACAAACTCCAGAAGGCCAGAGGCAGGCTGCTGAGAGAGAGAAAGA
GTTTCGTTGCTCGAGTGCAGCCAGCTCCAGGGTATCGGGTGGTTTTCTGAAGACAGCTCAAAGAAAAG
AATCTTGTATCCTGGAAAGCCAAACACAGCCTCAAGTGCAAGTCCAAGATGAAGAAATCACAGAGGATG
ACCTCCGACTCATTATGAGAGAGAGTCTTCAATCAGGCAGCTGGAAGCTGATATTATGGACATTAATGA
AATATTTAAAGACTTGGGGATGATGATTCATGAACAAGGCATATGATTGACAGCATAGAAGCCAATGTA
GAAAGTGCAGGAAAGTTCACGTTCAACAGGCAAACAGCAGCTGTCAAGAGCGGCAGACTATCAGCGCAAT
CCAGGAAAACCTCTGCATTATCATTATCTCTCGTGGTCCGGAATCGTGATCATCTGTCTCATCGTATG
GGACTGAAAGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG221576 representing NM_016797
Red=Cloning site Green=Tags(s)

MSYTPGIGGDSAQLAQRISSENIQKITQCSVEIQRTLNLGTPQDSPELRQQLQQKQYTNQLAKETDKYI
 KEFGSLPTTPSEQRQRKIQKDRLVAEFTTSLTNFQKAQRQAAEREKEFVARVRASSRVSGGFPESSKEK
 NLVSWESQTQPQVQVQDEEITEDDLRLIHERESSIRQLEADIMDINEIFKDLGMMIHEQGDMIDSIEANV
 ESAEVHVQQANQQLSRAADYQRKSRKTLCTIIIFILVVGIVIIICLIWGLKG

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_016797

ORF Size: 783 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016797.5
RefSeq Size:	2500 bp
RefSeq ORF:	786 bp
Locus ID:	53331
Cytogenetics:	10 A4