

## Product datasheet for **MG221589**

### **Klk6 (NM\_001164697) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Klk6 (NM\_001164697) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Klk6  
**Synonyms:** AI849898; BSP; Bssp; Klk29; MSP; neurosin; Prss9; Prss18  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG221589 representing NM\_001164697  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCCATGAAGATGCTGACAAATGAAGATGCTGGCCCTGTGCTTGGTTCTTGCTAAATCAGCCTGGTCGG  
AGGAACAGGAGAAGGTGGTTCATGGAGGCCGTGTTTGAAGGACTCCACCCTTCCAGGCTGCCCTCTA  
CACCTCAGGTCACCTGTGTGGTGGGGTCTCATTGACCCACAGTGGGTGCTGACAGCTGCCACTGC  
AAAAAACCGAATCTGCAGGTGATCTTGGGAAACACAACCTACGGCAAACAGAGACTTCCAAAGGCAAA  
TCTCAGTGACAGGACTATTGTCCATCCCCGCTACAACCCTGAAACCCACGACAATGACATCATGATGGT  
GCATCTGAAAAATCCAGTCAAATCTCTAAAAAGATCCAGCCTCTGCCCTTGAAGAATGACTGCTCTGAG  
GAGAATCCCAACTGCCAGATCCTGGGCTGGGGCAAGATGGAAAATGGTGACTTCCAGATACCATTCAGT  
GTGCTGATGTCCATCTGGTGCCCGGGAGCAGTGTGAGCGTGCCTACCCTGGCAAGATCACCCAGAGCAT  
GGTGTGCGCAGGCGACATGAAAGAAGGCAACGATTCTGTGAGGGTGATTCTGGAGGTCCCTAGTATGT  
GGGGTTCGCCTCCGAGGGCTCGTGTGATGGGGTGACATGCCCTGTGGATCAAAGGAGAAGCCAGGAGTTT  
ACACCGATGTCTGCACTCATATCAGATGGATCCAAAACATCCTCAGAAACAAGTGCGT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG221589 representing NM\_001164697  
Red=Cloning site Green=Tags(s)

MPMKMLTMKMLALCLVLAKSAWSEEQEKVVHGGPCLKDSHPFQAALYTSGHLLCGGVLIDPQWVLTAAHC  
 KKP NLQVILGKHNL RQTETFRQRI SVDRTIVHPRYNPETHDNDIMMVHLKNPVKFSKKIQPLPLKND CSE  
 ENPNCQILGWGMENGD FPD TIQCADVHLVPREQCERAYPGKITQSMVCAGDMKEGNDSCQGDSSGGPLVC  
 GGRLRGLVSWGDMPCGSKEKPGVYTDVCTHIRWIQNILRNKWL

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001164697

**ORF Size:** 759 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001164697.1</a> , <a href="#">NP_001158169.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 1346 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 762 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 19144                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Cytogenetics:</b>          | 7 28.28 cM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |