

## Product datasheet for **MG221692**

### Slc26a3 (NM\_021353) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Slc26a3 (NM_021353) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Slc26a3
Synonyms:	9030623B18Rik; 9130013M11Rik; AV376035; D; Dra
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG221692 representing NM\_021353  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGATCGAAGCCATAGGGAATCAGTATGTGGTGGCCAGACCGGTGTACTCCACAAAGACTTTTGGAGAGG  
AGTTTAAAGAAGACGCACAGACATCAAAAGACATTCCTGGACCATCTCAAAGGGTGTGTAGCTGCTCCTC  
ACAGAAGGCCAAGAAAATCGCCCTGTCTTTGTTCCCATAGCATCTTGGTTGCCAGCCTACAAGATAAAG  
GAATGGCTTCTCAGTGACATCGTCTCTGGCATCAGCACCGGGCTGGTGGCTGTACTGCAAGGTTTAGCAT  
TTGCTCTGCTGGTCAACATCCCTCCAGCCTACGGGTTGTACGCAGCCTTTTTCCAGTCATAACCTACTT  
TTTCTGGGCACATCTAGACACATATCTGTGGTCCATTTCCAGTTCTGAGTATGATGGTGGGAGTTGTC  
GTTACAAGAGTGGTCTCGGACCCCAATGCTTCTCAGAGTTGCCTCAAGCTCCACAGAAAATGACTCGT  
TCATAGAGGAGAAAGTAATGGTGGCTGCATCAGTAACTGTTCTTCTGGAATCATTAGTTGCTCCTGGG  
GGTTCTGCAGGTCGGCTTTGTGGTGATACCTATCCGAGTCCCTAATCAGTGGCTTACCACGGCTGCT  
GCCATTACGTTCTGGTTTCTCAACTGAAATTCATGCTACAGCTGCCAGTCCCAGCATACAGCGACCCAT  
TCTCTATATTCAAAGTCCTGGAGTCGGTCTTACACAAAATTCAGAAGACGAACATCGCAGACCTTGTGAC  
ATCTGTGATTATCTAGTGGTCGTGTTGTCTTTAAAGAAATCAATCAGCGTTACAGAAGCAAACCTCCG  
GTGCCATCCCATCGAACTCATGACTGTAATTGCCACAGGTGTGTCCTATGGCTGTAACCTTGAAG  
ACAGGTTTGGCGTGGCCGTGGTTGGGAACATGAGTCTTGGATTTCAACCCCTATTACACCCAGCGTGGA  
AGTTTTCAAGACACCATAGGAGACAGCTTCGGCATCGCGATCGTCGGCTTTGCCGTGGCCTTTCCGTT  
GCCAGCGTGTACTCCCTCAAATACGATTATCCGATTGATGGCAATCAGGAGTTAATTGCCTGGGAGTGA  
GCAACATATTCAGTGGACATTCAAAGGATTTGCAGGGAGCACAGCCCTTCCAGATCAGGGGTTTCAGGA  
GAGCACGGGAGGCAAAACACAGGTTGCGGGGCTTCTCTCAGCTGTCATCGTGTGATCGTCATAGTTGCC  
ATTGGGTTTCTTCTGCAGCCGCTACAAAAGTCTGTCTGGCAGCGCTAGCACTCGAAAATTAAGGGAA  
TGCTGATGCAGTTTGTGAAATAGGCAGATTGTGGAAAAAGGATAAATATGATTGTTAATTTGGATCAT  
GACCTTCATCTTCGCCATTGTCTTGGACTCGGATTAGGCCTGGCGGCAAGTGTAGCATTCAACTCCTA  
ACTATTGTGTTCCAGGACCAATTTCAAAGTGCAGCACATTGGCTAATGTTGGAAGGAGCAACATCTACA  
AGAATAAAAAGAATTATGCCGAAGTGTATGAGCCTGAAGGAGTAAAAATTTTAGATGTCCATCTCCAAT  
CTACTTCGAAACATTGGTTTCTTTAAGCAGAACTTATTGATGCTGTTGGCTTTAGCCCACTTCGAATT  
CTACGGAAGCGCAACAAAGCTTTGAAGAAAATCCGAAAACGCAGAAACGAGGTCTGATACAGATGACAC  
CTAAAGGTTTATATGCACCTCTGACGGTTTCAAAGACTCCGACGAAGAATTGGACAACAACCAGATAGA  
AGAAGTTGATCAGCCGATCAATACCACAGATCTGCCCTTTGACATAGACTGGAACGGGGATCTTCCCTC  
AACATCACCATCCCCAAAATCAGCCTCCACAGCCTCATCTCGATTTTTAGCAGTGTCTTTCCTTGACG  
TCTCCTCAATGAGGGGCTCAGAACGATTTTACAAGAAATTTAGGATCAAGGTGGACGTATATATTGT  
GGGAACCGATGATGATTTCAATGACAACTTGTCTGGTGTGAATCTTCGATGATGAAGTCACAGACTCA  
ATATTTTCTTAAACAATCCATGATGCTATTTTGCATATTTGATGAAGAAGGACTACAGCACTTCAAAGT  
TTAATTCAGCCAGGAAAAAGAAAGAAATTTGACTTTACCATAAATACAAATGGAGGACTACGTAATCG  
GGAATGTCAGGTACCAGTTGAAACGAAATTC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG221692 representing NM\_021353  
 Red=Cloning site Green=Tags(s)


MIEAIGNQYVVARPVYSTKTFGEFEKFKTHRHHTFLDHLKGCCSCSSQKAKKIALSLFPIASWLPAYKIK  
 EWLLSDIVSGISTGLVAVLQGLAFALLVNIPPAYGLYAAFFPVITYFFLGTSRHISVGPFPVLSMMVGVV  
 VTRVSDPNASSELSSSSTENDSFIEEKVMAASVTVLSGIIQLLLGVLQVGFVVIYLSSELISGFTTAA  
 AIHVLVSQKFMQLQLPVPAYSDFPFIFKVLSEVFTQIQKTNIAIDLVTSVIILVVVFVFKINQRYRSKLP  
 VPIPIELIMTVIATGVSYGCFEDRFGVAVVGNMSLGFQPPITPSVEVFQDTIGDSFGIAIVGFVAVFSV  
 ASVYSLKYDYPIDGNQELIALGVSNIIFTGAFKGFAGSTALSRSGVQESTGGKTQVAGLLSAVIVLIVIVA  
 IGFLQPLQKSVLAALALGNLKGMLMQFAEIGRLWKKDKYDCLIWIMTFIFAIVLGLGLAASVAFQLL  
 TIVFRTQFPKCSTLANVGRSNIYKNKNYAIEVYEPEGVKIFRCPSPIYFANIGFFKQKLIDAVGFSPLRI  
 LRKRKALKKKIRKLQKRGLIQMTPKGFICTSDGFKDSDEELDNNQIEELDQPINTDLPFDIDWNGDPL  
 NITIPKISLHSLILDFSAVSFLDVSSMRGLRTILQEFIRIKVDVYIVGTDDDFIDKLARCEFFDDEVTD  
 SIFFLTIDHAILHILMKDYSTSKFNSSQEKERKDFDFTINTNGGLRNRECQVPVETKF

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_021353

ORF Size: 2908 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021353.3](#), [NP\\_067328.1](#)

**RefSeq ORF:** 2274 bp

**Locus ID:** 13487

**UniProt ID:** [Q9WVC8](#)

**Cytogenetics:** 12 13.45 cM

**Gene Summary:** This gene encodes a member of the solute carrier/sulfate transporter family. The encoded protein is predominantly expressed in the intestine where it is essential for chloride absorption. Disruption of this gene results in chloride-rich diarrhea and compensatory up-regulation of ion-absorbing transporters. [provided by RefSeq, Dec 2012]