

Product datasheet for **MG221884**

Lrrfip1 (NM_001111311) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrfip1 (NM_001111311) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Lrrfip1
Synonyms:	AU024550; Fliiap1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG221884 representing NM_001111311
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACCAGCCAGAGGGCGCGCAAAACAAGGAAATCGACTGTCTGAGTCCAGAGGCGCAAAGACTGGCGG
AGGCCAGGCTCGCAGCGAAGCGTGCAGCCCGAGCCGAGGCTCGAGAGATCCGGATGAAGGAGCTGGAACG
GCAGCAGAAAGAGGTAGAAGAAAGACCAGACAAAGACTTTGCTGAGAAGGGCTCCCGTAACATGCCTAGC
CTGTCTGCAGCCACACTGGCATCTCTGGGTGGGACTTCCTCCCGAGAGGGAGTGGAGATACCTCCATCT
CCATGGACACCGAAGCATCCATCAGGAAATCAAGGACTCGTGGCCGAAGTTGAGGAGAAGTACAAGAA
AGCTATGGTCTCCAACGCCAGCTAGACAATGAGAAGACCAACTTCATGTACCAGGTGGACACACTGAAG
GACATGCTGCTGGAGCTGGAAGAGCAGCTGGCCGAGTCTCAGCGCAATACGAGGAGAAGAACAAGGAGT
TCGAGAGGGAGAAGCATGCCACAGCATCTGCAGTTCAGTTTGCCGAAGTGAAGAGGGCCCTGCGGCA
AAGGGAAGAGATGCTCGAGAAACATGGAATAATCCTAAATTCAGAAATAGCCACCAACGGAGAGACTTCA
GACACAGTAAATGACGTTGGGTATCAAGCCCTACGAAGATAACAAAAGAAGAGCTGAATGCCCTCAAGT
CGGCCGGGGAGGGGACCCTAGGAAAAGCCAAAAGAGTGGAGGTGAAAAAGGAAATTTGGGAGAAAGTGGG
GCAAAGAGAAACCTTACAGAATTCTGAGCAAGAACAGCCCAAACCAACACAGGAAAGGATTGTGTGGAC
AGAGGGGTGTCACATCCTGGTGAGAAGGCTGAGAACAGAGACCCGCTGAAGACAGTGCCTGTCCCCGG
GACCACTAGCAGGGGCCAAATGTGAGCAACAGGTTGAGAGCAAGATCAAGAGAATACTTCTGACCTCAA
AAACTCAGAGCAGATCGAGTCACACAAGGTCACAAACAAGTCAGATAGTAGGGCCAGTAACTCCCCTGAA
CAGTCCAGCTGCCTGGAAGGCTAGACAGTGAAGTGCCAGGGCCCACTGAGGATTTGAAAACAGATTTGG
GGAAAGGTAGCTTTGAACCATGTCCAGACTACATTTTGGGGCAAACAGCTGAAATTTGACAAGGTGACCTG
TACAGACTCAAGAGGACACAGGTGGAACCAAAAGGAGGATGAGGTACAGGCAGGGGACACCACAGTTGAG
GACCAAGTAGGGACAGTGGCTTCTGGACCTGCAAAGCAAAGCAAAGGCACCGAGAATCATGGTGAAAGTT
GTCTGAAGGATGGGCTGGGACAGAGCTCAGAGAGAGAACTCACCAAGAAGTAGCTGAGCCTGAGGAAGC
CATCGTTCAGATCCCACAGGCCGGTGGAGAGAACTATTACCAAGGCAGATGATGCCGAAGGAAGGGAT
GAGAAGCCATCCAGGCAGAAGCCAGGCAAGCCCTGGAGCTCCCATAAATCAGAGCGGCCATCAGGACA
CAACAGGACCAGGCAGCACCATGCCACAGCGCACACCCTCATGCAAAAGAGCGCAAGAAACAAGGGAA
GAGCGAGCAGCAGGCAGAGGCCCTGGATTACCCGCAAGAAGACGAAGAACAAGAAGAAGAACAAG
AAGAAGAAGCTGCAACTCCCAGGAAACCTGCAGAGATGCCAATGAAGAGTTGAACTGTCAGGACCCAG
ATGTGGGTGACATGGAGGAGGAAGAGCGGCTCCAGGCTACTGACAAGAAACAGGCATCAGGAAGCCAGA
GCAGAAAATCAGAGCGGGAAGCAGAGAACCTGTGGAAGATCCTCAAAGTGGGTCGAGTGGAAAACAGAAC
AAAGTGGAAAGAGGATGGTCCCACGGAAGGTCCCACAGACATCTTGGATCAGAACAGTCCCAGTGTGAAG
ACAGGGAGATCTCCCAGTGGGAGAGAAAGGTCCCCAATGCGATACCTCTCAGATAGGTAGTGAAGAGGG
ACATGTGACCTCCAACACGGTGGTCAGGCAGTTGAGAACCACAACCTAGATAACAGTGACCTGTCAGGG
CAGCTGGAGGGCTTCAATTCAGAAAGTGGAGGGCAGGCAAGAGAGGAGTTGGGAATAGCAAGAGTAAGG
AGGACTGCACCATGTCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG221884 representing NM_001111311
 Red=Cloning site Green=Tags(s)


MTSPEGAQNKEIDCLSPEAQRLEARLAAKRAARAEAREIRMKELERQQKEVEERPKDFAEKGSRNMP
 LSAATLASLGGTSSRRGSGDTSISMDEASIREIKDSLAEVEEKYKKAMVSNAQLDNEKTNFMYQVDTLK
 DMLLEEEQLAESQRQYEEKNEFEREKHAHSILQFQFAEVKEALRQREEMLEKHGIIILNSEIATNGETS
 DTVNDVGYQAPTKITKEELNALKSAGEGTLGKAKEVEVKKEIVEKVGQRETLQNSEQEQPKPNTGKDCVD
 RGVSHPGKAENQRPAEDSALSPLAGAKCEQQVQSQDQENTSDLKNSEQIESHKVTNKSDSRASNSPE
 QSSCLEGLDSEVPGPTEDLKTDLGKGFEPDYLILGQTAEIDKVTCTDSRGTGGNQREDEVQAGDTTVE
 DQVGTVASGPAKQSKGTENHGESCLKDGLGQSSERELTQEVAPEEAIIVQIPQAGGENTITKADDAEGRD
 EKPIQAEAQASPGAPINQSGHQDTTGPGSTDAQRTPPHAKERKKQKSEQQAEALDSPQKKTKNKKKKNK
 KKAATPAETCRDANEELNCQDPDVGDMEEERLQVTDKKQASGSPEQKIRAGSREPVEDPQSGSSGKQN
 KVEEDGPTGPTDILDQNSPQCEDREISPVGEGKGPQCDTSQIGSEEGHVT SQHGGQAVENHNLDNSDL SG
 QLEGFNSEGGQAREEVGNSKSKEDCTMS

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2138_b12.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:




ACCN: NM_001111311

ORF Size: 2187 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001111311.1</u> , <u>NP_001104781.1</u>
RefSeq Size:	3662 bp
RefSeq ORF:	2190 bp
Locus ID:	16978
UniProt ID:	<u>Q3UZ39</u>
Cytogenetics:	1 D
Gene Summary:	Transcriptional repressor which preferentially binds to the GC-rich consensus sequence (5'-AGCCCCGGCG-3') and may regulate expression of TNF, EGFR and PDGFA. May control smooth muscle cells proliferation following artery injury through PDGFA repression. May also bind double-stranded RNA (By similarity). Positively regulates Toll-like receptor (TLR) signaling in response to agonist probably by competing with the negative FLII regulator for MYD88-binding (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG221884