

Product datasheet for **MG221922**

Gal3st1 (NM_016922) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gal3st1 (NM_016922) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gal3st1
Synonyms:	Cst; Gcst
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG221922 representing NM_016922 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTCTGCTGCCAAAGAAGCCCTGCAAGTCCAAGGCCAAGGGTTACTGCTGGGAGCTCTTTACCA
GTTTCCTATTGCTGCTGTACTCCTATGTGGTACCTCCACTGTATCCCAACATGGCCTTCACGACCTCAGA
GGCTGCAGCACCTGCTCCCTATTCCCAATGAGCCAGTGGCAGCCACTCCCGCCAACGGCTCAGCAGGA
GGGTGTCAACCTCGCGAGATATTGTGTTTCATGAAGACGCACAAGACCGCCAGCAGCACACTGCTCAACA
TCCTGTTCCGCTTCGGCCAGAAGCACGAGCTCAAGTTCGCTTTCCCAATGGCCGAAACGACTTCCACTA
CCCGTCGTAATTTGCACGAAGCCTGGTGCAGGACTACCGCCTGGGGCGTGTTCACATCATCTGCAAC
CACATGCGCTTCCACTATGAGGAAGTGCAGCGGGCTGGTGCAGCGGGTGCACCTTCATCACCGTATCC
GCGACCCTGCTCGTCTCTTCGAGTCTTCCCTTCCACTACTTTGGATCCGTGGTGCCGCTCACCTGGAAGCT
GTCGAGCCGCGACAAGCTGGCCGAGTTCCTGCAGGACCCGATCGCTACTACGACCCGAGCAGCTACAAC
GCGCACTACCTCCGCAACCTGCTTCTTCGATCTGGGCTATGATAGTAGCCTGGACCCGGCCAGCCCGC
GCGTGCAGGAGCACATCCTGGAGGTGGAGCGCCGCTTCCACTTGGTGCTACTGCAGGAGTACTTCGACGA
ATCCCTGGTGTGCTCCGAGAGCTGCTGTGCTGGGACCTGGAGGACGTGCTGATTTCAAGCTCAATGCC
CGCGCGACTCGCCGGTTCGGCGCTCTCTGGCAGCTGTACCAGCCGAGCCACCGCTGGAACCTGCTGG
ACGTGCGCCTCTACCGTCACTTCAACGCCAGCTTCTGGCGAAGGTGGAAGCTTTCGGGCGTGAGCGCAT
GGCGCGGAGGTAGCCGAGCTGCGCCAAGCGAAGCAGCAGCATGCGCCACATCTGCATAGATGGCGGCCAA
GCAGTTGGTGCTGAGGCCATCCAGGACTCGGCCATGCAGCCCTGGCAACCCCTGGGCATCAAGTCCATCC
TGGGTTACAACCTCAAGAAGAGCATCGGGCCACAACACGAGCAGCTCTGCCGGCGCATGCTCACGCCGA
AATCCAGTACCTGTCTGACCTCGGTGCCAATCTCTGGGTACCAAGCTCTGGAAGTTCCTTAGGGACTTT
CTAAGGTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >MG221922 representing NM_016922
Red=Cloning site Green=Tags(s)

MTLLPKKPCKSKAKGLLLGALFTSFLLLLYSYVPPLYPNMAFTTSEAAAPCSPIPNEPVAATPANGSAG
 GCQPRRDIVFMKTHKTASSTLLNILFRFGQKHELKFAFPNGRNDFHYPYSYFARSLVQDYRPGACFNICN
 HMRFHYYEEVRLVRPGATFITVIRDPARLFESSHFYFGSVVPLTWKLSRDKLAEFLQDPDRYYDPSSYN
 AHYLRNLLFFDLGYDSSLDPASPRVQEHILEVERRFHLVLLQEYFDESLVLLRELLCWDLEDVLYFKLNA
 RRDSPVPRLSGELYRRATAWNLLDVRLYRHFNASFWRKVEAFGRERMAREVAELRQANEHMRHICIDGGQ
 AVGAEAIQDSAMQPWQPLGIKSLILGYNLKKSIGPQHEQLCRRMLTPEIQYLSDLGANLWVTKLWKFLRDF
 LRW

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_016922


ORF Size: 1269 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016922.3 , NP_058618.2
RefSeq Size:	1969 bp
RefSeq ORF:	1272 bp
Locus ID:	53897
UniProt ID:	Q9JHE4
Cytogenetics:	11 A1
Gene Summary:	Catalyzes the sulfation of membrane glycolipids. Seems to prefer beta-glycosides at the non-reducing termini of sugar chains attached to a lipid moiety. Catalyzes the synthesis of HSO3-3-galactosylceramide (sulfatide), a major lipid component of the myelin sheath and of HSO3-3-monogalactosylalkylacylglycerol (seminolipid), present in spermatocytes. Also acts on lactosylceramide, galactosyl 1-alkyl-2-sn-glycerol and galactosyl diacylglycerol (in vitro). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG221922