

Product datasheet for **MG221987**

Aamp (NM_146110) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aamp (NM_146110) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Aamp
Synonyms:	Aamp-rs; AU040907
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG221987 representing NM_146110 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCCGAATCGGAGAGCGGAGCGGCCGCTGACACCCCGCCTCTGGAGACTCTGAGCTTTCACGGAG
ATGAAGAGATCATTGAGGTGGTAGAGCTGGATCCAGGCCCGCCGACCCGGATGATCTGGCCAGGAGAT
GGAAGACGTGGACTTTGAGGAAGAGGAAGAAGAAGAGGGCAATGACGAGGGCTGGGTCTTGAACCC
CAGGAAGGGTGGTGGTAGCATGGAGGCCCGGATGATAGTGAGGTACCTTTGCATTGCACTCAGCGT
CTGTGTTCTGTGTGAGCCTGGACCCAAAACACACCTTGGCAGTGACAGGGGGCAAGACGACAAGGC
CTTCGTGTGGAGGCTCGGTGATGGAGAGCTGCTCTTTGAGTGTGCAGGCCACAAAGACTCTGTGACATGC
GCTGGTTTCAGTCACGATTCCACCTAGTGCCACAGGGGACATGAGTGGTCTTTTAAAAGTATGGCAGG
TAGACACCAAGGAAGAAGTCTGGTCTTTTGAAGCAGGAGACCTGGAGTGGATGGAGTGGCACCCCTCGGGC
TCCTGTCTGTTGGCGGCACAGCCGATGGCAACACTTGGATGTGGAAGGTCCCGAATGGTACTGTAAG
ACCTTTCAGGGCCCAACTGCCCTGCAACCTGTGGCCCGCTCCTCCTGATGGGAAGAGAGCTGTGGTAG
GGTATGAAGATGGTACCATCAGGATCTGGGACCTCAAGCAAGGAAACCCATCCATGTAATAAAGGGAC
TGAGGGTCATCAGGGCCCTCTGACCTGTGTAGCCACCAACCAGGATGGCAGTCTGATCCTGACTGGCTCT
GTAGACTGCCAGGCCAAGCTGGTCACTGAGCCACCAGTGGCAAGGTGGTGGGTGTTCAGACTGAGACAG
TCGCCTCCAGCCAGCCTGGGAGAGGGGAGGAGAGTGAAGTCTAACTCGTGGAGTCTTGGGCTTCTG
CAGTGTATGCCTCTGGCCGCTGTTGGCTATCTGGACGGAACCTTGGCCATCTATGACCTGTCTACGCAG
ACACTCCGGCACCAATGTCAGCACCAGTCGGGCATTGTTGAGCTGCTGTGGGAGGCTGGCACTGCTGTGG
TGTACACCTGTAGCCTGGATGGTGTGTACGTCTCTGGGATGCCCGGACTGGCCGCTGCTTACTGACTA
CCGGGGCCACACTGCTGAGATCCTGGACTTTGCCCTCAGCAAAGATGCTTCCCTGGTGGTGACCACATCA
GGAGACCACAAAGCAAAGTATTTGTGTCCAGAGACCTGACCGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG221987 representing NM_146110
 Red=Cloning site Green=Tags(s)


```
MESESESGAAADTPPLETLSFHGDEEIEVVVELDPGPPDPDDLAQEMEDVDFEEEEEEEGNDEGWLEP
QEGVVGSMEGPDDSEVTFALHSASVFCVSLDPKTNLAVTGGEDDKAFVWRLGDGELLFECAGHKDSVTC
AGF SHDSTLVATGDMSSLKLVWQVDTKEEVWSFEAGDLEWMEWHPRAPVLLAGTADGNTWMWKVPNGDCK
TFQGPNCPATCGRVLPDGKRAVVGVEDGTIRIWDLKQGNPIHVLKGTTEGHQGPLTCVATNQDGSILTGS
VDCQAKLVSAATTKVVGVFRPETVASQPSLGEGESESNVSLGFCVSMPLAAVGYLDGTLAIYDLSTQ
TLRHQCQHQSGIVQLLWEAGTAVVYTCSLDGVVRLWDARTGRLLTDYRGHTAEILDFALSKDASLVVTTT
GDHKAKVFCVQRPDR
```

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_146110

ORF Size:	1305 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_146110.3 , NP_666222.2
RefSeq Size:	1817 bp
RefSeq ORF:	1308 bp
Locus ID:	227290
Cytogenetics:	1 38.53 cM